

Manual No.	21
Revision	March 17, 2023
Product No.	121
Lot No.	A17-2 B16-2 DR43-3

WORKING INSTRUCTION

HLA-A*/-B*/-DRB1* CTS-PCR-SSP TRAY KIT

LOCUS- AND LOT-SPECIFIC MANUAL

To be applied to the following product:

Product No.	Description
121	HLA-A*/-B*/-DRB1* CTS-PCR-SSP TRAY KIT  IVD CE 0197

1. Main differences

- **Between lot A17-2 (the current lot) and lot A17-1:**
The kit was updated to cover new alleles included in the IMGT/HLA Sequence Database of January 2023. Deleted and renamed alleles were taken into consideration.
- **Between lot B16-2 (the current lot) and lot B16-1:**
The kit was updated to cover new alleles included in the IMGT/HLA Sequence Database of January 2023. Deleted and renamed alleles were taken into consideration.
- **Between lot DR43-3 (the current lot) and the previous lot DR43-2:**
The kit was updated to cover new alleles included in the IMGT/HLA Sequence Database of January 2023. Deleted and renamed alleles were taken into consideration.

2. Introduction

- Intended use: This kit reveals a low/intermediate resolution typing of HLA-A*, -B* and -DRB1* (besides HLA-DRB3*/4*/5* low resolution) by the PCR-SSP method.
- Allele coverage: IMGT/HLA Sequence Database Release 3.51.0, January 2023 for HLA-A*, -B* and -DRB1*, except:
 - **HLA-A*** A*01:01:65-01:01:66/01:01:131/01:38/01:72/01:102/01:139/01:167/01:221/01:267/01:284/01:390/01:399/01:418N, A*02:644, A*11:01:35/11:01:57/11:01:104/11:23/11:25:03/11:31/11:35/11:40/11:45/11:50Q/11:78N/11:98/11:111-11:112/11:176/11:183/11:191/11:211/11:223/11:226/11:229/11:250/11:258/11:269/11:290/11:310N/11:346/11:380/11:389/11:392/11:396-11:397N/11:432N/11:435, A*24:129, A*25:01:09/25:01:17/25:81, A*26:01:59/26:01:65/26:60N/26:85/26:125/26:131/26:143-26:144/26:204:01:01-26:204:01:02, A*29:01:07/29:02:36/29:60/29:111, A*30:12:02/30:26/30:82/30:84/30:184Q/30:189/30:206, A*31:01:07/31:01:47/31:21-31:22/31:75/31:90/31:104/31:109/31:202N/31:204/31:210, A*32:01:07/32:01:12/32:01:48/32:28/32:33:02/32:34/32:66/32:71/32:75/32:94/32:158, A*33:01:05/33:03:26-33:03:27/33:03:31/33:03:56/33:03:58/33:18:02-33:19/33:30/33:32:0133:32:02/33:35/33:66/33:81/33:89/33:122/33:130/33:135/33:140N/33:152/33:158-33:159/33:175Q/33:188/33:215, A*34:24, A*68:01:11/68:01:25/68:02:07/68:48-68:50:02/68:137/68:184/68:218/68:230/68:234/68:236N/68:251N/68:258/68:286/68:293Q, A*74:16:01, A*80:02/80:04
 - **HLA-B*** B*08:01:32/08:12:03/08:26:01-08:26:03/08:50/08:94/08:112/08:170/08:181/08:232/08:263/08:272/08:297/08:299, B*13:22:01-13:22:02/13:58/13:86/13:97, B*14:01:12/14:02:03/14:02:06/14:08:01-14:08:02/14:21/14:42:01/14:53/14:55, B*15:162/15:253/15:396/15:553, B*27:05:23/27:75/27:140/27:189/27:239/27:250, B*37:01:12/37:05/37:22:01-37:22:03/37:29/37:95, B*38:73, B*39:01:11/39:06:05/39:24:02/39:30/39:34/39:36/39:80/39:82/39:107/39:128/39:136/39:156, B*44:02:13/44:03:37/44:03:51/44:83/44:134/44:201/44:445, B*51:45/51:81/51:148/51:220, B*53:38/53:44, B*55:58/55:121, B*57:82, B*58:27, B*59:08
 - **HLA-DRB1*** DRB1*01:01:31/01:02:12/01:18:01-01:18:02/01:22/01:33N/01:44:01-01:44:02/01:60/01:76/01:134, DRB1*04:11:04/04:156, DRB1*08:01:06/08:57/08:80, DRB1*09:09/09:23, DRB1*10:12, DRB1*11:02:03/11:04:13/11:86/11:94/11:110/11:123/11:156/11:167/11:172/11:185/11:187/11:192/11:208/11:294N, DRB1*12:31N/12:42/12:57/12:83, DRB1*13:195, DRB1*14:231N/14:247, DRB1*15:01:22/15:34, DRB1*16:50, DRB3*01:68/01:95/01:104/01:109, DRB3*02:26-02:27/02:47/02:58/02:66/02:91/02:111/02:131/02:143/02:160, DRB4*01:14N/01:62/01:88, DRB4*03:01N, DRB5*01:01:10/01:51/01:73/01:79Q/01:87/01:97/01:110These alleles are considered to be rare.
- This manual is only valid for **Lot No. A17-2 B16-2 DR43-3**.
- This manual should be used together with the Main Manual (General Information) which is the 'Working instruction for the CTS-PCR-SSP **TRAY and MINITRAY KITS**' (Manual No. 100A).

3. Content

1. Main differences	1
2. Introduction	2
3. Content.....	3
4. Kit Composition	4
5. Materials, Reagents and Equipment not supplied	4
6. Sample Requirements, PCR and Gel Electrophoresis	4
7. Result Evaluation	4
8. Interpretation Hints	4
9. Special notes	4
10. Troubleshooting	5
11. Precaution	5
12. Contact	6
13. Appendix	7
14. Certificates of Analysis	76
 Table 1: Allele specificities and sizes of the PCR products of each HLA-A* CTS-PCR-SSP primer mix (Lot No. A17-2) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023	7
 Table 2: Allele specificities and sizes of the PCR products of each HLA-B* CTS-PCR-SSP primer mix (Lot No. B16-2) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023	14
 Table 3: Allele specificities and sizes of the PCR products of each HLA-DRB1* CTS-PCR-SSP primer mix (Lot No. DR43-3) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023	30
 Table 4: Amplification patterns for all HLA-A* specificities detected by the HLA-A* CTS-PCR-SSP primer mixes (Lot No. A17-2) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023	34
 Table 5: Amplification patterns for all HLA-B* alleles detected by the HLA-B* CTS-PCR-SSP reagents (Lot No. B16-2) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023	45
 Table 6: Amplification patterns of HLA-DRB1* and other HLA-DRB* alleles detected by the HLA-DRB1* CTS-PCR-SSP primer mixes (Lot No. DR43-3) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023	68

4. Kit Composition

- Number of PCR primer mixes per test: 96:
 - 24 HLA-specific mixes for HLA-A typing
 - 48 HLA-specific mixes for HLA-B typing
 - 23 HLA-specific mixes for HLA-DRB1 typing
 - 1 negative control mix
- Number of tests per tray: 1
- Number of trays per kit: 10
- The primer mixes are aliquoted and dried in thin-walled, blue PCR-Trays.
- PCR buffer: 3.0 ml of Mastermix SSP (without Taq polymerase).

For storage condition, please refer to Section 1 of the 'Working instruction for the CTS-PCR-SSP TRAY and MINITRAY KITS' (Manual No. 100A) supplied along with this product.

5. Materials, Reagents and Equipment not supplied

Please refer to Section 2 of the 'Working instruction for the CTS-PCR-SSP TRAY and MINITRAY KITS' (Manual No. 100A) supplied along with this product.

6. Sample Requirements, PCR and Gel Electrophoresis

Please refer to Section 3 to 6 of the 'Working instruction for the CTS-PCR-SSP TRAY and MINITRAY KITS' (Manual No. 100A) supplied along with this product.

7. Result Evaluation

- Check the approximate size of the PCR product against the Primer Mix Specificity Tables (Appendix / Table 1, 2 and 3) to confirm the correct product size.
- Use the Amplification Pattern Tables (Appendix / Table 4, 5 and 6) to make the allele assignments or use the SCORE Software for detailed result interpretation.

8. Interpretation Hints

- Weak or false positive reactions can occur if you use a different Taq polymerase. In addition, the quality and quantity of DNA is a crucial factor and can affect the mix reactivities. Under suboptimal test conditions, some mixes could give rise to false positive reactions (if there are any potentially false positive reactions when using this specific lot, they are indicated in Section 9).
- Alleles that are known to amplify weakly are listed with "w" (= weak) in the tables (Appendix).
- Please also refer to Sections 7 and 8 of the 'Working instruction for the CTS-PCR-SSP TRAY and MINITRAY KITS' (Manual No. 100A) supplied along with this product.

9. Special notes

- **HLA-A* locus:**
 - The following mixes will give **short** PCR products (<200 bp): Mixes **5, 17, 20, 23 and 24**.
 - These PCR products might be difficult to distinguish from the 90 bp amplification control (internal positive control)!
 - In general, these specific products will give a much stronger signal than that of the controls. Furthermore, they will not have migrated as far into the gel as the control bands. If you are not sure whether the strong signal is due to a specific or a generic PCR product, you might let the gel run for additional 15 minutes at a lower voltage. By this way the specific band will be separated from the control band, and you will be able to clearly see a double signal at this position: a very strong specific band, and a weaker, shorter amplification control band.
 - Since mix **5** corresponds to A*23 and mix **23** to A*02 (both alleles/allele groups are common), it is recommended to start the interpretation by checking these reactions very carefully first!
 - Some mixes contain specific primers which may give rise to PCR fragments of **two or more** different sizes simultaneously (see Table 1).

- **Mix 4:** 225 bp and 235 bp (these two bands tend to merge into a single band). An additional weak band of 740 bp may occur, but is usually invisible.
 - The following mixes could give false positive reactions: Mixes **2, 7, 8, 21**.
- **HLA-B* locus:**
 - The following mixes will give short allele-specific PCR products (<200 bp): mixes 5, 9, 20, 43, and 47. These PCR products might be difficult to be distinguished from the 90 bp amplification control (internal positive control) PCR product. In general, the allele-specific product will give a much stronger signal than that of the amplification control and will not migrate as far into the gel as the amplification control band. If you are not sure whether the (strong) signal represents an allele-specific product or just the amplification control, you might let the gel run for an additional 15 minutes at a lower voltage. By this way, the allele-specific band can be separated from the amplification control band and you will be able to clearly see a double signal at this position: a strong allele-specific band and - compared to that - a weaker shorter amplification control band. Since three of these mixes correspond to quite frequent alleles (B*18 and B*35 = mix 5; B*27 = mix 9; B*15(62) = mix 43), it is recommended to start your evaluation by checking these four reactions carefully first!
 - Some HLA-B* mixes will amplify two or even three fragments of different sizes (e. g. mixes 9, 15, 19, 25 and 35). Mixes with 2 fragments of almost the same size are not mentioned here.
 - **Mix 9**
 - The most common B*27 alleles (B* 27:02, *27:04 and *27:05) will generate two fragments: a short fragment (155 bp) + a long fragment (525 bp). In this case the long fragment can be weak.
 - Amplification of B*27:06 will result only in a short fragment (155 bp).
 - B*47:01:01:01-47:10 will only give rise to a long fragment (525bp).
 - **Special notes to HLA-B*07**
 - The common HLA-B*07 alleles react positively with mix 1 and mix 28. Some HLA-B*07 alleles may cause additional weakly positive reactions with mix 24, especially if high amounts of Taq polymerase or high DNA-concentrations are used.
 - **Special notes to HLA-B*13**
 - B*13:01 reacts positively with mix 3, whereas B*13:02 is amplified by mix 3 and mix 34.
 - The rare alleles B*13:03 and B*13:04 can only be detected by mix 34.
 - Mix 34 may show different PCR products for B*13, for example:
Amplification of B*13:04 results in two fragments (430 bp and 490 bp).
Amplification of B*13:02:01/03/08Q/09/14-16/18/19 shows only the smaller fragment (430 bp).
Amplification of B*13:10 shows only the larger fragment (490 bp).
 - The following mixes could give false positive reactions: Mixes **31** and **37**.
 - **HLA-DRB1* locus**
 - Mix **5** and mix **6** detect (among other alleles) the most common HLA-DRB1*03 alleles which correspond either to serological split HLA-DR17 or HLA-DR18, respectively.
 - The following mix could give a false positive reaction: Mix **7**.

10. Troubleshooting

Please refer to Section 8 of the 'Working instruction for the CTS-PCR-SSP TRAY and MINITRAY KITS' (Manual No. 100A) supplied along with this product.

11. Precaution

Please refer to Material Safety Data Sheet for the CTS-PCR-SSP TRAY and MINITRAY KITS (Manual No. 100B) supplied along with this product.

12. Contact

If you have any particular questions concerning this kit, which are not answered in this or the Main Manual, please do not hesitate to contact us at:

Phone: +49 6221 564013

Fax: +49 6221 564200

E-mail: dna.labor@med.uni-heidelberg.de

13. Appendix

Table 1: Allele specificities and sizes of the PCR products of each **HLA-A*** CTS-PCR-SSP primer mix (**Lot No. A17-2**) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023

Position	Mix	Allele	Serology	Size
H1	Mix 1	A*01:01:01:01:01:01:64/01:01:67-01:01:130/01:01:131?/01:01:132-01:04:01:02N/01:06-01:33/01:35-01:37:01:02/01:39-01:71/01:73-01:101/01:103-01:138/01:140-01:166/01:168-01:199/01:201-01:220/01:222-01:243/01:245-01:266:02/01:267w/01:268-01:283/01:285N-01:379N/01:381-01:389/01:391-01:398/01:399?/01:400N-01:417N/01:419-01:426, A*11:27/11:38-11:39/11:94/11:209/11:271/11:326, A*23:31/23:45/23:55/23:128, A*24:15/24:41/24:51/24:92/24:235/24:499/24:504/24:551/24:557, A*25:36/25:59/25:65, A*26:29/26:49/26:205, A*32:13, A*66:10, B*07:64, B*15:12:01-15:12:03/15:19/15:270/15:298/15:304N/15:649	A1, Null, -, A26(10), B76(15)	225 bp
G1	Mix 2	A*01:200/01:244/01:380:01:01-01:380:01:02/01:427, A*02:609w/02:1022, A*03:107/03:164/03:413, A*11:17, A*26:72, A*32:24, A*36:01:01:01-36:13	-, A36	580 bp
F1	Mix 3	A*02:01:01:01-02:01:07/02:01:09-02:01:15/02:01:18-02:01:19/02:01:21-02:01:25/02:01:27-02:01:51/02:01:53-02:01:63/02:01:65w/02:01:66-02:01:68/02:01:70-02:01:73/02:01:75w/02:01:76-02:01:80/02:01:83-02:01:94/02:01:96-02:03:06/02:03:08-02:05:02/02:05:04-02:05:08/02:05:09w/02:05:10-02:06:11/02:06:13-02:06:18/02:06:20-02:06:30/02:06:31?/02:07:01:01-02:07:14/02:07:16-02:08:01/02:09:01:01-02:22:02/02:24:01-02:35:01/02:35:03-02:47/02:49:01-02:49:02/02:51/02:53N-02:54/02:56:01/02:57-02:59/02:60:01w-02:60:02w/02:61-02:63/02:66-02:71/02:73-02:79:02/02:81-02:87/02:89:01-02:97:02/02:99:01-02:99:02/02:101:01-02:102/02:104/02:106-02:113:02N/02:115-02:116/02:118/02:120-02:121/02:123-02:128/02:130-02:134/02:136-02:151/02:153:01-02:168/02:170-02:177/02:179-02:184/02:186-02:187/02:189-02:192/02:193w/02:194/02:196-02:222N/02:224-02:228/02:230-02:245/02:247-02:253/02:254w/02:255-02:274/02:276-02:278/02:280-02:285/02:286w/02:287/02:289:01/02:290-02:297/02:299-02:300/02:302-02:303/02:305N-02:307/02:308w/02:310/02:311w/02:312-02:338/02:340-02:347/02:349-02:376/02:378-02:382/02:384-02:405/02:407-02:411/02:413-02:453/02:455-02:487/02:488w/02:489-02:490N/02:492-02:503/02:505-02:509/02:511-02:526/02:528:01-02:570:02/02:572-02:581Q/02:583-02:590/02:591:02-02:593/02:594w/02:595-02:600/02:602-02:621/02:623-02:640/02:642/02:645-02:648:02/02:649w/02:650-02:661/02:663-02:675N/02:677/02:679-02:680/02:682-02:710N/02:711w/02:712-02:780/02:782-02:788N/02:790-02:825/02:827Q/02:830:01-02:830:02/02:832N-02:871N/02:873-02:875/02:877-02:922/02:924-02:926/02:928-02:943/02:945N-02:984/02:986-02:991/02:993-02:1020/02:1021?/02:1022/02:1023?/02:1024-02:1027/02:1031N-02:1032/02:1033?/02:1034-02:1045Q/02:1047-02:1086/02:1087?/02:1088-02:1089, A*23:01:21, A*24:02:76/24:310:02, B*07:02:80, B*15:10:07	A2, Low A2, -, A203, A210, Null	455 bp
E1	Mix 4	A*01:13/01:176/01:194/01:427, A*02:34w-02:35:03w/02:56:01w-02:56:02w/02:62w/02:103w/02:237w/02:580w/02:741w/02:809w/02:907w/02:927w, A*03:01:01:01-03:193/03:195-03:199/03:201-03:259/03:261-03:458Q, A*11:130/11:199:01/11:199:02w-11:199:03w/11:222, A*30:55/30:89/30:186, A*32:04/32:17w/32:52, A*34:08, A*36:02, A*68:71, A*69:07, A*74:13/74:42?	-, A2, A3, Null	see below
		A*01:13/01:176/01:194, A*02:34w-02:35:03w/02:56:01w-02:56:02w/02:62w/02:103w/02:580w/02:741w/02:907w, A*03:01:01:01-03:01:07/03:01:09-03:01:11/03:01:13-03:01:22/03:01:24-03:01:68/03:01:70-03:01:82/03:01:84-03:01:93/03:01:95-03:07:02/03:10-03:11N/03:13-03:31/03:33-03:35/03:37-03:40/03:42-03:56:01:02/03:58/03:60-03:71/03:73-03:87/03:90-03:106/03:110/03:112-03:127/03:129N-03:141/03:143-03:148/03:150-03:151/03:153-03:171/03:174-03:175/03:177/03:181-03:193/03:195/03:197:01N-03:197:02N/03:201-03:202/03:204/03:206-03:210/03:212-03:215/03:217-03:218/03:220/03:222-03:242/03:244/03:246-03:251/03:253-03:258/03:261-03:266N/03:268-03:272/03:274-03:299/03:301-03:318/03:320-03:329N/03:331-03:381N/03:383-03:395/03:397/03:399-03:412/03:414-03:428/03:430-03:458Q, A*11:199:01/11:199:02w-11:199:03w/11:222, A*30:55/30:186, A*32:17w, A*34:08, A*68:71, A*69:07, A*74:13/74:42?	-, A2, A3, Null	225 bp

Position	Mix	Allele	Serology	Size
		A*01:427, A*02:237w/02:809w/02:927w, A*03:01:01:01-03:01:11/03:01:13-03:01:29/03:01:31-03:01:34/03:01:36-03:01:49/03:01:51-03:01:89/03:01:91-03:01:105/03:01:107-03:04:03/03:07:01:01-03:09/03:11N/03:13-03:17:02/03:19-03:39/03:41/03:43-03:58/03:60-03:71/03:73-03:87/03:89:02-03:96/03:99-03:104/03:106-03:121/03:123:01-03:134/03:136-03:148/03:150-03:166/03:168N-03:178N/03:181-03:186/03:188-03:193/03:195/03:197:01N-03:199/03:201-03:204/03:206-03:207/03:209-03:220/03:222-03:224/03:226-03:242/03:244/03:246-03:259/03:261-03:288/03:290-03:299/03:301-03:318/03:320-03:325/03:327-03:331/03:333-03:395/03:397-03:451/03:453-03:458Q, A*32:04/32:52, A*36:02	-, A3, Null	740 bp
		A*01:427, A*03:01:01:01-03:01:05/03:01:07-03:01:29/03:01:31-03:01:34/03:01:36-03:01:48/03:01:51-03:01:89/03:01:91-03:01:105/03:01:107-03:04:03/03:07:01:01-03:09/03:11N-03:17:02/03:19-03:39/03:41/03:43-03:74/03:76-03:83/03:85-03:94/03:96/03:99-03:104/03:107-03:121/03:123:01-03:134/03:136-03:166/03:168N-03:176/03:178N-03:186/03:188-03:193/03:195-03:199/03:201-03:203/03:205-03:207/03:209-03:214/03:216-03:224/03:226-03:234Q/03:236-03:248/03:250-03:259/03:261-03:287/03:290-03:321/03:323N-03:325/03:327-03:331/03:333-03:342N/03:344:01-03:346/03:348-03:376/03:378-03:413/03:415-03:417/03:419-03:451/03:453-03:458Q, A*11:130, A*30:89, A*32:04, A*36:02	-, A3, Null	235 bp
D1	Mix 5	A*23:01:01:01-23:17:02/23:17:03?/23:18-23:23/23:25-23:33/23:35-23:56/23:58-23:65/23:67-23:68/23:70-23:103N/23:105-23:112/23:114-23:127/23:129, A*24:24/24:71/24:315/24:392/24:527, A*29:07/29:49, A*31:29, A*32:170, A*33:187, B*18:27	A23(9), -, Null	195 bp
C1	Mix 6	A*02:46/02:48/02:70/02:129/02:571/02:1050, A*03:30/03:152/03:273, A*11:325, A*23:01:01:01-23:02/23:04-23:08N/23:10:01-23:22/23:24-23:27/23:29-23:36/23:38N-23:50/23:52-23:60/23:62-23:68/23:70-23:79/23:81-23:82/23:84N-23:99/23:101-23:112/23:114-23:127/23:129, A*24:02:01:01-24:02:16/24:02:18-24:02:48/24:02:50-24:02:63/24:02:64w/24:02:65/24:02:66w/24:02:67-24:02:88/24:02:90-24:07:04/24:09N-24:11N/24:13:01-24:15/24:17:01:01-24:20:02/24:22-24:23/24:25-24:28/24:30/24:32-24:41/24:43-24:64/24:66/24:68-24:76/24:78-24:88/24:90:01N-24:128/24:130-24:155N/24:157/24:159-24:190/24:192-24:199/24:201-24:210/24:212-24:236/24:238-24:272/24:274-24:289/24:291-24:308/24:311-24:325/24:327-24:341/24:343w/24:344-24:358/24:360-24:369/24:371-24:388N/24:390-24:391/24:393-24:414/24:415w/24:416-24:427N/24:428Nw/24:429N-24:470/24:472-24:488/24:490-24:505/24:507-24:521/24:522?/24:523-24:561:02/24:563-24:589, A*25:77/25:82, A*31:08, A*32:79, A*33:21/33:53, A*68:249	A2, -, A23(9), Null, A24(9), Low A24(9), A2403, A9, A24(9)/A3	see below
		A*03:30/03:152/03:273, A*11:325, A*23:31/23:45/23:53/23:70, A*24:17:01:01-24:17:01:02/24:41/24:62/24:106/24:208:01-24:208:02:02/24:296/24:330/24:488/24:499/24:503/24:535/24:551/24:555/24:557, A*25:77/25:82, A*31:08, A*32:79, A*33:21/33:53, A*68:249	-, A24(9)	435 bp
		A*02:46/02:48/02:70/02:129/02:571/02:1050, A*23:01:01:01-23:02/23:04-23:08N/23:10:01-23:22/23:24-23:27/23:29-23:36/23:38N-23:50/23:52-23:60/23:62-23:68/23:71-23:79/23:81-23:82/23:84N-23:99/23:101-23:112/23:114-23:127/23:129, A*24:02:01:01-24:02:16/24:02:18-24:02:48/24:02:50-24:02:63/24:02:64w/24:02:65/24:02:66w/24:02:67-24:02:88/24:02:90-24:07:04/24:09N-24:11N/24:13:01-24:15/24:17:01:01-24:20:02/24:22-24:23/24:25-24:28/24:30/24:32-24:41/24:43-24:64/24:66/24:68-24:76/24:78-24:88/24:90:01N-24:128/24:130-24:155N/24:157/24:159-24:190/24:192-24:199/24:201-24:207:02/24:209-24:210/24:212-24:236/24:238-24:272/24:274-24:289/24:291-24:308/24:311-24:325/24:327-24:341/24:343w/24:344-24:358/24:360-24:369/24:371-24:388N/24:390-24:391/24:393-24:414/24:415w/24:416-24:427N/24:428Nw/24:429N-24:470/24:472-24:488/24:490-24:498/24:500-24:505/24:507-24:521/24:522?/24:523-24:550/24:552-24:554/24:556N-24:561:02/24:563-24:589, A*68:249	A2, -, A23(9), Null, A24(9), Low A24(9), A2403, A9, A24(9)/A3	470 bp

Position	Mix	Allele	Serology	Size
B1	Mix 7	A*11:178/11:190, A*23:105, A*24:02:01:01-24:02:73/24:02:75-24:02:145/24:02:147-24:03:04/24:05:01-24:05:02/24:07:01:01-24:11N/24:14:01:01-24:15/24:17:01:01-24:17:01:02/24:20:01:01-24:21:03/24:23/24:25-24:27/24:29-24:43/24:45N-24:64/24:66-24:86N/24:88/24:90:01N-24:93/24:95-24:106/24:108/24:110-24:128/24:130-24:137/24:139-24:166/24:168-24:187/24:189-24:203/24:205-24:206/24:208:01-24:210/24:212/24:214-24:226:02/24:229-24:284/24:286-24:288/24:291-24:298/24:300-24:314/24:316-24:338/24:340-24:354/24:356-24:372/24:374-24:376/24:377w/24:378Q/24:380/24:382-24:391/24:393-24:404/24:407-24:423/24:425N-24:445N/24:447Q-24:470/24:472-24:477/24:479Q-24:514N/24:516-24:555/24:556N?/24:557-24:559/24:560?/24:561:01-24:568N/24:570-24:571/24:573-24:589, A*25:04/25:44, A*32:02/32:22, B*51:55, B*57:129	-, A24(9), Low A24(9), A2403, A9, Null, A32(19)	510 bp
A1	Mix 8	A*02:06:31/02:52/02:135/02:309/02:454/02:992, A*25:01:01:01-25:01:01:12/25:01:02w/25:01:03-25:01:08/25:01:10-25:01:16/25:01:18-25:10/25:11w/25:12N-25:23/25:25-25:32/25:34-25:42N/25:44-25:58/25:60-25:74/25:75w/25:76:01-25:80/25:83-25:84, A*30:01:01:01-30:04:01:03/30:04:03-30:04:04/30:06-30:12:01/30:13-30:20/30:22-30:25/30:27N-30:81/30:83/30:85-30:89/30:91-30:95/30:97-30:123N/30:125-30:151/30:153-30:183/30:185-30:188/30:189?/30:190-30:205/30:206?, A*32:101Q, B*07:02:74w/07:225w/07:260w, B*57:148?, 04:01:49, 12:328, 15:222	-, A25(10), Null, A30(19)	see below
		A*02:52/02:992, A*30:01:01:01-30:02:09/30:02:11-30:04:01:03/30:04:03-30:04:04/30:06-30:12:01/30:13-30:20/30:22-30:25/30:27N-30:81/30:83/30:85-30:89/30:91-30:95/30:97-30:123N/30:125-30:151/30:153-30:183/30:185-30:188/30:189?/30:190-30:205/30:206?, B*57:148?, C*12:328, 15:222	-, A30(19), Null	550 bp
		A*02:06:31/02:135/02:309/02:454, A*30:02:10, C*04:01:49	-	560 bp
		A*32:101Q, B*07:225w/07:260w	-	395 bp
		A*25:01:01:01-25:01:01:12/25:01:02w/25:01:03-25:01:08/25:01:10-25:01:16/25:01:18-25:10/25:11w/25:12N-25:23/25:25-25:32/25:34-25:42N/25:44-25:58/25:60-25:74/25:75w/25:76:01-25:80/25:83-25:84, B*07:02:74w	A25(10), -, Null	405 bp
H2	Mix 9	A*02:06:31/02:135/02:309/02:454, A*25:01:01:01-25:01:07/25:01:08w/25:01:10-25:01:16/25:01:18-25:13/25:15-25:38/25:39w/25:40-25:42N/25:44-25:58/25:60-25:76:02/25:78-25:80/25:83-25:84, A*26:01:01:01-26:01:26/26:01:28-26:01:37/26:01:38w/26:01:39-26:01:58/26:01:60-26:01:64/26:01:66-26:03:01:02/26:03:03-26:18/26:20:01-26:30/26:31w/26:32-26:43:02/26:45-26:59/26:61-26:71N/26:73-26:75/26:76w/26:77-26:83/26:86-26:124/26:126-26:130/26:132-26:142/26:143w/26:145N-26:169/26:170w/26:171-26:183/26:185-26:203/26:205-26:229, A*34:01:01:01-34:11/34:13-34:18/34:20-34:22/34:23w-34:24w/34:25-34:30N, A*66:01:01-66:22/66:24-66:29/66:31-66:40/66:41w/66:42-66:47, A*68:294, B*07:02:74, B*18:02:02, B*38:02:04	-, A25(10), Null, A26(10), A10, A34(10), A66(10)	see below
		A*02:06:31/02:135/02:309/02:454, A*25:46/25:83, A*26:07:01-26:07:02/26:92, A*66:42	-, A26(10)	450 bp
		A*25:01:01:01-25:01:07/25:01:08w/25:01:10-25:01:16/25:01:18-25:13/25:15-25:38/25:39w/25:40-25:42N/25:44-25:45/25:47-25:58/25:60-25:76:02/25:78-25:80/25:84, A*26:01:01:01-26:01:26/26:01:28-26:01:37/26:01:38w/26:01:39-26:01:58/26:01:60-26:01:64/26:01:66-26:03:01:02/26:03:03-26:06/26:08:01:01-26:18/26:20:01-26:30/26:31w/26:32-26:43:02/26:45-26:59/26:61-26:71N/26:73-26:75/26:76w/26:77-26:83/26:86-26:91/26:93-26:124/26:126-26:130/26:132-26:142/26:143w/26:145N-26:169/26:170w/26:171-26:183/26:185-26:203/26:205-26:229, A*34:01:01:01-34:11/34:13-34:18/34:20-34:22/34:23w-34:24w/34:25-34:30N, A*66:01:01:01-66:22/66:24-66:29/66:31-66:40/66:41w/66:43-66:47, A*68:294, B*07:02:74, B*18:02:02, B*38:02:04	A25(10), -, Null, A26(10), A10, A34(10), A66(10)	445 bp
G2	Mix 10	A*01:01:56, A*02:06:31/02:135/02:309/02:454, A*25:11/25:81?25:82?, A*26:01:01:01-26:01:01:53/26:01:03-26:01:12/26:01:14-26:01:58/26:01:60-26:01:62/26:01:64/26:01:66-26:03:01:02/26:03:03-26:32/26:34-26:43:02/26:45-26:59/26:61-26:84/26:86-26:109/26:111-26:124/26:126-26:130/26:132-26:142/26:143w/26:145N-26:203/26:205-26:224/26:225?/26:226/26:228-26:229, A*30:02:10, A*43:01-43:03	-, A26(10), Null, A10, A43	see below
		A*02:06:31/02:135/02:309/02:454, A*25:11, A*26:03:01:01-26:03:01:02/26:03:03/26:06/26:21/26:30/26:36/26:78/26:92/26:111/26:146/26:177/26:193/26:229	-, A26(10)	415 bp

Position	Mix	Allele	Serology	Size
		A*01:01:56, A*25:81?-25:82?, A*26:01:01:01-26:01:01:53/26:01:03-26:01:12/26:01:14-26:01:58/26:01:60-26:01:62/26:01:64/26:01:66-26:02:02/26:04-26:05/26:07:01-26:20:02/26:22-26:29/26:31-26:32/26:34-26:35/26:37-26:43:02/26:45-26:59/26:61-26:77/26:79-26:84/26:86-26:91/26:93-26:109/26:112:01-26:124/26:126-26:130/26:132-26:142/26:143w/26:145N/26:147-26:176/26:178-26:192/26:194-26:203/26:205-26:224/26:225?/26:226/26:228, A*30:02:10, A*43:01-43:03	-, A26(10), Null, A10, A43	405 bp
F2	Mix 11	A*02:741, A*26:48/26:69, A*33:183, A*34:01:01:01-34:01:07/34:05/34:11-34:12/34:14/34:16-34:18/34:23/34:27	-, A34(10)	580 bp
E2	Mix 12	A*03:103:01-03:103:02/03:282, A*11:116/11:140/11:372, A*25:02/25:75, A*26:13/26:19/26:33/26:226, A*29:32, A*31:89/31:115, A*32:17, A*33:125/33:131, A*34:01:01:01-34:01:07/34:05-34:06/34:11-34:12/34:14/34:16-34:19/34:23/34:27, A*66:01:01:01-66:02:01:03/66:04/66:06-66:14/66:16-66:27N/66:29-66:38/66:40-66:47, A*68:262/68:294, A*74:13/74:42	-, A25(10), A34(10), A66(10), Null	415 bp
D2	Mix 13	A*01:51, A*03:24, A*11:10, A*25:43, A*30:163, A*33:49/33:144, A*34:02:01:01-34:04/34:07-34:10N/34:13/34:15/34:20-34:22/34:25-34:26N/34:28?/34:29N-34:30N, A*68:103:01-68:103:02, B*35:08:14w/35:14:03w	-, A3, A11, A34(10), Null	385 bp
C2	Mix 14	A*02:12-02:13/02:19/02:27/02:37-02:38/02:44/02:49:01/02:54/02:142/02:226N/02:271/02:280/02:417/02:543/02:630/02:662/02:682-02:683/02:723/02:804/02:840/02:948, A*03:02:01:01w-03:02:06w/03:10/03:31w-03:32w/03:65w/03:69Nw/03:73w/03:76w/03:82w/03:90w/03:106w/03:113w/03:160w/03:167/03:198w/03:218w/03:223w/03:225/03:236w-03:237w/03:242w/03:244w/03:253w/03:274w-03:275Nw/03:281w/03:285w/03:309w/03:314w/03:318w-03:319w/03:333w/03:348w/03:356w/03:385w/03:453w/03:458Qw, A*11:01:01:01-11:01:34/11:01:36-11:01:56/11:01:58-11:01:103/11:01:104w/11:01:105-11:16/11:18-11:22/11:23w/11:26-11:27/11:29-11:30/11:32:01-11:34/11:36-11:39/11:41-11:44/11:46-11:49/11:51-11:52Q/11:54:01-11:77/11:79-11:93/11:95-11:97/11:98w/11:99N-11:110/11:113-11:129/11:130w/11:131-11:157/11:159-11:175/11:177/11:179-11:182Q/11:184-11:189/11:192-11:210N/11:212-11:222/11:224-11:225/11:227-11:228/11:229w/11:230-11:249/11:250w/11:251N-11:257/11:258w/11:259-11:268/11:270/11:272Q-11:289/11:291-11:309/11:311-11:325/11:327-11:345/11:346w/11:347N-11:371/11:373-11:379/11:381-11:388N/11:389?/11:390-11:391/11:392?/11:393-11:395/11:396?-11:397N?/11:398-11:431/11:432N?/11:433N-11:434/11:435?/11:436N-11:438, A*24:19/24:28/24:44/24:89/24:290/24:379/24:406/24:424, A*30:164, A*31:24w/31:97w/31:136w, A*33:19w/33:59w/33:102w/33:150w/33:177w/33:190w, A*66:19, A*68:09/68:26/68:28/68:129/68:134/68:171:02N/68:204, A*69:03	A2, -, Null, A3, A11, A9, A24(9), A28	520 bp
B2	Mix 15	A*03:01:39, A*29:01:01:01-29:01:06/29:01:08-29:02:35/29:02:37-29:59/29:61-29:110/29:111w/29:112N-29:172, A*32:30:01-32:30:02/32:32, A*68:01:28/68:01:36/68:130:02, A*74:10/74:41-74:42	-, A29(19), Null	see below
		A*29:01:01:01-29:01:06/29:01:08-29:02:35/29:02:37-29:12/29:14-29:18/29:20-29:38/29:40-29:47/29:49-29:59/29:61-29:97/29:99-29:110/29:111w/29:112N-29:158/29:160-29:171, A*74:10/74:41-74:42	A29(19), Null, -	160 bp
		A*03:01:39, A*29:01:01:01-29:01:06/29:01:08/29:01:10-29:02:35/29:02:37-29:31/29:33-29:59/29:61-29:110/29:111w/29:112N-29:172, A*32:30:01-32:30:02/32:32, A*68:01:28/68:01:36/68:130:02	-, A29(19), Null	515 bp
A2	Mix 16	A*30:01:01:01-30:01:07/30:01:08w/30:01:09-30:02:09/30:02:11-30:04:04/30:06-30:11:02/30:13-30:17/30:19-30:20/30:22-30:25/30:27N-30:40/30:42-30:54/30:56-30:81/30:83/30:85-30:151/30:153-30:165/30:167-30:183/30:185-30:188/30:189?/30:190-30:205/30:206?	A30(19), -, Null	570 bp
H3	Mix 17	A*29:14, A*30:152, A*31:01:02:01-31:01:06/31:01:08-31:01:46/31:01:47?/31:01:48-31:07/31:09-31:20/31:23-31:74/31:76-31:88/31:91-31:103/31:104w/31:105-31:108/31:110-31:201N/31:202N?/31:203/31:204?/31:205-31:209/31:211-31:219	-, A31(19), Null	185 bp
G3	Mix 18	A*02:829, A*03:152/03:219, A*23:03:01/23:83, A*24:18/24:21:03/24:204/24:208:01/24:213, A*29:13, A*31:07-31:08/31:10, A*32:01:01:01-32:01:06/32:01:08-32:01:11/32:01:13-32:01:47/32:01:49-32:27N/32:29-32:30:01/32:31-32:33:01/32:33:03/32:35-32:65/32:67-32:70/32:72-32:74/32:75w/32:76-32:93/32:94w/32:95-32:100/32:102-32:157/32:159-32:173	-, A24(9)/A3, A31(19), A32(19), Null	see below

Position	Mix	Allele	Serology	Size
		A*02:829, A*23:03:01/23:83, A*24:21:03/24:208:01, A*29:13, A*31:07-31:08/31:10, A*32:01:01:01-32:01:06/32:01:08-32:01:11/32:01:13-32:01:47/32:01:49-32:03:01:02/32:05-32:27N/32:29-32:30:01/32:31-32:33:01/32:33:03/32:35-32:65/32:67-32:70/32:72-32:74/32:75w/32:76-32:93/32:94w/32:95-32:100/32:102-32:157/32:159-32:173	-, A31(19), A32(19), Null	410 bp
		A*03:152/03:219, A*24:18/24:204/24:213, A*32:04/32:52	-, A24(9)/A3	520 bp
F3	Mix 19	A*02:24:02/02:65/02:152/02:507, A*29:48/29:105/29:164, A*32:15/32:93, A*33:01:01:01-33:01:04/33:01:06-33:01:18/33:03:01:01-33:03:25/33:03:26w/33:03:28-33:03:30/33:03:31w/33:03:32-33:03:55/33:03:56?/33:03:57/33:03:59-33:18:01/33:20/33:22-33:29/33:31/33:33-33:34/33:36-33:37/33:39-33:52/33:54-33:65/33:66w/33:67-33:80N/33:82-33:88/33:90-33:121/33:123N-33:129N/33:131-33:134/33:136-33:139/33:140Nw/33:141-33:151/33:153-33:157N/33:158w/33:160-33:174N/33:176N-33:187/33:188w/33:189-33:214/33:215?/33:216-33:234, A*74:04/74:21/74:43, B*07:02:40, B*08:01:07, B*15:02:07, B*55:02:10	-, A2/19, A33(19), A19, Null, B8	see below
		A*02:24:02/02:65/02:152/02:507, A*29:48, A*32:93, A*33:08-33:09, A*74:04/74:21/74:43	-, A2/19, A33(19), A19	475 bp
		A*29:105/29:164, A*32:15, A*33:01:01:01-33:01:04/33:01:06-33:01:18/33:03:01:01-33:03:25/33:03:26w/33:03:28-33:03:30/33:03:31w/33:03:32-33:03:55/33:03:56?/33:03:57/33:03:59-33:07/33:10-33:18:01/33:20/33:22-33:29/33:31/33:33-33:34/33:36-33:37/33:39-33:52/33:54-33:65/33:66w/33:67-33:80N/33:82-33:88/33:90-33:121/33:123N-33:129N/33:131-33:134/33:136-33:139/33:140Nw/33:141-33:151/33:153-33:157N/33:158w/33:160-33:174N/33:176N-33:187/33:188w/33:189-33:214/33:215?/33:216-33:234, B*07:02:40, B*08:01:07, B*15:02:07, B*55:02:10	-, A33(19), Null, B8	470 bp
E3	Mix 20	A*02:01:165, A*03:01:39, A*29:19/29:39/29:48/29:159, A*68:01:28/68:01:36/68:130:02, A*74:01:01:01-74:09/74:11-74:15/74:16:02-74:40/74:43-74:44	-, A74(19), Null	160 bp
D3	Mix 21	A*02:01:09/02:05:05/02:06:07/02:07:12/02:11:02/02:35:03/02:50/02:76:02-02:76:03/02:122/02:243:02-02:243:03/02:528:02/02:570:01/02:591:01, A*03:171/03:271/03:373/03:428/03:442, A*11:57/11:158/11:359, A*24:94/24:138/24:188/24:228/24:293/24:515, A*29:51/29:69/29:73, A*32:109, A*68:01:01:01-68:01:10/68:01:12-68:01:24/68:01:26-68:02:06/68:02:08-68:47/68:51-68:128/68:130:01-68:136/68:138-68:183/68:185-68:217/68:219-68:229/68:231-68:233/68:235/68:237-68:250/68:252-68:257/68:258?/68:259-68:285/68:287-68:292N/68:294-68:298, A*69:01:01:01-69:06/69:08-69:09, B*44:497	A2, -, A68(28), A28, Null, A69(28)	645 bp
C3	Mix 22	A*02:34-02:35:03/02:56:01-02:56:02/02:62/02:78/02:103/02:741/02:907, A*03:241/03:343, A*69:01:01:01-69:09, 07:1034?	A2, -, A69(28)	390 bp
B3	Mix 23	A*02:01:01:01-02:01:07/02:01:10-02:01:11/02:01:12w/02:01:13-02:01:15/02:01:17-02:01:19/02:01:21-02:01:32/02:01:33w/02:01:34-02:01:39/02:01:41-02:01:81/02:01:83-02:03:10/02:04:01-02:05:04/02:05:06-02:06:06/02:06:08-02:07:11/02:07:13-02:11:01:06/02:11:03-02:22:02/02:24:01-02:35:01/02:35:02w/02:36-02:42:01/02:43N-02:76:01/02:77/02:79:01-02:97:02/02:99:01-02:99:02/02:101:01-02:121/02:123-02:140/02:142-02:243:01/02:244-02:311/02:313-02:469/02:471-02:528:01/02:529-02:569/02:571-02:590/02:591:02-02:635/02:636w/02:637-02:643N/02:645-02:673/02:674w/02:675N-02:740/02:741w/02:742-02:795Q/02:797N-02:802/02:804-02:814/02:815w/02:816-02:828/02:830:01-02:958/02:959?/02:960-02:967/02:968?/02:969-02:981/02:982?/02:983-02:1009/02:1010?/02:1011Q-02:1089, A*03:89:01, A*24:02:06/24:340, A*30:07w/30:13w/30:16w/30:22w/30:44w/30:46w, B*07:02:80w, B*35:01:60w, B*51:01:34w, B*58:01:24, 06:247w, 07:02:27w/07:577w	A2, Low A2, -, A203, A210, Null, A2/19, A24(9)	see below
		A*02:556, B*51:01:34w	-	175 bp
		A*02:50/02:556	A2, -	205 bp
		A*02:50/02:556	A2, -	225 bp

Position	Mix	Allele	Serology	Size
		A*02:01:01:01-02:01:07/02:01:10-02:01:11/02:01:12w/02:01:13-02:01:15/02:01:18-02:01:19/02:01:21-02:01:32/02:01:33w/02:01:34-02:01:39/02:01:41-02:01:50/02:01:53-02:01:62/02:01:64-02:01:73/02:01:75-02:01:76/02:01:78-02:01:81/02:01:83-02:01:113/02:01:115-02:01:152/02:01:154-02:03:10/02:04:01-02:05:02/02:05:04/02:05:06-02:06:06/02:06:08-02:06:11/02:06:13/02:06:15-02:06:19/02:06:21-02:07:11/02:07:13-02:07:14/02:07:16/02:07:18-02:11:01:06/02:11:03/02:11:05-02:11:08/02:11:10-02:22:01:02/02:24:01-02:35:01/02:36-02:42:01/02:43N-02:44/02:47/02:49:01-02:49:02/02:51-02:54/02:57-02:60:02/02:62-02:76:01/02:77/02:79:01-02:91/02:93:01-02:97:02/02:99:01-02:99:02/02:101:01-02:102/02:104-02:112/02:114-02:121/02:123-02:128/02:130-02:140/02:142/02:145-02:168/02:170-02:175/02:177-02:184/02:186-02:194/02:196-02:200/02:202:01-02:204/02:207-02:216/02:218-02:243:01/02:244-02:248/02:250N-02:311/02:313-02:314:02N/02:316-02:320/02:322-02:357/02:359-02:368/02:370-02:371/02:373N-02:377/02:379-02:409/02:411-02:418/02:420-02:469/02:471-02:481/02:483-02:503/02:505-02:513/02:515-02:528:01/02:529-02:543/02:545-02:569/02:572-02:590/02:591:02-02:600/02:602-02:635/02:636w/02:637-02:643N/02:645-02:673/02:674w/02:675N-02:696:02N/02:698-02:740/02:741w/02:742-02:745/02:747-02:765/02:767-02:795Q/02:797N-02:802/02:804-02:814/02:815w/02:816-02:828/02:830:01-02:872/02:874-02:880N/02:882-02:958/02:959?/02:960-02:967/02:968?/02:969-02:981/02:982?/02:983-02:1009/02:1010?/02:1011Q-02:1021/02:1023-02:1026/02:1029-02:1030/02:1032-02:1089, A*03:89:01, A*30:22w, B*07:02:80w, B*35:01:60w, B*51:01:34w, B*58:01:24, 07:02:27w/07:577w	A2, Low A2, -, A203, A210, Null, A2/19	135 bp
		A*02:01:01:01-02:01:07/02:01:10-02:01:11/02:01:12w/02:01:13-02:01:15/02:01:17-02:01:19/02:01:21-02:01:30/02:01:32/02:01:33w/02:01:34-02:01:39/02:01:41-02:01:66/02:01:68-02:01:81/02:01:83-02:01:121/02:01:123-02:01:126/02:01:128-02:03:10/02:04:01-02:05:04/02:05:06-02:06:06/02:06:08-02:07:05/02:07:07-02:07:11/02:07:13-02:07:23/02:09:01:01-02:11:01:06/02:11:03-02:19/02:21-02:22:02/02:24:01-02:33/02:35:02w/02:36-02:42:01/02:43N-02:49:02/02:51-02:54/02:57-02:61/02:63-02:76:01/02:77/02:79:01-02:97:02/02:99:01-02:99:02/02:101:01-02:102/02:104-02:121/02:123-02:127/02:129-02:138/02:140/02:142-02:157:02/02:159-02:168/02:170-02:194/02:196-02:207/02:210-02:243:01/02:244-02:249/02:251-02:257/02:259-02:271/02:273-02:291/02:293Q-02:311/02:313-02:316/02:318-02:350N/02:352-02:385/02:387-02:414/02:416-02:469/02:472-02:492/02:494-02:528:01/02:529-02:569/02:571-02:579/02:581Q-02:584/02:586-02:590/02:591:02-02:608N/02:610:01-02:633/02:635/02:636w/02:637-02:643N/02:645-02:662/02:664-02:673/02:674w/02:675N-02:676/02:678-02:701/02:703-02:737/02:739-02:740/02:742-02:786/02:788N-02:795Q/02:797N-02:802/02:804/02:806N-02:814/02:816-02:828/02:830:01-02:853/02:855-02:906/02:908-02:944/02:946N-02:957/02:959?/02:960-02:967/02:968?/02:969-02:981/02:982?/02:983-02:1000/02:1002-02:1009/02:1010?/02:1011Q-02:1021/02:1023-02:1027/02:1029-02:1089, A*03:89:01, A*24:02:06/24:340, A*30:07w/30:22w	A2, Low A2, -, A203, Null, A2/19, A24(9)	165 bp
		A*02:01:01:01-02:01:02:02/02:01:04-02:01:07/02:01:10-02:01:11/02:01:12w/02:01:13-02:01:15/02:01:17-02:01:19/02:01:21-02:01:32/02:01:33w/02:01:34-02:01:39/02:01:41-02:01:48/02:01:50-02:01:81/02:01:83-02:01:107/02:01:109-02:01:154/02:01:156-02:03:10/02:04:01-02:05:04/02:05:06-02:05:12/02:06:01:01-02:06:01:16/02:06:03-02:06:06/02:06:08-02:07:11/02:07:13-02:10/02:12-02:22:02/02:24:01-02:34/02:36-02:42:01/02:43N-02:47/02:49:01-02:49:02/02:51-02:68/02:70-02:76:01/02:77/02:79:01-02:80/02:82N-02:86:02/02:88N-02:89:02/02:91-02:97:02/02:99:01-02:99:02/02:101:01-02:111/02:113:01N-02:115/02:117-02:121/02:123/02:125N-02:127/02:130-02:135/02:137-02:140/02:142-02:145/02:147-02:183/02:185-02:243:01/02:244-02:246/02:248-02:260/02:262-02:296:02/02:299-02:307/02:309-02:311/02:313-02:330/02:332-02:365/02:367-02:380/02:382-02:393/02:395N-02:442/02:444-02:456/02:458-02:469/02:471-02:489/02:491-02:493/02:495-02:496/02:498-02:502/02:504-02:513/02:516N-02:528:01/02:529-02:533/02:535-02:569/02:572-02:576/02:578/02:581Q-02:583/02:585-02:588/02:590/02:591:02-02:626/02:628-02:635/02:636w/02:637-02:643N/02:646-02:673/02:674w/02:675N-02:700/02:702-02:704/02:706-02:713/02:715N-02:731/02:733-02:737/02:739-02:740/02:742-02:746/02:748N/02:750/02:752-02:753/02:755-02:756/02:758-02:772/02:774-02:795Q/02:797N-02:802/02:804-02:807N/02:809-02:814/02:816-02:820/02:822-02:828/02:830:01-02:892/02:894-02:905/02:908-02:922/02:924-02:925/02:927-02:931/02:933/02:935-02:939/02:942-02:958/02:960-02:963/02:965-02:967/02:968?/02:969-02:981/02:982?/02:983-02:985/02:987-02:990/02:992-02:995/02:997Q-02:1009/02:1010?/02:1011Q-02:1030/02:1033-02:1089, A*03:89:01, A*30:13w/30:16w/30:44w/30:46w, 06:247w	A2, Low A2, -, A203, A210, Null, A2/19	190 bp
A3	Mix 24	A*43:01-43:03, A*80:01:01-80:01:02/80:03/80:05-80:09N	A43, Null, -, A80	see below

Position	Mix	Allele	Serology	Size
		A*80:01:01:01-80:01:02/80:03/80:05-80:09N	A80, -, Null	170 bp
		A*43:01-43:03	A43, Null, -	445 bp

Amplification control (internal positive control): 90 base pairs (bp)

w = weak

? = nucleotide sequence information not available for the primer matching sequence or alleles with unknown reactivities

Bold: mixes which result in PCR fragments of different sizes (the specificities are first indicated all in one row, then split into several groups in the subsequent rows depending on the fragment size)

Table 2: Allele specificities and sizes of the PCR products of each **HLA-B*** CTS-PCR-SSP primer mix (**Lot No. B16-2**) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023

Position	Mix	Allele	Serology	Size
H4	Mix 1	B*07:02:01:01-07:02:05/07:02:07-07:02:32/07:02:33w/07:02:34-07:02:42/07:02:43w/07:02:44-07:02:84/07:02:86-07:05:05/ 07:05:06w/07:05:07-07:06:02/07:06:03w/07:06:04-07:18:02/07:20-07:32/07:34-07:39/07:41-07:52/07:55-07:59/07:61-07:93/07:95- 07:99/07:101-07:121/07:124-07:138/07:140-07:164/07:166-07:218/07:220-07:234/07:236-07:242/07:244-07:289/07:291-07:342/ 07:344-07:358/07:360Q-07:385/07:387/07:389-07:399/07:401-07:441/07:443-07:452/07:454-07:467Q/07:469N-07:470/07:471?/ 07:472N-07:474, B*08:79:01-08:79:02/08:156/08:187, B*15:297, B*18:35w/18:110w, B*35:66, B*38:63, B*40:25/40:106/40:166, B*42:05:01/42:05:02w, B*46:06/46:87, B*48:06, B*81:01:01:01-81:10, C*02:23, C*03:260, C*04:77, C*07:38:01-C*07:38:02	B7, -, B703, B7/ B42, Null, B81, Cw7	see below
		B*07:02:01:01-07:02:05/07:02:07-07:02:32/07:02:33w/07:02:34-07:02:42/07:02:43w/07:02:44-07:02:84/07:02:86-07:05:05/07:05:06w/ 07:05:07-07:06:02/07:06:03w/07:06:04-07:12/07:14:01-07:18:02/07:20-07:32/07:34-07:39/07:41-07:52/07:55-07:59/07:61-07:93/07:95- 07:99/07:101-07:121/07:124-07:138/07:140-07:164/07:166-07:218/07:220-07:234/07:236-07:242/07:244-07:289/07:291-07:342/07:344- 07:358/07:360Q-07:385/07:387/07:389-07:399/07:401-07:441/07:443-07:452/07:454-07:467Q/07:469N-07:470/07:471?/07:472N- 07:474, B*08:79:01-08:79:02/08:156/08:187, B*15:297, B*18:35w/18:110w, B*35:66, B*38:63, B*40:25/40:106/40:166, B*42:05:01/ 42:05:02w, B*48:06, B*81:01:01:01-81:10	B7, -, B703, B7/ B42, Null, B81	630 bp
		B*07:13, B*46:06/46:87, C*02:23, C*03:260, C*04:77, C*07:38:01-07:38:02	-, Cw7	615 bp
G4	Mix 2	B*07:29/07:186/07:383, B*08:01:01:01-08:01:31/08:01:33-08:05/08:07-08:12:02/08:13-08:25/08:27-08:48/08:51-08:61/08:63-08:93/ 08:95-08:106/08:108-08:111/08:113-08:169/08:171-08:180/08:182-08:225/08:227-08:231/08:233-08:238/08:240-08:241/08:243- 08:262/08:264-08:271/08:273-08:296/08:298/08:300Q-08:305/08:306?, B*13:93, B*15:142/15:180/15:429, B*38:84, B*39:12:01:01- 39:12:01:02/39:178, B*40:328, B*44:49/44:156, B*49:25, B*51:68/51:176, B*52:77, C*02:14:01-C*02:14:02/02:107/02:164, C*06:02:37/06:106:01/06:127:02/06:264, C*07:01:07/07:04:03/07:14/07:27:02/07:31:02/07:43:02/07:184/07:199:02/07:257:02/ 07:257:04/07:294/07:463/07:526:01/07:681/07:758/07:850, C*12:16:01-12:16:02/12:217	-, B8, Null, B39(16), Cw7	see below
		B*08:09/08:84, B*13:93, B*39:178, B*49:25, B*51:68/51:176, B*52:77, C*07:463/07:758	B8, -	545 bp
		B*07:29/07:186/07:383, B*08:01:01:01-08:01:31/08:01:33-08:05/08:07-08:08N/08:10-08:12:02/08:13-08:25/08:27-08:48/08:51-08:61/ 08:63-08:83/08:85-08:93/08:95-08:106/08:108-08:111/08:113-08:169/08:171-08:180/08:182-08:225/08:227-08:231/08:233-08:238/ 08:240-08:241/08:243-08:262/08:264-08:271/08:273-08:296/08:298/08:300Q-08:305/08:306?, B*15:142/15:180/15:429, B*38:84, B*39:12:01:01-39:12:01:02, B*40:328, B*44:49/44:156, C*02:14:01-02:14:02/02:107/02:164, C*06:02:37/06:106:01/06:127:02/06:264, C*07:01:07/07:04:03/07:14/07:27:02/07:31:02/07:43:02/07:184/07:199:02/07:257:02/07:257:04/07:294/07:526:01/07:681/07:850, C*12:16:01-12:16:02/12:217	-, B8, Null, B39(16), Cw7	560 bp
F4	Mix 3	B*13:01:01:01-13:02:32/13:02:33?/13:02:34-13:02:38/13:06-13:21/13:23/13:25-13:36/13:37w/13:38-13:39/13:41-13:44/13:47/13:49N- 13:57/13:59-13:61/13:63N-13:85/13:87-13:96/13:98-13:127/13:129-13:142/13:143w/13:144/13:146-13:152/13:154-13:158/13:160- 13:173/13:175-13:176	B13, -, Null	585 bp
E4	Mix 4	B*15:64:02, B*18:01:01:01-18:01:08/18:01:10-18:01:11/18:01:13-18:01:23/18:01:25-18:07:01/18:08-18:12:01/18:13-18:15/18:17N- 18:25/18:27-18:28/18:30-18:36/18:37:02-18:40/18:42-18:47/18:49-18:54/18:56-18:60/18:62-18:81/18:83-18:96/18:98:01-18:98:02/ 18:99w/18:100-18:106Q/18:107w/18:108-18:121/18:123-18:164/18:166-18:167/18:169-18:177/18:179-18:188/18:190/18:192-18:226, B*35:01:57/35:20:02, B*44:03:12w, B*53:01:26, C*06:02:101?, C*12:301	-, B18, Null, B35	465 bp

Position	Mix	Allele	Serology	Size
D4	Mix 5	B*07:65/07:134/07:400, B*08:32/08:133/08:180/08:239, B*15:202/15:239/15:308/15:376, B*18:01:01:01-18:03:01:03/18:04:01-18:08/18:10-18:15/18:17N-18:36/18:38-18:47/18:50-18:53/18:55-18:65/18:68-18:78/18:80-18:106Q/18:109-18:124/18:126-18:131:01:02/18:133-18:135:01:02/18:137-18:141/18:143-18:177/18:179-18:193/18:196-18:197/18:199-18:226, B*35:01:01:01-35:03:02/35:03:04-35:08:01:13/35:08:03-35:08:04/35:08:06-35:18/35:20:01-35:24:02/35:28:01:01-35:43:01/35:43:03-35:45/35:48/35:50-35:55/35:57-35:62/35:64:01-35:72/35:74-35:153/35:155-35:197/35:199-35:202/35:204-35:206/35:208:01-35:209/35:211-35:216N/35:218-35:253/35:255-35:270/35:272-35:273/35:275-35:298/35:300-35:310/35:312-35:316/35:318-35:354/35:356-35:360/35:362-35:449/35:451-35:470/35:472-35:521/35:522?/35:523-35:524/35:525?/35:526-35:549/35:552-35:569/35:571-35:575, B*37:11/37:14/37:37/37:67, B*39:19:01-39:19:02/39:186, B*40:185/40:209/40:308/40:350/40:365, B*45:22, B*46:40:01-46:40:02, B*55:130, B*56:06/56:45, B*58:64, B*78:01:01:02-78:05/78:07-78:08, C*07:1011	-, B18, Null, B35, B35/22, B39(16), B78	135 bp
C4	Mix 6	B*08:01:65/08:107, B*35:01:01:01-35:01:16/35:01:18-35:01:55/35:01:56w/35:01:57-35:01:66/35:01:68-35:02:05/35:02:06w/35:02:07-35:03:06/35:03:08-35:04:02/35:05:01:01-35:09:04/35:11:01-35:11:03/35:12:04/35:14:01-35:15:02/35:17:01:01-35:19/35:21/35:23-35:25/35:27/35:29:01-35:30:01:02/35:32:01/35:33-35:38/35:40N-35:42:02/35:45-35:50/35:51:02-35:58/35:61:02-35:62/35:64:01-35:66/35:68:01-35:68:02/35:70-35:71/35:74-35:75/35:77-35:78/35:81/35:84/35:89-35:94/35:96-35:98/35:100-35:101:02/35:103-35:111/35:113-35:116/35:119-35:134N/35:136-35:148/35:150:01-35:163/35:165N-35:175/35:177-35:184/35:187-35:190/35:192-35:198/35:200/35:202w/35:203-35:206/35:209-35:211/35:213-35:214/35:216N-35:225/35:227-35:229/35:231-35:232/35:235-35:241/35:243-35:250/35:252-35:260/35:262-35:264/35:267-35:269/35:271-35:275/35:277-35:293/35:295/35:297-35:316/35:319-35:360/35:362-35:372/35:373w/35:374-35:381N/35:383-35:413/35:415-35:424/35:426-35:431/35:433-35:438/35:440-35:462/35:464-35:470/35:472-35:473/35:478-35:533/35:535-35:540/35:542N-35:546/35:548/35:550-35:559/35:560?/35:562-35:569/35:572-35:574N/35:575?, B*49:67, B*51:56:02/51:372, B*53:01:01:01-53:06/53:08:01-53:08:02/53:10-53:16/53:18-53:27/53:29-53:37/53:39/53:41-53:43:02/53:45-53:61/53:63-53:72, C*06:147	-, B35, Null, B53	400 bp

Position	Mix	Allele	Serology	Size
B4	Mix 7	B*07:12/07:14:01w-07:14:02w/07:18:01/07:18:02w/07:77/07:137/07:162-07:163/07:193/07:235/07:248/07:359w, B*08:49/08:129w/08:255w, B*13:36, B*14:10w/14:102, B*15:02:01:01-15:02:04/15:02:05w/15:02:06-15:02:12/15:13:01-15:13:03/15:21:01:01-15:21:01:02/15:44/15:80w/15:88-15:89:02/15:112w/15:121/15:139/15:144/15:170/15:194/15:213-15:214/15:223/15:265/15:283/15:291/15:297/15:301-15:302N/15:308/15:319/15:325/15:328/15:330/15:341/15:345/15:357-15:358/15:374/15:378/15:394/15:399/15:404/15:418/15:420/15:425/15:437-15:438/15:454N/15:460w/15:491/15:495/15:505/15:537/15:557/15:567/15:569/15:574/15:592/15:596N/15:615/15:623/15:627/15:630/15:646/15:648, B*18:18:01:01-18:18:01:02/18:105/18:149w-18:150w/18:204, B*35:01:01:01-35:01:03/35:01:04w/35:01:05-35:01:72/35:01:73w/35:02:01:01-35:02:16/35:02:17w/35:02:18-35:03:19/35:03:20w/35:03:21-35:04:03/35:06-35:08:02/35:08:03w/35:08:04-35:09:04/35:11:01-35:11:03/35:12:01:01w-35:12:01:04w/35:12:02/35:12:03w-35:12:04w/35:14:01-35:15:02/35:18-35:21/35:23w/35:24:01/35:24:02w/35:25-35:27/35:29:01-35:29:03/35:33-35:36/35:38/35:39w/35:40N-35:42:02/35:45-35:50/35:52-35:57/35:59:01w/35:59:02/35:61:01w-35:61:02w/35:62/35:64:01-35:65Q/35:70-35:71/35:74/35:76-35:78/35:81-35:85/35:88/35:90-35:96/35:98/35:100-35:101:02/35:103-35:112:02/35:115-35:116/35:119/35:120w/35:121-35:124/35:126-35:134N/35:136-35:148/35:149w/35:150:01-35:150:02/35:152-35:175/35:177-35:184/35:187-35:198/35:200/35:201w/35:203-35:206/35:209-35:220/35:222-35:225/35:227-35:231/35:233w/35:236-35:241/35:243-35:250/35:252-35:257/35:258w/35:259-35:264/35:266w/35:267-35:269/35:271-35:272/35:274-35:275/35:277-35:285/35:287/35:288w/35:289-35:291/35:293-35:316/35:319-35:360/35:362-35:366/35:369-35:372/35:374-35:388/35:390N-35:394/35:395w/35:396-35:415/35:417-35:424/35:425w/35:426/35:428Q-35:439:01:02/35:441-35:458/35:460Q-35:470/35:473-35:474/35:478-35:480/35:482-35:497/35:500-35:508N/35:510-35:518/35:520/35:522/35:524/35:526-35:533/35:536-35:540/35:542N-35:548/35:550-35:569/35:571w/35:573-35:575, B*37:08/37:67, B*38:20:01w/38:20:02, B*39:09:01:01-39:09:01:06/39:09:03-39:09:05/39:42w/39:88/39:112/39:127/39:155/39:158/39:170-39:171, B*40:68, B*42:13, B*44:06w/44:12w/44:232w/44:263/44:356w/44:365w, B*51:04:01/51:04:02w/51:42w/51:46w/51:56:01w/51:56:02-51:56:03/51:139/51:200/51:227w/51:281w/51:359w, B*53:01:01:01-53:13/53:15-53:16/53:18-53:27/53:29-53:37/53:39-53:43:02/53:45-53:46/53:48N-53:66/53:68-53:71/53:72w, B*56:09/56:11/56:12w/56:60, B*83:01w, C*15:27	B7, -, B75(15), B77(15), B70, B15, Null, B35, B35/22, B37, B39(16), B44(12), B51(5), B53, B56(22), B55(22)	see below

Position	Mix	Allele	Serology	Size
		B*07:12/07:14:01w-07:14:02w/07:18:01/07:18:02w/07:137/07:162/07:235/07:359w, B*08:129w/08:255w, B*13:36, B*14:10w/14:102, B*15:02:01:01-15:02:04/15:02:05w/15:02:06-15:02:12/15:13:01-15:13:03/15:21:01:01-15:21:01:02/15:44/15:80w/15:88-15:89:02/15:112w/15:121/15:139/15:144/15:170/15:194/15:213-15:214/15:223/15:265/15:283/15:291/15:297/15:301-15:302N/15:308/15:319/15:325/15:328/15:330/15:341/15:345/15:357-15:358/15:374/15:378/15:394/15:399/15:404/15:418/15:420/15:425/15:437-15:438/15:454N/15:460w/15:491/15:495/15:505/15:537/15:557/15:567/15:569/15:574/15:592/15:596N/15:615/15:623/15:627/15:630/15:646/15:648, B*18:105/18:149w-18:150w, B*35:01:01:01-35:01:03/35:01:04w/35:01:05-35:01:72/35:01:73w/35:02:01:01-35:02:16/35:02:17w/35:02:18-35:03:19/35:03:20w/35:03:21-35:04:03/35:06-35:08:02/35:08:03w/35:08:04-35:09:04/35:11:01-35:11:03/35:12:01:01w-35:12:01:04w/35:12:02/35:12:03w-35:12:04w/35:14:01-35:15:02/35:18-35:21/35:23w/35:24:01/35:24:02w/35:25-35:27/35:29:01-35:29:03/35:33-35:36/35:38/35:39w/35:40N-35:42:02/35:45-35:50/35:52-35:57/35:59:01w/35:59:02/35:61:01w-35:61:02w/35:62/35:64:01-35:65Q/35:70-35:71/35:74/35:76-35:78/35:81-35:85/35:88/35:90-35:96/35:98/35:100-35:101:02/35:103-35:112:02/35:115-35:116/35:119/35:120w/35:121-35:124/35:126-35:134N/35:136-35:148/35:149w/35:150:01-35:150:02/35:152-35:175/35:177-35:184/35:187-35:198/35:200/35:201w/35:203-35:206/35:209-35:220/35:222-35:225/35:227-35:231/35:233w/35:236-35:241/35:243-35:250/35:252-35:257/35:258w/35:259-35:264/35:266w/35:267-35:269/35:271-35:272/35:274-35:275/35:277-35:285/35:287/35:288w/35:289-35:291/35:293-35:316/35:319-35:360/35:362-35:366/35:369-35:372/35:374-35:388/35:390N-35:394/35:395w/35:396-35:415/35:417-35:424/35:425w/35:426/35:428Q-35:432/35:434-35:439:01:02/35:441-35:458/35:460Q-35:470/35:473-35:474/35:478-35:480/35:482-35:497/35:500-35:508N/35:510-35:518/35:520/35:522/35:524/35:526-35:533/35:536-35:540/35:542N-35:548/35:550-35:569/35:571w/35:573-35:575, B*38:20:01w/38:20:02, B*39:42w, B*40:68, B*42:13, B*44:06w/44:12w/44:232w/44:263/44:356w/44:365w, B*51:04:01/51:04:02w/51:42w/51:46w/51:56:01w/51:56:02-51:56:03/51:139/51:227w/51:281w/51:359w, B*53:01:01:01-53:06/53:08:01-53:13/53:15-53:16/53:18-53:27/53:29-53:37/53:39-53:43:02/53:45-53:46/53:48N-53:64/53:66/53:68-53:71/53:72w, B*56:09/56:11/56:12w/56:60, B*83:01w, C*15:27	B7, -, B75(15), B77(15), B70, B15, Null, B35, B35/22, B44(12), B51(5), B53, B56(22), B55(22)	380 bp
		B*07:77/07:163/07:193/07:248, B*08:49, B*18:18:01:01-18:18:01:02/18:204, B*35:48/35:433, B*37:08/37:67, B*39:09:01:01-39:09:01:06/39:09:03-39:09:05/39:88/39:112/39:127/39:155/39:158/39:170-39:171, B*51:200, B*53:07/53:65	-, B37, B39(16), B53	395 bp
A4	Mix 8	B*15:13:01-15:13:03/15:418, B*27:30w, B*38:20:01w/38:20:02, B*44:06w/44:25w/44:50:01w/44:50:03/44:95w, B*49:61w, B*51:04:01/51:04:02w/51:42w/51:46w/51:56:01w/51:56:02-51:56:03/51:139/51:227w/51:281w/51:310/51:359w, B*52:88w, B*53:01:01:01-53:02/53:04-53:06/53:08:01-53:08:02/53:10/53:15-53:30/53:32-53:35/53:37/53:39-53:43:02/53:45-53:46/53:48N-53:54/53:56-53:61/53:62w/53:63-53:64/53:66/53:68-53:71/53:72w, B*57:05/57:131, B*58:01:01:01-58:01:15/58:01:17-58:01:19/58:01:20w/58:01:21-58:01:44/58:04-58:05/58:09-58:13/58:17N/58:19/58:21-58:24/58:28:01-58:29/58:31N-58:37/58:39N-58:42/58:44-58:45:02/58:47-58:59:02/58:61-58:63/58:65-58:74/58:75w/58:76-58:92/58:93Nw/58:94N-58:103/58:105-58:113/58:114w/58:115-58:120/58:123-58:131/58:133N-58:136/58:138-58:140, C*15:71	B77(15), -, B27, B44(12), B51(5), B53, Null, B58(17)	330 bp
H5	Mix 9	B*27:01-27:05:22/27:05:24-27:21:02/27:23-27:74/27:76-27:139/27:141-27:152/27:154-27:188/27:190-27:214/27:216-27:238/27:240-27:249/27:251-27:265, B*40:124:01-40:124:02/40:180/40:429, B*44:97, B*47:01:01:02-47:12, B*50:20	B27, -, B2708, Null, B47	see below
		B*27:02:01:01-27:05:10/27:05:12-27:05:22/27:05:24-27:11/27:13:01-27:15/27:17/27:19:01:01-27:21:02/27:24-27:25/27:27-27:28/27:30/27:32-27:38/27:40-27:43/27:45-27:58/27:60-27:74/27:76/27:78-27:91/27:93-27:100/27:102-27:118/27:120-27:128/27:130/27:132-27:136/27:138-27:139/27:141-27:152/27:154-27:156/27:158-27:174/27:176N-27:180/27:182-27:188/27:190-27:203/27:205-27:214/27:216-27:238/27:240-27:241/27:243N-27:249/27:251-27:256/27:258-27:265, B*44:97	B27, -, B2708, Null	155 bp

Position	Mix	Allele	Serology	Size
		B*27:01-27:05:22/27:05:24-27:05:31/27:05:33-27:05:35/27:05:37-27:05:59/27:08-27:10/27:12:01:01-27:19:01:02/27:23/27:25-27:26/27:28-27:32/27:36-27:40/27:42/27:44-27:45/27:47-27:74/27:77/27:79-27:90:04/27:92-27:101/27:103-27:105/27:108-27:124/27:126/27:128-27:129/27:131-27:134/27:137/27:139/27:141-27:150/27:152/27:155-27:163/27:165/27:167/27:169-27:179/27:181-27:185Q/27:188/27:190-27:191/27:193/27:195-27:207/27:209-27:214/27:216-27:217/27:220-27:221/27:223N-27:226/27:228-27:232/27:234-27:238/27:240-27:249/27:252-27:261N/27:263N-27:265, B*40:124:01-40:124:02/40:180/40:429, B*47:01:01:02-47:12, B*50:20	B27, -, B2708, Null, B47	530 bp
G5	Mix 10	B*07:27/07:236, B*08:126/08:131/08:139, B*18:54/18:195, B*37:01:01:01-37:01:11/37:01:13-37:04:02/37:06:01-37:09/37:12-37:13/37:15-37:21/37:23-37:28/37:30N-37:33N/37:35-37:36/37:38-37:66/37:68-37:94/37:96-37:104/37:105?/37:106-37:107, B*38:17, A*11:133, A*33:60, A*68:234	-, B37, Null	200 bp
F5	Mix 11	B*13:46/13:51/13:79/13:140, B*14:93, B*35:471, B*37:10/37:34, B*38:01:01:01-38:02:03/38:02:04w/38:02:05-38:02:09/38:02:10w/38:02:11-38:16/38:18-38:51/38:53-38:72/38:74-38:87/38:89-38:104/38:160-38:173Q, B*39:104, B*44:117/44:251, B*49:50, B*51:13:01-51:13:02/51:92:01-51:92:02/51:172, B*52:14:01-52:14:02/52:75, B*53:04/53:07/53:31/53:33, B*57:76/57:118, B*58:20/58:34, B*59:03, A*24:96/24:256	-, B37, B38(16), B16, Null, B53	395 bp
E5	Mix 12	B*07:45w/07:420w, B*08:55, B*15:48w/15:69/15:108w/15:136w/15:186/15:235/15:249w/15:265/15:449/15:494, B*18:168, B*35:35/35:236w/35:259w/35:276/35:468, B*39:01:01:01-39:01:01:29/39:01:03:01-39:01:10/39:01:12-39:06:04/39:06:06-39:20/39:22-39:24:01/39:24:03-39:29/39:31:01:01-39:31:01:02/39:32w/39:35/39:37:01:01-39:42/39:44-39:46/39:47w/39:48-39:49/39:50w/39:51-39:73/39:74w/39:75-39:79/39:81/39:83:01-39:101/39:102w/39:103/39:105-39:106/39:108-39:111/39:113-39:127/39:129-39:135/39:137-39:155/39:157N-39:173/39:174w/39:175N-39:184/39:185?/39:186-39:194Q, B*40:184/40:424, B*46:19w/46:93w, B*50:49, B*55:123, B*67:01:01-67:07, 03:215w, C*07:137:01-07:137:02/07:868, C*14:63	-, B62(15), B71(70), B35, B3901, Low B39, B3902, B39(16), Null, B67	520 bp
D5	Mix 13	B*07:45w/07:420w, B*38:26, B*55:123, B*67:01:01-67:01:03/67:03-67:07	-, B67	545 bp
C5	Mix 14	B*07:133, B*40:01:01-40:01:55/40:01:56w/40:01:57-40:04:02/40:06:01:01-40:11:02/40:13-40:16:01:02/40:18-40:20:01:02/40:22N-40:25/40:27:01-40:27:03/40:29-40:38/40:40/40:42-40:45/40:47-40:50:01:02/40:52-40:62/40:64:01:01-40:70:02/40:72:01-40:91/40:94-40:136/40:138-40:162/40:164-40:170/40:171w/40:172-40:173/40:175-40:176/40:177w/40:178-40:183/40:186:01-40:198/40:200-40:207/40:210-40:264/40:266-40:293/40:294w/40:295-40:307/40:309-40:314/40:315w/40:316-40:323:02/40:325-40:342/40:344-40:349/40:351-40:364/40:366-40:377/40:378w/40:379:01:01-40:423/40:425-40:429/40:431-40:457/40:459-40:495/40:497-40:522/40:524-40:528, B*41:19, B*44:16/44:21/44:31/44:64:01-44:64:02/44:124/44:150/44:347, B*45:10, B*47:01:01:02-47:07/47:09-47:12, B*48:41, B*49:07, B*50:31	-, B60(40), B61(40), Null, B40, B47, B44(12)	640 bp
B5	Mix 15	B*07:27w/07:50w/07:183w/07:400/07:453, B*08:04:01w-08:04:02w/08:17/08:54:01-08:54:02/08:110w/08:285, B*13:21/13:71/13:108/13:166/13:174, B*15:07:01:01-15:07:03/15:33/15:45/15:55w/15:68/15:78:01:01-15:78:04/15:107/15:116/15:126/15:141/15:150/15:188/15:207/15:248/15:287/15:324/15:331/15:352/15:360/15:383/15:405/15:421/15:431/15:450/15:524/15:614:01:01-15:614:01:02, B*18:14w/18:21w-18:22w/18:40w/18:61w/18:189w/18:217w, B*35:72w/35:534/35:541, B*37:07/37:09/37:12/37:81, B*40:01:01-40:01:08/40:01:10-40:06:31/40:09-40:12:02/40:14:01-40:16:01:02/40:18-40:24/40:25w/40:27:01-40:27:03/40:29-40:40/40:42-40:60/40:62-40:67/40:69-40:74/40:76/40:78-40:86/40:88-40:94/40:96-40:105/40:107-40:120/40:122-40:157/40:159-40:161/40:163-40:165/40:167-40:183/40:185-40:195/40:197-40:219/40:221/40:223-40:231/40:233-40:236/40:238-40:299/40:301-40:312/40:313w/40:314-40:321/40:323:01-40:325/40:326w/40:327-40:362/40:364-40:367/40:370-40:421Q/40:423-40:428N/40:430-40:457/40:459-40:480/40:482-40:512/40:514N-40:516/40:518-40:519/40:521-40:528, B*41:01:01:01-41:45N/41:47-41:79, B*42:11, B*44:31/44:53:01-44:54/44:106/44:110/44:135/44:158/44:166/44:182/44:184/44:192:01-44:192:04/44:211/44:213/44:343/44:506, B*45:27, B*48:01:01:01-48:01:12/48:03:01:01-48:05/48:07-48:24/48:27-48:30/48:32-48:57, B*49:40, B*50:60, B*52:54, C*07:937, C*12:301w, C*16:85	-, B62(15), B15, B35, B37, B60(40), B61(40), B4005, Null, B41, B44(12), B48, B48/40	see below

Position	Mix	Allele	Serology	Size
		B*07:27w/07:50w/07:183w/07:400/07:453, B*08:04:01w-08:04:02w/08:17/08:54:01-08:54:02/08:110w/08:285, B*13:71/13:108/13:174, B*15:07:01:01-15:07:03/15:45/15:55w/15:68/15:126/15:207/15:324/15:331/15:360/15:383/15:405/15:421/15:431/15:450/15:524, B*18:14w/18:22w, B*35:72w/35:541, B*37:07/37:09, B*40:02:01:01-40:02:08/40:02:10-40:03:02/40:05:01:01-40:05:02/40:09/40:15-40:16:01:02/40:18-40:19/40:27:01-40:27:03/40:29/40:31-40:32/40:35:01-40:35:03/40:37/40:39-40:40/40:45/40:50:01:01-40:50:01:02/40:56-40:58/40:71/40:78/40:80/40:82/40:85/40:89-40:91/40:94/40:97-40:98/40:104-40:105/40:107/40:111/40:115/40:119/40:122/40:133Q/40:142N/40:144N-40:145/40:157/40:164/40:169/40:174/40:176/40:180-40:181/40:189/40:200-40:202/40:205-40:206/40:209/40:211/40:214/40:219/40:224-40:225:01/40:226/40:229/40:246/40:248/40:255/40:266-40:267/40:271/40:274/40:276/40:283-40:284/40:287/40:289-40:293/40:296-40:297/40:302-40:305/40:309/40:314/40:320/40:331/40:334:01-40:335/40:345N/40:356/40:359/40:371/40:375/40:378/40:384/40:390/40:399N/40:405/40:409/40:412-40:415/40:420/40:432/40:435/40:437-40:438N/40:440-40:442/40:444/40:455/40:457/40:465/40:467/40:470/40:472/40:489/40:493/40:496/40:500/40:504-40:505/40:511N/40:514N/40:518/40:521-40:522, B*41:02:01:01-41:02:11/41:04/41:10-41:11/41:13/41:15/41:18-41:19/41:23-41:24/41:27/41:31/41:36/41:38-41:43/41:45N/41:47/41:49/41:51-41:52/41:55/41:58/41:62/41:65/41:68-41:69/41:71-41:72/41:74, B*42:11, B*44:54/44:106/44:110/44:135/44:158/44:166/44:184/44:211/44:213, B*48:01:01:01-48:01:11/48:04:01:01-48:05/48:07-48:16/48:18-48:20/48:22/48:24/48:27-48:30/48:32-48:39/48:41-48:56	-, B62(15), B15, B35, B37, B61(40), B4005, B60(40), B41, Null, B44(12), B48	375 bp
		B*07:27w/07:50w/07:183w/07:400/07:453, B*08:04:01w/08:17/08:110w, B*13:21/13:71/13:166/13:174, B*15:33/15:78:01:01-15:78:04/15:107/15:116/15:141/15:150/15:188/15:248/15:287/15:352/15:421/15:524/15:614:01:01-15:614:01:02, B*18:21w/18:40w/18:61w/18:189w/18:217w, B*35:534, B*37:07/37:09/37:12/37:81, B*40:01:01-40:01:08/40:01:10-40:02:13/40:02:15-40:02:19/40:02:22-40:02:25/40:02:27-40:06:31/40:09-40:12:02/40:14:01-40:16:01:02/40:18-40:24/40:25w/40:27:01-40:27:03/40:29-40:34/40:36-40:40/40:42-40:60/40:62-40:67/40:69-40:70:02/40:72:01-40:74/40:76/40:78-40:81/40:83-40:86/40:88-40:94/40:96-40:105/40:107-40:120/40:122-40:157/40:159-40:161/40:163-40:165/40:167-40:183/40:185-40:195/40:197-40:210/40:212-40:219/40:221/40:223-40:231/40:233-40:236/40:238-40:254/40:256N-40:299/40:301-40:312/40:313w/40:314-40:321/40:323:01-40:325/40:326w/40:327-40:362/40:364-40:367/40:370-40:421Q/40:423-40:428N/40:430-40:457/40:459-40:480/40:482-40:512/40:514N-40:516/40:518-40:519/40:521-40:528, B*41:01:01:01-41:02:01:01/41:03:01-41:09/41:12/41:14/41:16-41:17/41:20-41:23/41:25-41:26/41:28-41:35/41:37/41:44/41:48/41:50/41:53-41:54/41:56Q-41:68/41:70-41:79, B*44:31/44:53:01-44:53:02/44:182/44:192:01-44:192:04/44:343/44:506, B*45:27, B*48:01:01:01-48:01:12/48:03:01:01-48:05/48:07-48:11/48:13/48:15-48:24/48:27-48:28/48:30/48:33-48:49/48:51-48:57, B*49:40, B*50:60, B*52:54, C*07:937, C*12:301w, C*16:85	-, B15, B62(15), B37, B60(40), B61(40), B4005, Null, B41, B44(12), B48, B48/40	475 bp
A5	Mix 16	B*27:24/27:36/27:125/27:173, B*40:01:01-40:01:18/40:01:19w/40:01:20-40:01:80/40:07/40:22N-40:23/40:25/40:30-40:31/40:33-40:34/40:36/40:38/40:42-40:43/40:45-40:48/40:51-40:52/40:54-40:55/40:59-40:63/40:65-40:67/40:69/40:73-40:74/40:76/40:79-40:81/40:84/40:87:01-40:88/40:92/40:100-40:102/40:106/40:108/40:112-40:114:02/40:116-40:118N/40:123-40:126/40:128/40:130:01-40:130:02/40:132/40:134-40:135/40:138-40:141/40:146-40:147/40:149-40:156/40:160:01-40:160:02/40:163/40:166/40:168/40:170-40:172/40:175/40:178-40:179/40:182-40:183/40:185-40:188/40:191w/40:192-40:193/40:195-40:196/40:197w/40:198-40:199/40:204/40:207-40:208/40:212-40:213/40:215-40:218/40:221-40:223/40:227-40:228/40:233-40:236/40:238-40:242/40:245:01-40:245:02/40:247/40:249-40:253/40:257-40:265N/40:272-40:273/40:277-40:282/40:285-40:286N/40:288/40:298:01-40:301/40:308/40:310/40:312/40:315-40:316/40:319/40:321/40:323:01-40:324/40:326-40:327/40:329/40:332-40:333/40:336/40:338N-40:339/40:344/40:346/40:352-40:353:02/40:357-40:358/40:360/40:367-40:368/40:370/40:372N/40:377/40:379:01:01-40:379:02/40:381:01-40:383/40:385-40:386/40:391-40:398/40:400-40:404/40:406-40:407/40:416-40:419/40:421Q-40:422/40:424/40:426N/40:428N/40:430-40:431/40:433-40:434/40:436/40:439/40:443/40:446-40:448/40:450/40:452/40:456/40:461/40:463/40:466/40:469/40:474-40:475:02/40:479-40:480/40:483N/40:485-40:486/40:490-40:491/40:495/40:497-40:499/40:501-40:503/40:506N-40:510N/40:512-40:513/40:516-40:517/40:523-40:525/40:528, B*48:22/48:30	-, B60(40), Null, B40	690 bp

Position	Mix	Allele	Serology	Size
H6	Mix 17	B*07:133, B*08:155, B*15:46/15:53/15:106/15:143/15:212/15:525, B*18:48/18:107, B*35:19/35:47/35:63/35:154/35:217/35:274/35:317/35:450, B*37:06:01-37:06:02, B*39:69, B*40:02:01:01-40:06:31/40:08:01:01-40:09/40:11:01:01-40:11:02/40:13-40:16:01:02/40:18-40:20:01:02/40:24/40:26-40:29/40:32/40:35:01-40:35:03/40:37/40:39-40:40/40:44/40:50:01:01-40:50:01:02/40:53/40:56-40:58/40:64:01:01-40:64:02/40:68/40:70:01-40:72:03/40:75/40:77-40:78/40:82-40:83/40:86/40:89-40:91/40:94-40:99/40:103-40:105/40:107/40:109-40:111/40:115/40:119-40:122/40:127/40:129/40:131/40:133Q/40:136/40:142N-40:145/40:148/40:157-40:159/40:161-40:162/40:164-40:165/40:167/40:169/40:173-40:174/40:176/40:177w/40:180-40:181/40:184/40:189-40:190/40:200-40:203/40:205-40:206/40:211/40:214/40:219-40:220/40:224-40:226/40:229-40:230/40:232/40:237/40:243-40:244/40:246/40:248/40:254-40:256N/40:266-40:271/40:274-40:276/40:283-40:284/40:287/40:289-40:293/40:294w/40:295-40:297/40:302-40:307/40:309/40:311/40:313-40:314/40:317-40:318/40:320/40:322/40:330-40:331/40:334:01-40:335/40:337N/40:340-40:343/40:345N/40:347-40:349/40:351/40:354/40:356/40:359/40:361N-40:364/40:366/40:369/40:371w/40:373-40:376/40:378w/40:384/40:387-40:390/40:399N/40:405/40:408-40:415/40:420/40:423/40:425/40:427/40:429/40:432/40:435/40:437-40:438N/40:440-40:442/40:444-40:445/40:449/40:451/40:453-40:455/40:457-40:460/40:464-40:465/40:467-40:468/40:470-40:473/40:476/40:478/40:481N-40:482/40:484/40:487N-40:489/40:492-40:494/40:496/40:500/40:504-40:505/40:511N/40:514N-40:515/40:518-40:522/40:526-40:527, B*41:01:01:01-41:02:09/41:02:10?/41:02:11-41:79, B*44:02:01:01-44:02:12/44:02:14-44:03:36/44:03:38-44:03:50/44:03:52-44:05:06/44:07/44:09-44:30/44:32-44:39/44:41:01-44:43:02/44:45-44:51/44:53:01-44:56N/44:58N-44:59:02/44:61N-44:82/44:84:01-44:129/44:131-44:133/44:135-44:144/44:145w/44:146-44:155/44:157-44:168/44:170-44:181/44:183-44:200/44:202-44:209/44:211-44:232/44:234-44:245:02/44:247-44:365/44:437-44:444/44:446-44:524/44:526-44:529/44:530w/44:531-44:543/44:544N?/44:545-44:546, B*45:01:01:01-45:21/45:23-45:25/45:27-45:30, B*47:01:01:02-47:12, B*49:01:01:01-49:01:02/49:01:04-49:01:17/49:01:18w/49:01:19-49:02/49:04:01-49:14/49:16-49:25/49:27-49:72/49:74-49:78, B*50:01:01:01-50:02:01:02/50:04:01-50:05/50:07-50:22/50:31-50:74/50:76-50:79, B*51:112, B*53:30/53:39, B*56:67, B*57:45/57:51/57:69, B*58:111, C*16:85	-, B72(70), B35, B61(40), B4005, B60(40), B21, B41, Null, B44(12), B45(12), B12, B47, B49(21), B50(21)	580 bp
G6	Mix 18	B*15:339/15:342, B*18:165w, B*27:18/27:29, B*37:02/37:102, B*40:124:01-40:124:02/40:180/40:429, B*47:01:01:02-47:12, B*48:26, B*50:20	-, B27, B47	425 bp
F6	Mix 19	B*07:12/07:14:01w-07:14:02w/07:18:01/07:18:02w/07:137/07:162/07:174/07:202/07:222/07:235/07:248w/07:279/07:359w, B*08:129w/08:213:01-08:213:02/08:255w, B*14:10w/14:102, B*15:62/15:80w/15:124, B*27:19:01:01w-27:19:01:02w/27:24/27:30w/27:36/27:125/27:127w/27:173, B*35:25-35:26/35:49:01-35:49:02/35:198/35:254/35:361, B*38:20:01w/38:20:02/38:63, B*39:42w, B*40:12:01-40:12:02/40:46/40:430, B*42:13, B*44:40w/44:44w/44:130w/44:156w/44:210:01w-44:210:02w, B*48:01:01:01-48:04:02/48:06-48:07/48:09/48:11/48:15-48:40/48:42-48:57, B*56:09/56:11/56:12w/56:60, B*58:44, B*81:01:01:01-81:10, B*83:01w, C*03:342, C*15:99	B7, -, B70, B27, B35, B44(12), B48, B48/40, B56(22), B55(22), B81, Null	see below
		B*07:12/07:14:01w-07:14:02w/07:18:01/07:18:02w/07:137/07:162/07:235/07:248w/07:359w, B*08:129w/08:255w, B*14:10w/14:102, B*15:62/15:80w, B*27:19:01:01w-27:19:01:02w/27:30w/27:127w, B*35:25-35:26/35:49:01-35:49:02/35:198/35:254/35:361, B*38:20:01w/38:20:02, B*39:42w, B*42:13, B*44:40w/44:44w/44:130w/44:156w/44:210:01w-44:210:02w, B*48:02:01-48:02:03/48:17/48:25, B*56:09/56:11/56:12w/56:60, B*58:44, B*83:01w, C*03:342, C*15:99	B7, -, B70, B27, B35, B44(12), B48, B48/40, B56(22), B55(22)	435 bp
		B*07:174/07:202/07:222/07:279, B*08:213:01-08:213:02, B*15:124, B*27:24/27:36/27:125/27:173, B*38:63, B*40:12:01-40:12:02/40:46/40:430, B*48:01:01:01-48:01:12/48:03:01:01-48:04:02/48:06-48:07/48:09/48:11/48:15-48:24/48:26-48:40/48:42-48:57, B*81:01:01:01-81:10	-, B48, B48/40, B81, Null	580 bp

Position	Mix	Allele	Serology	Size
E6	Mix 20	B*13:36, B*15:02:01:01-15:02:12/15:08:01:01-15:08:02/15:11:01-15:11:07/15:13:01-15:13:03/15:15:01:01-15:15:01:02/15:21:01:01-15:21:01:02/15:31:01:01-15:31:01:03/15:44/15:55/15:76/15:88-15:89:02/15:112/15:121/15:139/15:144/15:148/15:170w/15:189:01-15:189:02/15:191/15:194/15:209N/15:213-15:215/15:223/15:265/15:283-15:284/15:291/15:297/15:301-15:303/15:305/15:319/15:325/15:328/15:330/15:341/15:345/15:350/15:357-15:358/15:367/15:374/15:378/15:384:01-15:384:02/15:394/15:399/15:404/15:418/15:420/15:425/15:437-15:438/15:454N/15:460/15:479/15:486/15:491/15:495/15:505/15:508/15:513/15:537/15:557/15:567/15:569/15:574/15:592/15:596N/15:615/15:622-15:623/15:627/15:629-15:630/15:637/15:646/15:648, B*18:143w, B*35:46/35:311, B*51:104:01-51:104:02	-, B75(15), B77(15), B62(15), B15, Null, B35	140 bp
D6	Mix 21	B*07:09:01-07:09:02/07:11/07:17/07:162/07:235/07:237/07:255/07:287, B*08:28/08:35/08:37/08:89/08:107/08:291, B*14:56, B*15:03:01:01-15:03:13/15:18:01:01-15:18:02/15:18:04-15:18:10/15:23/15:29:01:01-15:29:01:02/15:47:01-15:47:02/15:49/15:51-15:52/15:54/15:61-15:62/15:64:01-15:64:02/15:68-15:69/15:72/15:74/15:80/15:91/15:93/15:98/15:103/15:108/15:114-15:115/15:119/15:123-15:124/15:127/15:131-15:132/15:134/15:151/15:153/15:156/15:158/15:161/15:173/15:176/15:180/15:186/15:197:01-15:198/15:200/15:210/15:220:01:01-15:221/15:226N/15:235/15:238/15:242:01-15:242:02/15:255/15:259/15:263/15:266/15:274-15:275:02/15:281-15:282/15:290/15:292-15:294N/15:306-15:307/15:311/15:314/15:335/15:337/15:351/15:354/15:369/15:379-15:380N/15:382/15:388-15:389/15:397/15:414/15:426/15:433/15:444/15:447-15:448/15:451/15:475/15:487N/15:498/15:502/15:504:01-15:504:02/15:509/15:511/15:517/15:526-15:527/15:539/15:548/15:560/15:570-15:571/15:573/15:582/15:588/15:607/15:610/15:612/15:614:01:01-15:614:01:02/15:616Q-15:617/15:628/15:643/15:650, B*18:37:01-18:37:02/18:125, B*27:41/27:78/27:107, B*35:25-35:26/35:49:01-35:49:02/35:254, B*38:52/38:88, B*39:07/39:43/39:91/39:134/39:169, B*42:09/42:32, B*48:02:01-48:02:03/48:14/48:23/48:25, B*55:14/55:23/55:32/55:68, B*56:03/56:09/56:18/56:31-56:32/56:50, B*58:44, C*07:46/07:356/07:546	B7, -, B72(70), B71(70), B15, B70, B35, Null, B48, B22	500 bp
C6	Mix 22	B*08:62/08:226/08:242, B*27:83/27:128, B*41:01:01:01-41:07/41:10-41:19/41:20w/41:21-41:32/41:34-41:46/41:48-41:61/41:63-41:79, B*42:02:01:01-42:02:02/42:09/42:17-42:18/42:26-42:27, B*49:76	-, B41, Null, B42	790 bp
B6	Mix 23	B*07:04:01-07:04:02/07:19/07:25/07:146/07:349, B*08:123/08:132/08:146, B*42:01:01:01-42:02:02/42:04-42:10/42:12-42:18/42:20-42:25/42:27-42:33, B*44:90/44:356?, B*45:06, B*55:20/55:56, B*56:13, B*82:01:01:01-82:04, B*83:01	B7, -, B42, B56(22), B82	see below
		B*07:04:01-07:04:02/07:19/07:25/07:146/07:349, B*08:123/08:132/08:146, B*42:01:01:01-42:02:02/42:04-42:10/42:12-42:18/42:20-42:25/42:27-42:33, B*44:90/44:356?, B*45:06, B*55:20/55:56, B*56:13, B*82:01:01:01-82:04, B*83:01	B7, -, B42, B56(22), B82	550 bp
		B*07:04:01-07:04:02/07:19/07:25/07:146/07:349, B*08:123/08:132, B*42:01:01:01-42:02:02/42:04-42:06/42:08-42:10/42:12-42:18/42:20-42:25/42:27-42:33, B*44:90/44:356?, B*45:06, B*55:20/55:56, B*56:13, B*82:01:01:01-82:04, B*83:01	B7, -, B42, B56(22), B82	540 bp

Position	Mix	Allele	Serology	Size
A6	Mix 24	B*07:02:01:01w-07:02:05w/07:02:07w-07:04:02w/07:07:01:01w-07:12w/07:14:01w-07:31:02w/07:33:01w-07:33:03w/07:35w-07:39w/07:41w-07:52w/07:54w-07:64w/07:66w-07:68:03w/07:70w-07:77w/07:79w/07:81w-07:89w/07:91w-07:96:02w/07:98w-07:104w/07:106w-07:111Nw/07:113w-07:114w/07:116w-07:122w/07:124w-07:133w/07:135Nw-07:136:02w/07:139:01w-07:139:02w/07:141w-07:161Nw/07:163w-07:174w/07:177w-07:181Nw/07:183w-07:196w/07:198w-07:200w/07:202w-07:205w/07:208w/07:211w-07:212w/07:215w-07:221w/07:223w-07:224w/07:226w-07:236w/07:238w-07:245w/07:248w/07:250w-07:257w/07:259w/07:261w/07:263w/07:265w-07:268w/07:271w-07:277w/07:279w-07:282w/07:285Nw-07:286w/07:288w-07:292w/07:294w-07:303:02w/07:305w-07:315Nw/07:318Nw-07:323w/07:325Nw-07:326w/07:328w-07:331w/07:334w-07:339w/07:341w-07:348w/07:350w-07:351Nw/07:354w-07:355w/07:359w-07:367w/07:369w-07:375w/07:377w/07:379w-07:386Nw/07:388w/07:390w-07:394w/07:396w-07:398w/07:401w-07:406w/07:410w-07:425Nw/07:427w-07:433w/07:435Nw-07:442w/07:444w-07:456w/07:458w-07:459w/07:461w-07:473w, B*08:07w/08:28w/08:37w/08:87w/08:89w/08:107w/08:111w/08:156w/08:255w, B*14:23w/14:56w, B*15:10:01:01w-15:10:08w/15:37w/15:99w/15:312w-15:313w/15:323w/15:338w/15:445w/15:455w/15:503w/15:506w/15:563w, B*27:06:01:01w-27:06:02w/27:21:01w-27:21:02w/27:34w/27:91w/27:106w/27:127w/27:136w/27:138w/27:154w/27:192w/27:208w/27:222w/27:227w/27:262w, B*39:29w, B*40:09w/40:18w/40:24w/40:31w-40:33w/40:42w/40:44w/40:98w/40:120w/40:414/40:458/40:498, B*41:04w/41:21w/41:63w, B*42:16:01w-42:16:02w, B*44:02:01:01-44:02:12/44:02:14-44:02:69/44:02:70w/44:02:71-44:03:36/44:03:38-44:03:50/44:03:52-44:04/44:05:01:01w-44:05:06w/44:07/44:09/44:11-44:13/44:16/44:19N-44:30/44:32-44:42/44:44-44:56N/44:58N-44:59:02/44:61N/44:63-44:69:02/44:70w/44:71-44:74/44:75w/44:76/44:78w/44:79:01-44:81/44:84:01-44:106/44:107w/44:108N-44:116/44:118-44:122/44:123w/44:124-44:133/44:135/44:136w/44:137-44:138Q/44:139w/44:141-44:143/44:145w/44:146-44:155/44:156w/44:157/44:158w/44:159-44:161/44:163:01-44:165/44:167-44:183/44:185-44:197/44:198Nw/44:199-44:200/44:202-44:209/44:210:01w-44:210:02w/44:211-44:212/44:214-44:216/44:217Nw/44:218-44:225/44:227-44:229/44:231-44:235/44:236w/44:237N-44:245:02/44:247-44:249/44:250w/44:252-44:255/44:257-44:262/44:263w-44:264w/44:265-44:275/44:277-44:287/44:288w/44:289-44:307/44:309N-44:329/44:330w/44:331-44:338/44:341N-44:354N/44:356-44:365/44:437-44:444/44:446-44:458/44:460-44:464:01:02/44:465w/44:466N-44:524/44:526-44:529/44:530w/44:531-44:540/44:542-44:546, B*48:08w/48:46w, B*49:61, B*55:27w/55:110w, B*83:01w, C*03:342w	B7, -, B703, B7/ B42, Null, B8, B71(70), B27, B61(40), B60(40), B44(12), B45(12), B47, B21	490 bp
H7	Mix 25	B*07:20/07:140/07:153/07:359, B*15:14/15:91/15:131/15:161/15:518/15:551/15:578/15:592/15:634, B*18:56/18:169, B*35:45/35:71/35:169/35:309/35:452/35:545, B*39:173, B*40:132/40:178/40:458/40:523, B*44:02:10/44:06/44:09/44:34:01/44:46:01-44:46:02/44:75/44:90/44:129/44:131/44:254/44:277/44:356/44:438N, B*45:01:01:01-45:03/45:05-45:07/45:09/45:11-45:23/45:25/45:27-45:30, B*46:17, B*50:02:01:01-50:02:01:02/50:21/50:56, B*51:23/51:42/51:108/51:359, B*55:122, B*56:74, B*82:01:01:01-82:03, B*83:01, C*08:93, C*14:61	-, B76(15), B44(12), B45(12), Null, B82	see below
		B*44:02:10/44:06/44:34:01, B*51:23/51:42/51:108/51:359	B44(12), -	715 bp
		B*07:20/07:140/07:153/07:359, B*15:14/15:91/15:131/15:161/15:518/15:551/15:578/15:592/15:634, B*18:56/18:169, B*35:45/35:71/35:169/35:309/35:452/35:545, B*39:173, B*40:132/40:178/40:458/40:523, B*44:09/44:46:01-44:46:02/44:75/44:90/44:129/44:131/44:254/44:277/44:356/44:438N, B*45:01:01:01-45:03/45:05-45:07/45:09/45:11-45:23/45:25/45:27-45:30, B*46:17, B*50:02:01:01-50:02:01:02/50:21/50:56, B*55:122, B*56:74, B*82:01:01:01-82:03, B*83:01, C*08:93, C*14:61	-, B76(15), B45(12), Null, B82	545 bp
G7	Mix 26	B*07:133, B*15:46/15:53/15:106/15:143/15:212/15:525, B*40:05:01:01-40:05:02/40:15-40:16:01:02/40:23/40:26/40:28/40:32/40:51/40:95/40:98/40:148/40:158/40:161/40:174/40:183/40:198/40:309/40:355/40:466/40:496, B*44:76/44:79:01-44:79:02/44:146/44:150/44:459, B*47:10, B*49:01:01:01-49:01:01:18/49:01:03-49:02/49:04:01-49:10/49:12-49:14/49:16-49:25/49:27-49:44/49:46-49:72/49:74-49:78, B*50:01:01:01-50:02:01:02/50:04:01-50:05/50:07-50:22/50:31-50:48/50:50-50:69/50:71-50:74/50:76-50:79, B*51:112	-, B72(70), B4005, B61(40), B21, B49(21), Null, B50(21), B45(12)	605 bp

Position	Mix	Allele	Serology	Size
F7	Mix 27	B*07:78/07:182N/07:424/07:473, B*13:09/13:39/13:67, B*14:37, B*15:18:03/15:42/15:73:01:01-15:73:01:02/15:86/15:224/15:303/15:485-15:486/15:552, B*18:191, B*35:60/35:291/35:432/35:446/35:545, B*39:17, B*40:48/40:71/40:377/40:522, B*45:01:01:01-45:01:12/45:03-45:08/45:10-45:25/45:27-45:28N/45:30, B*46:11/46:18, B*50:01:01:01-50:01:17/50:01:19-50:02:01:02/50:04:01-50:08/50:10-50:13/50:15/50:17-50:19/50:21-50:22/50:31-50:32/50:34-50:40/50:42-50:79, B*54:01:01:01-54:03/54:05N/54:07-54:08N/54:10/54:13/54:16-54:35/54:36w/54:37-54:40/54:41w/54:43-54:45/54:46w, B*55:01:01:01-55:01:12/55:01:14-55:03/55:05/55:07/55:09-55:12/55:15-55:16/55:18-55:19/55:21-55:22/55:24-55:26/55:29-55:31/55:33-55:48/55:50/55:52/55:53w/55:54-55:57/55:59-55:67/55:69-55:80/55:82-55:85/55:87-55:89N/55:91-55:102/55:104-55:109/55:111-55:120/55:122-55:134, B*56:01:01:01-56:02:03/56:04:01:01-56:04:04/56:08/56:10/56:13-56:14/56:16-56:17/56:19N-56:20:02/56:23-56:30/56:33-56:45/56:47-56:49/56:51-56:64/56:66-56:87/56:89-56:91, B*82:01:01:01-82:04, C*01:90/01:136/01:145:01N, C*03:19/03:102/03:307/03:318N, C*07:289/07:536/07:725, 08:207, C*12:03:66/12:15/12:208, C*15:07:01:01-15:07:01:02/15:43/15:85/15:144	-, Null, B71(70), B62(15), B35, B60(40), B45(12), B50(21), B54(22), B55(22), B22, B56(22), B82, Cw3	395 bp
E7	Mix 28	B*07:02:01:01-07:02:05/07:02:07-07:02:27/07:02:28w/07:02:29-07:02:86/07:02:87?/07:02:88-07:02:97/07:04:01-07:07:03/07:09:01w-07:09:02w/07:10/07:11w/07:12/07:15/07:18:01-07:26/07:29-07:31:02/07:33:01-07:36/07:39-07:49N/07:51-07:68:03/07:70-07:84/07:86-07:95/07:97-07:114/07:116-07:117/07:118w/07:119-07:125/07:126w/07:127-07:131/07:132w/07:133-07:179/07:181N-07:196/07:198-07:206/07:208-07:218/07:220-07:224/07:226-07:227/07:229:01-07:234/07:235w/07:236/07:237w/07:238-07:245/07:246w/07:247-07:254/07:255w/07:256-07:258/07:261-07:283/07:285N/07:287-07:301/07:302w/07:303:01-07:305/07:307-07:313/07:315N-07:326/07:327w/07:328-07:336/07:337w/07:338-07:350/07:352-07:358/07:360Q-07:375/07:376w/07:377-07:399/07:401-07:427/07:429-07:452/07:454-07:474, B*08:123/08:132, B*15:76w/15:101w/15:255w, B*35:76w/35:229w/35:368w/35:561w, B*38:26, B*40:166/40:513w, B*42:01:01:01-42:02:02/42:04-42:06/42:08/42:09w/42:10/42:12-42:16:01/42:16:02w/42:17-42:19/42:21-42:25/42:27-42:31/42:32w/42:33, B*45:06, B*54:01:01:01-54:05N/54:06w-54:07w/54:08N/54:09w/54:10-54:13/54:14w/54:16-54:29/54:31-54:35/54:36w/54:37-54:40/54:42-54:45, B*55:01:01:01-55:02:10/55:02:12-55:05/55:07-55:17/55:19-55:22/55:23w/55:24-55:26/55:27w/55:28:01-55:31/55:32w/55:33-55:52/55:53w/55:54-55:57/55:59-55:60/55:62-55:63/55:64w/55:65-55:67/55:68w/55:69-55:77/55:79-55:93/55:95-55:117N/55:119-55:120/55:122-55:134, B*56:01:01:01-56:02:03/56:03w/56:04:01:01-56:07/56:09w/56:10/56:11w/56:12-56:13/56:15-56:16/56:18w/56:19N-56:21/56:22w/56:24-56:30/56:31w-56:32w/56:34-56:41/56:43/56:45-56:49/56:50w/56:51-56:53/56:55:01:02-56:66/56:68-56:87/56:88w/56:90, B*67:01:01-67:01:03/67:03-67:07, B*81:01:01:01-81:10, B*82:01:01:01-82:04	B7, -, B7/B42, Null, B35/22, B54(22), B55(22), B22, B56(22), B78, B67, B81, B82	see below
		B*07:02:01:01-07:02:05/07:02:07-07:02:27/07:02:28w/07:02:29-07:02:86/07:02:87?/07:02:88-07:02:97/07:04:01-07:07:03/07:09:01w-07:09:02w/07:10/07:11w/07:12/07:15/07:18:01-07:26/07:29-07:31:02/07:33:01-07:36/07:39-07:49N/07:51-07:68:03/07:70-07:77/07:79-07:84/07:86-07:95/07:97-07:114/07:116-07:117/07:118w/07:119-07:125/07:126w/07:127-07:131/07:132w/07:133-07:179/07:181N/07:183-07:196/07:198-07:206/07:208-07:218/07:220-07:224/07:226-07:227/07:229:01-07:234/07:235w/07:236/07:237w/07:238-07:245/07:246w/07:247-07:254/07:255w/07:256-07:258/07:261-07:283/07:285N/07:287-07:301/07:302w/07:303:01-07:305/07:307-07:313/07:315N-07:326/07:327w/07:328-07:336/07:337w/07:338-07:350/07:352-07:358/07:360Q-07:375/07:376w/07:377-07:399/07:401-07:423/07:425N-07:427/07:429-07:452/07:454-07:472N/07:474, B*08:123/08:132, B*15:76w/15:101w/15:255w, B*35:76w/35:229w/35:368w/35:561w, B*38:26, B*40:166/40:513w, B*42:01:01:01-42:02:02/42:04-42:06/42:08/42:09w/42:10/42:12-42:16:01/42:16:02w/42:17-42:19/42:21-42:25/42:27-42:31/42:32w/42:33, B*54:04/54:06w-54:07w/54:09w/54:11/54:14w/54:25/54:31/54:42, B*55:01:04/55:01:17/55:04/55:08/55:13-55:14/55:17/55:20/55:23w/55:27w/55:28:01-55:28:02/55:32w/55:34/55:39/55:42/55:44/55:47/55:49/55:51/55:53w/55:64w/55:68w/55:81/55:86/55:110, B*56:03w/56:05:01-56:06/56:09w/56:11w/56:12/56:15/56:18w/56:21/56:22w/56:31w-56:32w/56:50w/56:53/56:88w, B*67:01:01-67:01:03/67:03-67:07, B*81:01:01:01-81:10	B7, -, B7/B42, Null, B35/22, B54(22), B55(22), B56(22), B22, B78, B67, B81	410 bp

Position	Mix	Allele	Serology	Size
		B*07:78/07:182N/07:424/07:473, B*45:06, B*54:01:01:01-54:03/54:05N/54:07w/54:08N/54:10/54:12-54:13/54:16-54:29/54:31-54:35/54:36w/54:37-54:40/54:43-54:45, B*55:01:01:01-55:02:10/55:02:12-55:03/55:05/55:07/55:09-55:12/55:15-55:16/55:19/55:21-55:22/55:24-55:26/55:29-55:31/55:33-55:48/55:50/55:52/55:53w/55:54-55:57/55:59-55:60/55:62-55:63/55:64w/55:65-55:67/55:69-55:77/55:79-55:80/55:82-55:85/55:87-55:93/55:95-55:109/55:111-55:117N/55:119-55:120/55:122-55:134, B*56:01:01:01-56:02:03/56:04:01:01-56:04:04/56:07/56:10/56:13/56:16/56:19N-56:20:02/56:24-56:30/56:34-56:41/56:43/56:45-56:49/56:51-56:53/56:55:01:02-56:66/56:68-56:87/56:90, B*82:01:01:01-82:04	-, Null, B54(22), B55(22), B22, B56(22), B82	425 bp
D7	Mix 29	B*54:01:01:01-54:01:06/54:01:08-54:05N/54:07-54:08N/54:10-54:37/54:39-54:46, B*55:07/55:124, B*59:06, C*03:464, C*07:463w/07:758, C*14:92	B54(22), -, Null	435 bp
C7	Mix 30	B*07:20w/07:24/07:60/07:100/07:131/07:220/07:223/07:235/07:358/07:359w/07:464, B*15:76/15:101/15:255/15:509, B*35:76/35:229w/35:368/35:561, B*44:90/44:356w, B*45:06w, B*54:03/54:06/54:33/54:40, B*55:08/55:51/55:134, B*56:01:01:01-56:01:12/56:01:14-56:07/56:09/56:11/56:13/56:15-56:16/56:20:01-56:22/56:24-56:30/56:34-56:41/56:43/56:45-56:53/56:55:01:02-56:59/56:61-56:65/56:68-56:73/56:74w/56:75-56:88/56:90, B*82:01:01:01w-82:02:02w/82:03-82:04, B*83:01w, C*03:542	-, B7, B35/22, B56(22), B22, B78, Null, B82	560 bp
B7	Mix 31	B*44:18:01:01-44:18:01:02/44:355, B*49:01:01:01-49:01:24/49:03/49:06-49:14/49:16-49:17/49:19N-49:49/49:51-49:60N/49:62-49:69/49:71-49:78, B*51:15/51:62/51:106:01-51:106:02/51:157/51:223/51:225/51:281/51:300, B*52:25:01-52:25:02, B*54:12, B*55:103, B*56:65, B*58:73, B*59:01:01:01-59:01:01:02/59:04-59:07/59:09-59:10N/59:12, A*24:90:02N, C*15:71	-, B49(21), Null, B59	395 bp
A7	Mix 32	B*07:143, B*15:71/15:175/15:204/15:210/15:212/15:225/15:344/15:381/15:426/15:559/15:571, B*18:15/18:19/18:21/18:30/18:57:01-18:57:02/18:101/18:152/18:164, B*27:04:01-27:04:06/27:06:01:01-27:06:02/27:10/27:15/27:18/27:20-27:21:02/27:24-27:25/27:40/27:54/27:63/27:66N/27:68-27:69/27:79/27:86/27:91-27:92/27:100/27:103/27:105-27:109/27:112-27:115/27:120/27:138/27:147/27:149/27:153/27:165/27:168/27:173/27:183/27:192/27:208/27:210-27:211/27:214/27:222/27:226/27:229/27:240/27:245/27:247/27:257/27:262, B*37:07/37:97, B*40:05:01:01-40:05:02/40:15-40:16:01:02/40:23/40:26/40:28/40:32/40:51/40:95/40:98/40:148/40:158/40:161/40:174/40:183/40:198/40:309/40:355/40:466/40:496, B*46:14, B*49:01:01:01-49:01:01:18/49:01:03-49:09/49:12-49:14/49:16-49:23/49:26-49:30/49:32-49:44/49:46-49:78, B*50:01:01:01-50:02:01:02/50:04:01-50:12/50:14-50:22/50:31-50:42/50:44-50:48/50:50-50:69/50:71-50:79, B*51:103/51:116/51:151, B*52:51, A*26:89, C*03:278/03:504, C*06:77	-, B62(15), B27, Null, B4005, B61(40), B21, B49(21), B50(21), B45(12)	715 bp
H8	Mix 33	B*35:186/35:570, B*46:01:01:01-46:01:04/46:01:05w/46:01:07-46:01:25/46:01:27-46:05/46:07N-46:30/46:31w/46:32:01-46:40:02/46:42-46:44/46:46-46:53/46:55N-46:78/46:80-46:98N, B*56:14, B*67:02:01:01-67:02:01:02, C*01:32:01-01:32:02, C*02:56, C*03:102/03:263:01-03:263:03/03:514, C*04:180:01-04:180:02, C*05:217, C*06:20, C*07:81/07:168/07:896, C*08:123/08:139, C*14:82/14:92, C*15:126, C*16:98/16:102/16:110	-, B46, Null	470 bp

Position	Mix	Allele	Serology	Size
G8	Mix 34	B*07:409, B*13:02:01:01-13:04/13:08-13:10/13:14-13:16/13:18-13:19/13:26:01-13:27/13:30-13:35/13:37-13:38/13:40-13:42/13:44-13:49N/13:53-13:56:02N/13:59/13:62/13:64-13:70/13:72/13:74-13:75/13:81-13:82/13:84-13:85/13:88-13:89/13:91/13:93-13:94/13:96/13:98-13:105/13:111-13:121/13:123Q-13:128/13:130-13:134/13:136-13:137N/13:139N-13:140/13:143/13:145/13:147/13:149-13:160/13:164-13:167/13:169/13:172w/13:174-13:176, B*15:01:01:01-15:01:04/15:01:06-15:01:46/15:01:47w/15:01:48-15:01:49/15:01:50w/15:01:51-15:02:12/15:04:01:01-15:08:02/15:11:01-15:16:05/15:19-15:21:01:02/15:24:01-15:28/15:31:01:01-15:36/15:38:01-15:40:02/15:42-15:45/15:50/15:55-15:57/15:60/15:63/15:65-15:67/15:70-15:71/15:75-15:79N/15:81-15:85/15:87-15:89:02/15:92/15:94N-15:97/15:101-15:102/15:104-15:105/15:107/15:109-15:113/15:116-15:118/15:120-15:122/15:125:01-15:126/15:128-15:129/15:135-15:142/15:144-15:149N/15:152/15:154-15:155/15:157/15:159-15:160:01:02/15:163-15:167/15:169/15:170w/15:171-15:172/15:174-15:175/15:178-15:179:02/15:181N-15:185/15:187/15:189:01-15:195/15:199/15:201/15:202w/15:203-15:207/15:209N/15:211/15:213-15:215/15:217-15:218Q/15:222-15:223/15:225/15:227-15:228/15:231-15:234/15:236-15:237/15:240-15:241/15:244-15:251:03/15:254/15:256-15:258N/15:260-15:262:02N/15:264-15:265/15:267-15:272N/15:276-15:280/15:283-15:289/15:291/15:295-15:305/15:309-15:310/15:315-15:322/15:325-15:328/15:330-15:333/15:336/15:340-15:341/15:343-15:350/15:352-15:353/15:355/15:357-15:360/15:362-15:368/15:370-15:375N/15:377Q-15:378/15:381/15:383-15:387/15:390-15:395/15:398-15:402/15:404-15:410/15:412-15:413/15:415/15:417-15:422/15:425/15:428-15:432/15:437-15:438/15:440-15:443/15:449-15:450/15:452/15:454N/15:456-15:459/15:461/15:463N-15:474/15:476-15:485/15:489-15:497/15:499/15:501/15:505/15:507-15:508/15:510/15:512-15:515/15:518-15:522/15:524/15:528N-15:529/15:531/15:533:01-15:538/15:540N/15:542-15:545/15:547/15:549N-15:552/15:554/15:556-15:559/15:561-15:562N/15:564-15:569/15:572/15:574/15:576-15:581/15:583-15:587/15:589/15:590?/15:591-15:602/15:604N-15:605/15:608-15:609/15:613/15:615/15:618/15:621-15:627/15:629-15:636/15:638-15:641/15:645/15:647-15:649/15:651-15:652, B*18:143w/18:178, B*35:46/35:207, B*44:08/44:57/44:60, B*46:01:01:01-46:01:26/46:01:27w/46:01:28-46:03/46:05-46:39/46:41N-46:46/46:48-46:60/46:62-46:69/46:71-46:73/46:75-46:76/46:78-46:85/46:87-46:98N, B*51:104:01-51:104:02, B*52:02:01-52:02:02, B*55:35, B*56:39, B*57:01:01:01-57:01:01:27/57:01:03-57:01:32/57:01:34-57:01:46/57:04:01-57:04:03/57:06/57:08/57:10-57:11/57:13-57:15/57:18-57:27/57:29-57:31:02/57:33-57:38/57:40-57:41/57:43-57:44/57:47-57:50/57:52/57:53w/57:54-57:56/57:58-57:62/57:64-57:65/57:67:01-57:68/57:71-57:79N/57:81/57:83/57:85-57:93/57:95/57:97-57:100/57:102-57:106/57:108-57:117/57:119/57:121-57:123/57:125-57:136/57:138-57:147/57:149-57:154/57:156-57:165, B*58:36, C*01:73, C*07:336w	-, B13, B15/B21, Null, B62(15), B75(15), B76(15), B77(15), B63(15), B35, B44(12), B46, B57(17)	see below
		B*13:02:01:01-13:04/13:08-13:09/13:14-13:16/13:18-13:19/13:27/13:30-13:35/13:37-13:38/13:40-13:42/13:44-13:49N/13:53-13:56:02N/13:59/13:62/13:64-13:70/13:72/13:74-13:75/13:81-13:82/13:84-13:85/13:88-13:89/13:91/13:93-13:94/13:96/13:98-13:105/13:111-13:121/13:123Q-13:128/13:130-13:134/13:136-13:137N/13:139N-13:140/13:143/13:145/13:147/13:149-13:160/13:164-13:167/13:169/13:172w/13:175-13:176, B*15:04:01:01-15:04:01:03/15:04:03-15:04:04/15:42/15:83:01:01-15:83:01:02/15:137/15:155/15:303/15:310/15:485/15:499/15:531/15:552, B*46:11/46:18/46:32:01, B*51:104:01-51:104:02, B*52:02:01-52:02:02, B*55:35, B*56:39	B13, -, B15/B21, Null, B62(15)	430 bp

Position	Mix	Allele	Serology	Size
		B*07:409, B*13:04/13:10/13:26:01-13:26:02/13:72/13:127/13:174, B*15:01:01:01-15:01:04/15:01:06-15:01:46/15:01:47w/15:01:48-15:01:49/15:01:50w/15:01:51-15:02:12/15:04:01:01-15:08:02/15:11:01-15:16:05/15:19-15:21:01:02/15:24:01-15:28/15:31:01:01-15:36/15:38:01-15:40:02/15:43-15:45/15:50/15:55-15:57/15:60/15:63/15:65-15:67/15:70-15:71/15:75-15:79N/15:81-15:82/15:84-15:85/15:87-15:89:02/15:92/15:94N-15:97/15:101-15:102/15:104-15:105/15:107/15:109-15:113/15:116-15:118/15:120-15:122/15:125:01-15:126/15:128-15:129/15:135-15:136/15:138-15:142/15:144-15:149N/15:152/15:154-15:155/15:157/15:159-15:160:01:02/15:163-15:167/15:169/15:170w/15:171-15:172/15:174-15:175/15:178-15:179:02/15:181N-15:185/15:187/15:189:01-15:195/15:199/15:201/15:202w/15:203-15:207/15:209N/15:211/15:213-15:215/15:217-15:218Q/15:222-15:223/15:225/15:227-15:228/15:231-15:234/15:236-15:237/15:240-15:241/15:244-15:251:03/15:254/15:256-15:258N/15:260-15:262:02N/15:264-15:265/15:267-15:272N/15:276-15:280/15:283-15:289/15:291/15:295-15:302N/15:304N-15:305/15:309-15:310/15:315-15:322/15:325-15:328/15:330-15:333/15:336/15:340-15:341/15:343-15:350/15:352-15:353/15:355/15:357-15:360/15:362-15:368/15:370-15:375N/15:377Q-15:378/15:381/15:383-15:387/15:390-15:395/15:398-15:402/15:404-15:410/15:412-15:413/15:415/15:417-15:422/15:425/15:428-15:432/15:437-15:438/15:440-15:443/15:449-15:450/15:452/15:454N/15:456-15:459/15:461/15:463N-15:474/15:476-15:484/15:489-15:497/15:499/15:501/15:505/15:507-15:508/15:510/15:512-15:515/15:518-15:522/15:524/15:528N-15:529/15:531/15:533:01-15:538/15:540N/15:542-15:545/15:547/15:549N-15:551/15:554/15:556-15:559/15:561-15:562N/15:564-15:569/15:572/15:574/15:576-15:581/15:583-15:587/15:589/15:590?/15:591-15:602/15:604N-15:605/15:608-15:609/15:613/15:615/15:618/15:621-15:627/15:629-15:636/15:638-15:641/15:645/15:647-15:649/15:651-15:652, B*18:143w/18:178, B*35:46/35:207, B*44:08/44:57/44:60, B*46:01:01:01-46:01:26/46:01:27w/46:01:28-46:03/46:05-46:10/46:12-46:17/46:19-46:39/46:41N-46:46/46:48-46:60/46:62-46:69/46:71-46:73/46:75-46:76/46:78-46:85/46:87-46:98N, B*57:01:01:01-57:01:01:27/57:01:03-57:01:32/57:01:34-57:01:46/57:04:01-57:04:03/57:06/57:08/57:10-57:11/57:13-57:15/57:18-57:27/57:29-57:31:02/57:33-57:38/57:40-57:41/57:43-57:44/57:47-57:50/57:52/57:53w/57:54-57:56/57:58-57:62/57:64-57:65/57:67:01-57:68/57:71-57:79N/57:81/57:83/57:85-57:93/57:95/57:97-57:100/57:102-57:106/57:108-57:117/57:119/57:121-57:123/57:125-57:136/57:138-57:147/57:149-57:154/57:156-57:165, B*58:36, C*01:73, C*07:336w	- , B15/B21, B62(15), Null, B75(15), B76(15), B77(15), B63(15), B35, B44(12), B46, B57(17)	490 bp
F8	Mix 35	B*57:01:01:01-57:10/57:12-57:81/57:84-57:117/57:119-57:130N/57:132-57:165, B*58:14	B57(17), -, Null	see below
		B*57:01:01:01-57:04:03/57:06-57:10/57:12-57:81/57:84-57:117/57:119-57:130N/57:132-57:165, B*58:14	B57(17), -, Null	360 bp
		B*57:04:01/57:05	B57(17), -	410 bp
E8	Mix 36	B*57:01:19/57:14:01, B*58:01:01:01-58:01:04/58:01:06-58:01:28/58:01:30-58:02:02/58:04-58:26/58:28:01-58:29/58:31N-58:76/58:78-58:81/58:83-58:115/58:117-58:118/58:120-58:134/58:135?/58:136-58:137/58:138?/58:139-58:140	-, B58(17), B17, Null	490 bp

Position	Mix	Allele	Serology	Size
D8	Mix 37	B*15:13:01-15:13:03/15:16:01:01w-15:16:05w/15:67w/15:95w/15:222w/15:254w/15:362w/15:408w/15:418/15:550w/15:602w/15:613w, B*27:30w, B*38:20:01w/38:20:02, B*40:109/40:340, B*41:60w, B*44:06w/44:18:01:01w-44:18:01:02w/44:25w/44:50:01w/44:50:03/44:95w/44:355w, B*49:01:01:01w-49:01:01:03w/49:01:01:04/49:01:01:05w-49:01:24w/49:03w-49:14w/49:16w-49:46w/49:48w-49:53w/49:54/49:55w-49:69w/49:71?-49:72?/49:73w-49:78w, B*51:01:01:01-51:04:01/51:04:02w/51:05:01:01-51:05:01:02/51:07:01-51:24:05/51:26-51:40:02/51:42w/51:43-51:44N/51:46w/51:48-51:53/51:55/51:56:01w/51:56:02-51:58/51:60-51:61:02/51:63:01-51:63:02/51:65-51:77/51:79-51:80/51:82-51:113/51:114w/51:115-51:134/51:136-51:147/51:149N-51:189/51:191-51:193:02/51:195-51:212/51:214-51:219/51:221-51:226/51:227w/51:228-51:246/51:248-51:257/51:259-51:280/51:281w/51:282-51:292/51:294-51:318N/51:320-51:340/51:342-51:351/51:354-51:358/51:359w/51:360-51:370/51:372w/51:373-51:377N, B*52:01:01:01-52:19/52:21:01-52:67/52:69-52:87/52:88w/52:89N-52:111, B*53:01:01:01-53:02/53:04-53:06/53:08:01-53:08:02/53:10/53:15-53:30/53:32-53:35/53:37/53:39-53:43:02/53:45-53:46/53:48N-53:54/53:56-53:64/53:66/53:68-53:71/53:72w, B*54:12, B*55:103, B*56:21/56:65, B*57:05/57:131, B*58:01:01:01-58:01:15/58:01:17-58:01:19/58:01:20w/58:01:21-58:01:44/58:04-58:05/58:08:01-58:13/58:17N/58:19/58:21-58:24/58:28:01-58:29/58:31N-58:37/58:39N-58:42/58:44-58:45:02/58:47-58:59:02/58:61-58:63/58:65-58:74/58:75w/58:76-58:92/58:93Nw/58:94N-58:103/58:105-58:113/58:114w/58:115-58:120/58:122-58:131/58:133N-58:136/58:138-58:140, B*59:01:01:01-59:07/59:09-59:12, C*15:71	B77(15), -, B63(15), B27, B44(12), B49(21), Null, B51(5), B5102, B5103, B52(5), B53, B58(17), B17, B59	330 bp
C8	Mix 38	B*08:274, B*35:01:10/35:03:32/35:04:02/35:251:01, B*51:01:01:01-51:01:03/51:01:05-51:01:10/51:01:11w/51:01:12/51:01:14-51:01:16/51:01:17w/51:01:18-51:02:03/51:02:05-51:06:04/51:08:01:01-51:08:04/51:08:06-51:09:03/51:11N-51:12/51:13:02-51:14/51:16-51:21/51:23/51:24:03-51:24:05/51:26-51:41N/51:43/51:46/51:48-51:61:02/51:63:01/51:64-51:80/51:82-51:91/51:92:02-51:106:01/51:107w/51:108-51:117/51:119-51:135/51:137-51:147/51:149N-51:156/51:158:01-51:167/51:169/51:170w/51:171-51:193:01/51:193:02?/51:194-51:196/51:198/51:200-51:213/51:215-51:219/51:221/51:223-51:224/51:226-51:244/51:246-51:250/51:252-51:294/51:296-51:299/51:301-51:309/51:311-51:326/51:328-51:337/51:339/51:341-51:357N/51:360-51:377N, B*78:01:01:02-78:02:02/78:04/78:07-78:10	-, B35, B51(5), B5102, B5103, Null, B78	450 bp
B8	Mix 39	B*07:02:11/07:05:03, B*08:274, B*15:09:01-15:09:02/15:252/15:435, B*18:205, B*35:01:10/35:03:32/35:04:02/35:251:01, B*51:01:01:01-51:01:03/51:01:05-51:01:12/51:01:14-51:01:37/51:01:39-51:02:03/51:02:05-51:08:04/51:08:06-51:09:03/51:11N-51:12/51:13:02-51:14/51:16-51:23/51:24:03-51:24:05/51:26-51:41N/51:43-51:44N/51:46/51:48-51:61:02/51:63:01/51:64-51:80/51:82-51:91/51:92:02-51:106:01/51:107w/51:108-51:135/51:137-51:147/51:149N-51:156/51:158:01-51:167/51:169/51:170w/51:171-51:177/51:179-51:193:01/51:193:02?/51:194-51:196/51:198/51:200-51:219/51:221-51:224/51:226-51:230:01:03/51:232:01-51:294/51:296-51:299/51:301-51:309/51:311-51:326/51:328-51:337/51:339-51:357N/51:360-51:377N, B*55:01:04/55:02:11/55:28:02, B*56:05:01-56:06/56:21, B*78:01:01:02-78:04/78:07-78:10	B7, -, B70, B35, B51(5), B5102, B5103, Null, B55(22), B56(22), B78	460 bp
A8	Mix 40	B*13:02:38/13:62, B*15:01:02/15:277-15:278, B*18:205w, B*40:01:06/40:06:30/40:26/40:28/40:95, B*41:02:09, B*44:03:40/44:62/44:459, B*48:01:08, B*49:18:01, B*50:01:21/50:14, B*51:07:01w-51:07:02w/51:214w/51:222w/51:251w/51:310/51:340w, B*52:01:01:01-52:01:04/52:01:06-52:01:14/52:01:16-52:05/52:07:01-52:13/52:14:02-52:24/52:25:02-52:72/52:73w/52:74-52:96N/52:98-52:111, B*53:28, B*78:05-78:06	-, B62(15), B60(40), B21, B51(5), B52(5), Null	450 bp
H9	Mix 41	B*07:02:11/07:05:03, B*08:274, B*15:01:02/15:09:01-15:09:02/15:243/15:252/15:277-15:278/15:435, B*18:205, B*27:153, B*35:01:10/35:03:32/35:04:02/35:251:01, B*40:01:06/40:06:30/40:26/40:28/40:95, B*41:02:09, B*48:01:08, B*50:01:21/50:14, B*55:01:04/55:02:11/55:28:02, B*56:05:01-56:06, B*78:01:01:02-78:10, C*12:03:31, C*14:02:37/14:92	B7, -, B62(15), B70, B35, B60(40), B21, B55(22), B56(22), B78	410 bp

Position	Mix	Allele	Serology	Size
G9	Mix 42	B*14:19/14:27, B*15:09:01-15:10:08/15:18:01:01-15:18:10/15:21:01:01-15:21:01:02/15:23/15:37/15:51/15:72/15:80/15:90/15:99/15:108/15:115/15:119/15:133-15:134/15:153/15:161w/15:176/15:189:01-15:189:02/15:197:01-15:198/15:200/15:221/15:226N/15:229/15:263/15:275:01-15:275:02/15:290/15:292-15:294N/15:306-15:307/15:311-15:314/15:329/15:335/15:337-15:338/15:351/15:354/15:380N/15:382/15:388/15:394/15:404/15:414/15:420/15:426-15:427/15:433-15:436/15:439/15:445/15:447/15:451/15:455/15:487N/15:498/15:503-15:504:02/15:506/15:526/15:530/15:548/15:563/15:570/15:573/15:607/15:611-15:612/15:615/15:620/15:623/15:643/15:646, B*35:26/35:82/35:85/35:135/35:294, B*38:30/38:65, B*39:32-39:33/39:74, B*44:263w, B*49:26, B*51:22, B*56:33, B*78:03	-, B70, B71(70), B75(15), Null	575 bp
F9	Mix 43	B*07:409, B*13:09/13:67, B*15:01:01:01-15:01:04/15:01:06-15:01:46/15:01:48/15:01:51-15:01:70/15:01:71?/15:01:72-15:01:74/15:01:76-15:02:12/15:04:01:01-15:08:02/15:11:01-15:12:03/15:14-15:15:01:02/15:19-15:21:01:02/15:25:01:01-15:28/15:30:01:01-15:35:02/15:38:01-15:40:02/15:42/15:44-15:45/15:48/15:50/15:55-15:56/15:58/15:60/15:63/15:65-15:66/15:70-15:71/15:73:01:01-15:73:01:02/15:75-15:79N/15:81-15:86/15:88/15:92/15:94N/15:96-15:97/15:101-15:102/15:104-15:105/15:107/15:109-15:113/15:116-15:118/15:120-15:122/15:125:01-15:126/15:128-15:129/15:135-15:142/15:144-15:150/15:152/15:154-15:155/15:159-15:160:01:02/15:163-15:167/15:169/15:170w/15:171-15:172/15:174-15:175/15:178-15:179:02/15:181N-15:185/15:187-15:195/15:199/15:201/15:202w/15:203-15:207/15:209N/15:211/15:213-15:215/15:217-15:219/15:223-15:225/15:227-15:228/15:231-15:234/15:236-15:237/15:240-15:241/15:244-15:251:03/15:257-15:258N/15:260-15:262:02N/15:264-15:265/15:267/15:269-15:272N/15:276-15:280/15:283-15:288/15:291/15:295-15:305/15:309-15:310/15:315-15:317/15:319-15:322/15:324-15:328/15:330-15:334/15:336/15:340-15:350/15:352-15:353/15:355/15:357-15:359/15:363:01-15:368/15:370-15:375N/15:377Q-15:378/15:381/15:383-15:387/15:390-15:395/15:398-15:402/15:404-15:407/15:409-15:410/15:412-15:413/15:415-15:417/15:419-15:422/15:425/15:428-15:432/15:437-15:438/15:440-15:441/15:443/15:449-15:450/15:452-15:454N/15:456-15:461/15:463N-15:474/15:476-15:486/15:488-15:497/15:499/15:501/15:507-15:508/15:510/15:512-15:515/15:518-15:522/15:524/15:528N-15:529/15:531/15:533:01-15:538/15:540N-15:545/15:547/15:549N/15:551-15:552/15:554/15:556-15:559/15:561-15:562N/15:564-15:569/15:572/15:574/15:576-15:578/15:580-15:581/15:583-15:585/15:587/15:589/15:590w/15:591-15:593/15:595-15:601/15:604N-15:606/15:608-15:609/15:615/15:618/15:621-15:627/15:629-15:641/15:645-15:649/15:651-15:652, B*18:143w/18:178, B*35:46/35:207/35:311, B*40:21/40:137, B*46:72, B*55:35, B*56:39, B*57:12	-, B62(15), Null, B75(15), B76(15), B15, B35	150 bp
E9	Mix 44	B*15:12:01/15:12:02w/15:12:03/15:14/15:19/15:270/15:298/15:304N/15:518/15:551/15:578/15:592/15:634/15:649, B*44:08/44:57/44:60, B*46:17, B*57:07/57:93/57:119	B76(15), -, Null, B44(12)	650 bp
D9	Mix 45	B*07:219, B*15:16:01:01-15:17:07/15:67/15:95/15:168/15:177/15:196/15:208/15:216/15:230/15:254/15:268/15:273/15:356/15:361-15:362/15:403/15:408/15:411/15:423-15:424/15:446/15:462/15:500/15:516/15:523/15:532/15:546Q/15:550/15:555/15:575N/15:602-15:603/15:613/15:619/15:642/15:644, B*57:13/57:22/57:57/57:63/57:161, B*58:06/58:08:01-58:08:02/58:19/58:46/58:78/58:91	-, B63(15), Null, B17	530 bp
C9	Mix 46	B*14:01:01:01-14:01:11/14:01:13-14:02:02/14:02:04-14:02:05/14:02:07-14:04/14:07N/14:09/14:11-14:12/14:14-14:20/14:22-14:36/14:38-14:41N/14:42:02-14:52/14:54/14:57-14:95/14:97-14:101N/14:103-14:116	B64(14), -, B65(14), B14, Null	400 bp
B9	Mix 47	B*14:02:01:01-14:02:02/14:02:04-14:02:05/14:02:07-14:03/14:05-14:06:02/14:09/14:11/14:13/14:15-14:18/14:20/14:22-14:23/14:24w/14:25/14:27/14:29-14:31/14:33-14:39/14:41N/14:42:02-14:45/14:48/14:50-14:52/14:56-14:64/14:66-14:70Q/14:72N-14:76N/14:79N/14:84-14:85N/14:88-14:94/14:96/14:98-14:101N/14:106-14:107/14:109/14:111/14:114-14:116, B*15:189:01-15:189:02, B*35:01:64/35:03:17/35:26/35:29:03, B*38:05/38:33, B*39:04/39:129, B*49:26, B*51:01:68	B65(14), -, B14, Null, B38(16), B39(16)	185 bp
A9	Mix 48	B*51:241, B*73:01:01:01-73:03, C*07:463w/07:758	-, B73	300 bp

Amplification control (internal positive control): 90 base pairs (bp)

w = weak

? = nucleotide sequence information not available for the primer matching sequence or alleles with unknown reactivities

Bold: mixes which result in PCR fragments of different sizes (the specificities are first indicated all in one row, then split into several groups in the subsequent rows depending on the fragment size)

Table 3: Allele specificities and sizes of the PCR products of each **HLA-DRB1*** CTS-PCR-SSP primer mix (**Lot No. DR43-3**) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023

Position	Mix	Allele	Serology	Size
H10	Mix 1	DRB1*01:01:01:01-01:01:30/01:01:32-01:01:44/01:01:45?/01:01:46-01:02:11/01:02:13-01:17/01:19-01:21/01:23-01:32/01:34-01:38/01:39Nw-01:40Nw/01:41-01:43/01:45-01:48/01:49w/01:50-01:59/01:61-01:64/01:65:01w-01:65:02w/01:66-01:75/01:77-01:133/01:135-01:141	DR1, -, DR103, Null	140 bp
G10	Mix 2	DRB1*15:01:01:01-15:01:21/15:01:23-15:33/15:35-15:112/15:113Nw/15:114-15:216, DRB1*16:09:01-16:10:02/16:33/16:36-16:37/16:58/16:66	DR15(2), -, DR2, Null, DR16(2)	140 bp/ 215 bp
F10	Mix 3	DRB1*15:21, DRB1*16:01:01:01-16:01:08/16:01:09w/16:01:10-16:05:02/16:07-16:49/16:51-16:73	-, DR16(2), DR2, Null	210 bp
E10	Mix 4	DRB1*03:01:01:01-03:01:01:25/03:01:02:01w-03:01:02:02w/03:01:03/03:01:05-03:01:13/03:01:14w/03:01:15-03:02:01:02/03:02:02w/03:02:03-03:05:02/03:05:03w/03:06-03:11:01/03:11:03-03:13:01/03:13:02w/03:14-03:15:01:02/03:15:02w/03:16-03:41:02/03:43-03:75/03:77-03:83/03:85-03:86/03:88-03:91/03:93-03:96/03:98-03:110/03:112-03:124/03:127w/03:128-03:129/03:130w/03:131-03:134/03:135w/03:136-03:144/03:145w/03:146-03:156N/03:158-03:161/03:163-03:173/03:175-03:176/03:178-03:189N/03:190w/03:191-03:205, DRB1*11:07:01/11:07:02w/11:53/11:103:01-11:103:02/11:105/11:107/11:125/11:136/11:173, DRB1*14:104w/14:111w/14:171, DRB3*01:14w	DR17(3), -, DR18(3), DR3, Null, DR11(5)	225 bp
D10	Mix 5	DRB1*03:01:01:01-03:01:41/03:04:01-03:04:02/03:06/03:08/03:10-03:11:01/03:11:03/03:12?/03:13:01-03:13:02/03:15:01:01-03:16/03:18-03:20/03:22-03:23/03:25:01-03:26/03:28/03:30-03:34/03:36-03:37/03:39/03:42-03:48/03:50-03:52:02/03:54w/03:55-03:56/03:58-03:62/03:64-03:73/03:75-03:80/03:82-03:87/03:89/03:91-03:96/03:98-03:99/03:100:01w-03:100:02w/03:104/03:106-03:112/03:114/03:116-03:118/03:121/03:123-03:129/03:132-03:153/03:156N-03:166/03:168-03:170/03:172-03:175/03:177/03:180-03:184/03:186-03:187/03:189N-03:204N/03:205w, DRB1*09:06, DRB1*11:59/11:83/11:135/11:182, DRB1*13:27:01-13:27:02/13:71/13:129/13:144/13:176/13:326, DRB1*14:82/14:95/14:132, DRB3*01:34, DRB3*03:33	DR17(3), -, DR3, Null, DR13(6)	220 bp
C10	Mix 6	DRB1*03:02:01:01-03:02:06/03:05:01-03:05:03/03:09/03:14/03:27/03:29/03:38/03:53/03:74/03:81/03:88/03:90/03:102-03:103/03:115/03:122/03:130-03:131/03:154/03:167/03:176/03:178-03:179/03:188, DRB1*11:09:01:01-11:09:01:02/11:20/11:28:02/11:87/11:122, DRB1*13:02:01:01-13:02:03/13:02:05-13:02:22/13:05:01:01-13:05:03/13:26:01-13:26:02/13:29:01:01-13:29:01:02/13:31/13:34/13:36/13:39/13:41/13:50:02-13:50:03/13:56/13:63/13:65/13:67/13:73-13:74/13:85/13:96:01-13:96:02/13:99/13:103-13:104/13:109/13:119-13:120/13:123-13:124/13:126/13:128/13:136/13:139-13:140/13:142N-13:143/13:145/13:147/13:155/13:157-13:159/13:165/13:168/13:170-13:171:02/13:179/13:189/13:196/13:198-13:199/13:202-13:203/13:207-13:212/13:220/13:225/13:235-13:237/13:239-13:241/13:243-13:244/13:255N/13:257/13:264-13:265/13:268N-13:269/13:275/13:277/13:280/13:282/13:285/13:288-13:289N/13:292-13:293/13:296-13:298N/13:301-13:302/13:306/13:310N/13:312/13:318/13:320/13:325/13:331-13:332/13:334, DRB1*14:02:01:01-14:02:07/14:02:09-14:03:02/14:09/14:13/14:19/14:24/14:27:01-14:27:02/14:30/14:46-14:48/14:51/14:63/14:67/14:85/14:89/14:98/14:109/14:115/14:134-14:135/14:144/14:174/14:179/14:194/14:198/14:200/14:209/14:218/14:222N-14:223/14:226/14:239/14:242:01-14:242:02, DRB1*16:08, DRB3*01:08/01:46/01:52/01:71, DRB3*02:06?/02:20	DR18(3), -, DR3, DR11(5), DR13(6), Null, DR14(6), DR1403	190 bp
B10	Mix 7	DRB1*03:49, DRB1*04:01:01:01-04:05:11/04:05:13-04:11:03/04:11:05-04:155/04:157N-04:361, DRB1*08:38, DRB1*11:22/11:102:01-11:102:02/11:175/11:275/11:278, DRB1*14:10:01:01-14:10:01:02/14:57/14:143/14:161/14:182	-, DR4, Null, DR14(6)	85 bp/ 145 bp/ 200 bp/ 260 bp

Position	Mix	Allele	Serology	Size
A10	Mix 8	DRB1*01:21w, DRB1*07:01:01:01-07:01:01:02/07:01:01:04-07:01:30/07:03-07:139/07:140?/07:141-07:146	-, DR7, Null	115 bp/ 190 bp
H11	Mix 9	DRB1*08:01:01:01-08:01:02/08:01:04-08:01:05/08:01:07-08:02:04/08:02:05?/08:03:02:01-08:04:01:05/08:04:03-08:19/08:21-08:30:03/08:32-08:56/08:58-08:66/08:68:01-08:79/08:81-08:119, DRB1*11:67/11:193:01-11:193:02, DRB1*12:16:01-12:16:03/12:22/12:39:01-12:39:02/12:46/12:96, DRB1*14:04:01:01-14:04:09/14:15/14:28/14:31/14:50:01-14:50:02/14:68:01-14:68:02/14:71/14:73/14:76/14:79/14:93/14:107/14:120/14:138/14:145/14:148/14:152N/14:180/14:196/14:201/14:205:01/14:206/14:210Q/14:220/14:224/14:229/14:233/14:241/14:249-14:250	DR8, -, Null, DR1404, DR14(6)	165 bp/ 215 bp/ 250 bp
G11	Mix 10	DRB1*09:01:02:01-09:05/09:07-09:08/09:10-09:22/09:24-09:52N	DR9, -, Null	220 bp
F11	Mix 11	DRB1*01:34w, DRB1*10:01:01:01/10:01:01:03-10:11/10:13-10:30/10:31w/10:32-10:45	-, DR10	205 bp
E11	Mix 12	DRB1*01:75, DRB1*03:08/03:65/03:140, DRB1*04:15/04:61/04:138/04:145/04:154/04:170/04:265/04:291/04:323/04:339, DRB1*08:31/08:41/08:63/08:75, DRB1*10:23, DRB1*11:01:01:01/11:01:01:03-11:02:02/11:02:04-11:04:12/11:04:14-11:11:01/11:11:03-11:21/11:22?/11:23:01-11:25/11:26?/11:27:01-11:39/11:40?/11:41-11:70/11:72-11:85/11:87-11:93/11:95-11:109/11:111-11:122/11:124:01-11:155/11:157-11:166/11:168-11:171/11:173-11:184/11:186/11:188-11:191/11:193:01-11:202/11:203w/11:204-11:207/11:209-11:210/11:211w/11:212-11:253/11:254w/11:255-11:263/11:265-11:292/11:295-11:312, DRB1*12:04/12:49, DRB1*14:11:01:01-14:11:01:02, DRB1*15:99/15:109, DRB1*16:16/16:27, DRB3*01:42, DRB3*02:18, DRB5*01:13/01:41	-, DR4, DR11(5), DR11(5)/DR13 (6), Null, DR14(6)	125 bp
D11	Mix 13	DRB1*08:25/08:32/08:47/08:53/08:87, DRB1*12:01:01:01-12:01:04/12:01:05w/12:01:06-12:19:02/12:20w/12:21-12:30/12:32-12:41/12:43-12:56/12:58-12:82/12:84-12:88/12:89?/12:90-12:103	-, DR12(5), Null	105 bp/ 210 bp
C11	Mix 14	DRB1*03:125, DRB1*11:09:01:01-11:09:01:02/11:16:01:01-11:16:01:02/11:20/11:28:02/11:40w/11:58:01/11:59w/11:83/11:87/11:113/11:115/11:122w/11:182, DRB1*13:01:01:01-13:02:22/13:05:01:01-13:06/13:10/13:15-13:16:01:02/13:18/13:26:01-13:28:02/13:31-13:32/13:34-13:36/13:39-13:42:02/13:50:02-13:51/13:53/13:59/13:61:01-13:61:02/13:63w/13:65/13:67/13:69/13:73-13:74/13:79-13:80/13:85/13:87/13:91-13:93/13:96:01-13:96:02/13:98-13:99/13:102-13:105/13:109-13:110/13:111w/13:112-13:113N/13:117:01:01-13:117:01:02/13:120-13:121/13:123-13:128/13:130-13:131/13:136-13:148/13:153/13:155-13:160/13:165-13:166/13:168/13:170-13:171:02/13:173/13:176-13:177/13:179/13:183-13:187/13:189-13:191/13:198-13:203/13:205/13:207-13:215/13:217-13:218/13:220-13:222/13:223w/13:225-13:226/13:233/13:234w/13:236-13:245/13:248-13:249N/13:251-13:252N/13:254-13:275/13:277/13:279-13:283/13:285/13:287-13:294/13:296w/13:297-13:306/13:308-13:311/13:313/13:315/13:318/13:320/13:322N-13:325/13:327-13:328/13:331-13:332/13:334, DRB1*14:27:01-14:27:02/14:63/14:67/14:78, DRB1*16:08	-, DR11(5), DR11(5)/DR13 (6), DR13(6), Null, DR14(6)	130 bp
B11	Mix 15	DRB1*03:12/03:155/03:158, DRB1*11:03:01-11:03:04/11:11:01/11:11:03/11:36/11:40-11:41/11:55/11:59/11:63:01-11:63:02/11:68/11:76/11:85/11:122/11:138-11:139/11:151/11:161/11:168/11:176/11:213/11:231/11:233/11:237/11:241/11:260/11:301N, DRB1*13:03:01:01-13:04/13:12:01-13:13/13:20:01:01-13:21:02/13:24:01:01-13:24:01:03/13:29:01:01-13:30/13:32-13:33:03/13:38/13:43/13:48-13:49:02/13:54-13:55/13:58/13:63/13:65-13:66:02/13:71/13:75/13:78/13:81/13:88-13:90/13:93-13:95/13:101/13:108/13:111/13:115/13:118/13:120/13:122/13:133-13:134/13:139/13:149/13:151-13:152/13:154/13:161/13:164/13:167/13:169/13:174/13:188/13:194/13:216/13:219/13:227-13:228/13:230/13:247/13:253/13:276/13:295N/13:307/13:317/13:319N/13:321/13:329N/13:333, DRB1*14:13/14:63/14:65/14:78/14:85/14:169-14:170	DR3, -, DR11(5), Null, DR13(6), DR14(6)	170 bp/ 200 bp
A11	Mix 16	DRB1*11:09:01:01-11:09:01:02/11:28:02/11:40w/11:58:01/11:59w/11:83/11:87/11:113/11:115/11:122w/11:182, DRB1*13:05:01:01-13:05:03/13:18/13:26:01-13:26:02/13:42:01-13:42:02/13:50:02-13:50:03/13:63w/13:111w/13:136/13:144/13:158:01-13:158:02/13:189/13:199/13:203/13:277/13:302, DRB1*14:27:01-14:27:02, DRB1*16:08	DR11(5), -, DR13(6), DR14(6)	130 bp

Position	Mix	Allele	Serology	Size
H12	Mix 17	DRB1*03:01:01:01-03:01:06/03:01:07w/03:01:08-03:01:13/03:01:14w/03:01:15-03:07:02/03:09-03:11:01/03:11:03-03:37/03:38w/03:39-03:41:02/03:43-03:45/03:46w/03:47-03:58/03:59w/03:60-03:64/03:66-03:78/03:79w/03:80-03:86/03:88-03:110/03:112-03:126/03:127w/03:128-03:129/03:130w/03:131-03:134/03:135w/03:136-03:139/03:141-03:155/03:157-03:161/03:163-03:186/03:187?/03:188-03:189N/03:190w/03:191-03:205, DRB1*08:01:01:01-08:01:02/08:01:04-08:01:05/08:01:07-08:02:04/08:02:05?/08:03:02:01-08:04:01:05/08:04:03-08:24:02/08:25w/08:26-08:30:03/08:32-08:33/08:34w/08:35-08:40/08:42-08:46/08:47w/08:48-08:51/08:52w/08:53-08:56/08:58-08:62/08:64-08:68:02/08:70-08:74/08:76-08:79/08:81-08:108/08:109w/08:110-08:111/08:112w-08:113w/08:114-08:119, DRB1*11:264/11:293, DRB1*12:01:01:01w-12:01:07w/12:01:09w-12:01:11w/12:01:12?/12:02:01:01w-12:02:05w/12:02:07w-12:03:03w/12:05w-12:08w/12:09/12:10w-12:17w/12:19:01w-12:23w/12:24N-12:25/12:26w-12:30w/12:32w/12:34w-12:41w/12:43w-12:45w/12:46/12:48/12:50w-12:53w/12:54/12:55w-12:56w/12:58w-12:61w/12:62/12:63w-12:66w/12:67/12:68w-12:75w/12:76/12:77w-12:82w/12:84w-12:88w/12:89?/12:90w-12:95w/12:96/12:97w-12:103w, DRB1*13:01:01:01-13:01:40/13:01:41w/13:02:01:01-13:02:11/13:02:12w/13:02:13-13:53/13:54w/13:55-13:66:02/13:68-13:76/13:77w/13:78-13:99/13:100w/13:101-13:102/13:103w/13:104-13:159/13:160w/13:161/13:162w/13:163-13:164/13:165w/13:166-13:182/13:184-13:194/13:196-13:198/13:199w/13:201-13:217/13:218w/13:219-13:227/13:229-13:247/13:250-13:298N/13:300-13:335, DRB1*14:01:01:01-14:09/14:12:01-14:38:02/14:40-14:45/14:47/14:48w/14:49-14:54:01:17/14:54:03-14:54:07/14:54:08w/14:54:09-14:56/14:58-14:60/14:62-14:65/14:67-14:124/14:125w/14:126:01-14:140/14:142/14:144-14:147/14:149-14:176/14:177w/14:178-14:181/14:183-14:194/14:195Nw/14:196-14:198/14:200-14:221/14:223-14:230/14:232-14:246/14:248?/14:249-14:251, DRB3*01:14w	DR17(3), -, DR18(3), DR3, Null, DR8, DR12(5), DR13(6), DR14(6), DR1403, DR1404, DR6	165 bp/ 170 bp/ 175 bp
G12	Mix 18	DRB1*03:10/03:97/03:126/03:191, DRB1*08:08/08:29, DRB1*11:17/11:52-11:54:02/11:89/11:150/11:159/11:202/11:267/11:283/11:293, DRB1*12:66/12:96, DRB1*13:43/13:45/13:76/13:159/13:171:01/13:179/13:188/13:191/13:197/13:206/13:217/13:302, DRB1*14:01:01:01-14:01:05/14:04:01:01-14:05:05/14:07:01-14:08/14:11:01:01-14:11:01:02/14:14/14:16/14:18/14:22-14:23:04/14:25:01-14:26/14:28/14:31-14:33/14:35-14:36/14:38:01-14:39/14:42-14:45/14:49-14:51/14:53-14:56/14:58-14:62/14:64/14:68:01-14:76/14:79/14:82/14:86-14:88/14:90-14:93/14:95-14:97/14:99/14:101/14:103-14:107/14:110-14:114/14:117-14:120/14:122/14:124/14:125w/14:126:01-14:126:02/14:128-14:134/14:137N-14:140/14:142/14:145-14:155/14:157-14:158/14:160-14:164/14:166N-14:169/14:171-14:173/14:175/14:178/14:180/14:184-14:193/14:196-14:197N/14:199-14:208/14:210Q-14:211/14:213-14:217/14:219-14:221/14:224-14:225/14:227-14:230/14:232-14:238/14:240-14:241/14:243-14:246/14:248?/14:249?/14:250-14:251	DR17(3), -, DR11(5), DR14(6), DR1404, DR6, DR13(6), Null	175 bp/ 225 bp
F12	Mix 19	DRB1*03:126/03:157w/03:177, DRB1*04:211, DRB1*11:09:01:01-11:09:01:02/11:28:02w/11:58:01w/11:83/11:87/11:113/11:115w/11:182/11:280, DRB1*12:58, DRB1*13:05:01:01-13:06/13:18/13:26:01-13:26:02/13:42:01-13:42:02/13:50:02w-13:50:03w/13:56/13:77/13:119/13:136/13:144/13:146/13:156-13:158:02/13:163/13:181/13:189w/13:193/13:199/13:203/13:235/13:277/13:302/13:331, DRB1*14:02:01:01-14:03:02/14:06:01-14:06:06/14:09/14:12:01-14:13/14:17-14:18/14:27:01-14:27:02/14:29-14:30/14:33/14:46-14:48/14:51-14:52/14:59/14:63-14:64/14:67/14:78/14:80-14:81/14:83/14:85/14:89/14:95/14:98w/14:106/14:108/14:115/14:121/14:132-14:135/14:144/14:154/14:165/14:170/14:174/14:176-14:177/14:179/14:194-14:195N/14:198/14:200/14:209/14:212/14:218/14:222N-14:223/14:239/14:242:01-14:242:02/14:246, DRB1*16:08w	-, DR11(5), DR13(6), DR14(6), DR1403, DR6, Null	145 bp/ 150 bp
E12	Mix 20	DRB1*09:07, DRB5*01:01:01:01-01:01:09/01:01:11-01:02:01/01:02:02w/01:03-01:25/01:26w/01:27N-01:47/01:48Nw/01:49N-01:50/01:52N-01:72/01:74w/01:75-01:78/01:80-01:86/01:88-01:96/01:98-01:109/01:111-01:112/01:113w/01:114-01:130, DRB5*02:02:01-02:37	-, DR51, Null	175 bp/ 210 bp
D12	Mix 21	DRB1*03:42/03:87, DRB1*14:141, DRB3*01:01:02:01-01:13/01:15-01:67/01:69-01:94/01:96-01:103/01:105N-01:106/01:107?/01:108/01:110-01:112, DRB3*02:01-02:25/02:28-02:46/02:48-02:57/02:59-02:65/02:67N-02:77/02:78w/02:79-02:90/02:92-02:108/02:109Nw/02:110/02:112-02:116/02:117w/02:118-02:123/02:124w/02:125N-02:130/02:132-02:142/02:144-02:159/02:161/02:162w/02:163-02:189/02:190?, DRB3*03:01:01:01-03:01:12/03:01:13w/03:01:14-03:21/03:23-03:45/03:46w/03:47-03:61	-, DR52, Null	215 bp

Position	Mix	Allele	Serology	Size
C12	Mix 22†	DRB4*01:01:01:01/01:01:01:03-01:01:01:04/01:02-01:03:01:01/01:03:01:03-01:03:01:12/01:03:01:14-01:03:02:02/01:152, DRB4*02:01N	DR53, -, Null	150 bp
B12	Mix 23	DRB4*01:03:01:02N/01:03:01:13N	Null	275 bp
A12	Mix 24	Negative Control		None (440 bp)

Amplification control (internal positive control): 440 base pairs (bp)

w = weak

? = nucleotide sequence information not available for the primer matching sequence or alleles with unknown reactivities

Mix 6: HLA-DRB1*03:81 serology unclear

HLA-DRB1*08:04:02 is not listed since the primer matching sequences are unknown

† = Alleles which are not sequenced in the primer binding sites of mixes 22 (according to the IMGT/HLA Database of January 2023) have been excluded from the list of allele specificities of this mix: HLA-DRB4*01:01:02-01:01:10/01:03:03-01:13/01:15-01:61N/01:63-01:87/01:89-01:151/01:153-01:167.

Table 4: Amplification patterns for all **HLA-A*** specificities detected by the HLA-A* CTS-PCR-SSP primer mixes (**Lot No. A17-2**) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*01:01:01:01-01:01:55/01:01:57-01:01:64/01:01:67-01:01:130/01:01:132-01:04:01:02N/01:06-01:12/01:14-01:33/01:35-01:37:01:02/01:39-01:50/01:52:01N-01:71/01:73-01:101/01:103-01:138/01:140:01-01:166/01:168-01:175/01:177-01:193/01:195-01:199/01:201-01:220/01:222-01:243/01:245-01:266:02/01:268-01:283/01:285N-01:379N/01:381-01:389/01:391-01:398/01:400N-01:417N/01:419-01:426, A*11:94/11:271/11:326, A*23:128, A*25:59, B*07:64, B*15:12:01-15:12:03/15:19/15:270/15:298/15:304N/15:649	A1, Null, -, B76(15)	1																							
A*01:01:56	-	1									10														
A*01:01:131/01:399	-	?																							
A*01:13/01:176/01:194	-	1			4																				
A*01:51	-	1												13											
A*01:200/01:244/01:380:01-01:380:01:02, A*11:17, A*36:01:01:01-36:01:02/36:03-36:13	-, A36		2																						
A*01:267	-	w																							
A*01:427, A*03:107/03:164/03:413, A*36:02	-		2		4																				

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*02:01:01:01-02:01:07/02:01:10-02:01:11/02:01:13-02:01:15/02:01:18-02:01:19/02:01:21-02:01:25/02:01:27-02:01:32/02:01:34-02:01:39/02:01:41-02:01:51/02:01:53-02:01:63/02:01:66-02:01:68/02:01:70-02:01:73/02:01:76-02:01:80/02:01:83-02:01:94/02:01:96-02:01:164/02:01:166-02:03:06/02:03:08-02:03:10/02:04:01-02:05:02/02:05:04/02:05:06-02:05:08/02:05:10-02:06:06/02:06:08-02:06:11/02:06:13-02:06:18/02:06:20-02:06:30/02:07:01:01-02:07:11/02:07:13-02:07:14/02:07:16-02:08:01/02:09:01:01-02:11:01:06/02:11:03-02:11:12/02:14-02:18/02:20:01-02:22:02/02:24:01/02:25-02:26/02:28-02:33/02:36/02:39-02:42:01/02:43N/02:45/02:47/02:49:02/02:51/02:53N/02:57-02:59/02:61/02:63/02:66-02:69/02:71/02:73-02:76:01/02:77/02:79:01-02:79:02/02:81-02:87/02:89:01-02:97:02/02:99:01-02:99:02/02:101:01-02:102/02:104/02:106-02:113:02N/02:115-02:116/02:118/02:120-02:121/02:123-02:128/02:130-02:134/02:136-02:140/02:143-02:151/02:153:01-02:168/02:170-02:177/02:179-02:184/02:186-02:187/02:189-02:192/02:194/02:196-02:222N/02:224-02:225N/02:227N-02:228/02:230-02:236/02:238-02:243:01/02:244-02:245/02:247-02:253/02:255-02:270/02:272-02:274/02:276-02:278/02:281-02:285/02:287/02:289:01/02:290-02:297/02:299-02:300/02:302-02:303/02:305N-02:307/02:310/02:313-02:338/02:340-02:347/02:349-02:376/02:378-02:382/02:384-02:405/02:407-02:411/02:413-02:416/02:418-02:453/02:455-02:469/02:471-02:487/02:489-02:490N/02:492-02:503/02:505-02:506N/02:508-02:509/02:511-02:526/02:528:01/02:529-02:542/02:544-02:569/02:572-02:579/02:581Q/02:583-02:590/02:591:02-02:593/02:595-02:600/02:602-02:608N/02:610:01-02:621/02:623-02:629/02:631-02:635/02:637-02:640/02:642/02:645-02:648:02/02:650-02:661/02:663-02:673/02:675N/02:677/02:679-02:680/02:684-02:710N/02:712-02:722/02:724-02:740/02:742-02:780/02:782-02:788N/02:790-02:795Q/02:797N-02:802/02:805Q-02:808/02:810-02:814/02:816-02:825/02:827Q/02:830:01-02:830:02/02:832N-02:839/02:841-02:871N/02:873-02:875/02:877-02:906/02:908-02:922/02:924-02:926/02:928-02:943/02:945N-02:947/02:949-02:958/02:960-02:967/02:969-02:981/02:983-02:984/02:986-02:991/02:993-02:1009/02:1011Q-02:1020/02:1024-02:1027/02:1031N-02:1032/02:1034-02:1045Q/02:1047-02:1049/02:1051-02:1086/02:1088-02:1089	A2, Low A2, -, A203, A210, Null			3																			23		
A*02:01:09/02:05:05/02:06:07/02:07:12/02:11:02/02:76:02-02:76:03/02:243:02-02:243:03/02:528:02/02:570:01	A2, -			3																		21			
A*02:01:12/02:01:33/02:636/02:674/02:815, B*07:02:80	A2, -			3																				w	

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*02:01:17/02:01:26/02:01:52/02:01:64/02:01:69/02:01:74/02:01:81/ 02:01:95/02:03:07/02:05:03/02:06:12/02:06:19/02:07:15/02:08:02/02:55/ 02:64:01-02:64:02/02:72/02:80/02:88N/02:105/02:114/02:117/02:119/ 02:169/02:178/02:185/02:188/02:195/02:223N/02:229/02:246/02:275/ 02:279/02:288/02:289:02/02:298/02:301N/02:304/02:339/02:348/02:377/ 02:383/02:406/02:412/02:491/02:504/02:510/02:527/02:582/02:601/ 02:622N/02:641/02:643N/02:676/02:678/02:681/02:781/02:789N/02:826Q/ 02:828/02:831N/02:872/02:876/02:923/02:944/02:985/02:1028-02:1030/ 02:1046, B*58:01:24	A2, -, Null																							23	
A*02:01:40/02:03:11/02:42:02/02:141/02:312/02:470/02:570:02/02:796N/ 02:803N, B*15:10:07	A2, -, Null			3																					
A*02:01:65/02:01:75/02:05:09/02:60:01-02:60:02/02:193/02:254/02:286/ 02:308/02:311/02:488/02:594/02:649/02:711	-			w																				23	
A*02:01:165	-			3																	20			23	
A*02:06:31	-			?					8	9	10													23	
A*02:12-02:13/02:19/02:27/02:37-02:38/02:44/02:49:01/02:54/02:142/ 02:226N/02:271/02:280/02:417/02:543/02:630/02:682-02:683/02:723/ 02:804/02:840/02:948	A2, -, Null			3											14									23	
A*02:24:02/02:507	-			3																19				23	
A*02:34-02:35:01/02:56:01/02:62/02:907	A2, -			3	w																		22	23	
A*02:35:02	A2				w																		22	w	
A*02:35:03	A2			3	w																	21	22		
A*02:46/02:70/02:1050	A2, -			3			6																	23	
A*02:48/02:129/02:571	-						6																	23	
A*02:50	A2																					21		23	
A*02:52/02:992	-								8															23	
A*02:56:02/02:103	-				w																		22	23	
A*02:65/02:152	A2/19, -																			19				23	
A*02:78	A2			3																			22		

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*02:122/02:591:01, A*11:158, A*68:01:01:01-68:01:10/68:01:12-68:01:24/68:01:26-68:01:27/68:01:29-68:01:35/68:01:37-68:02:06/68:02:08-68:08:02/68:10-68:25/68:27:01-68:27:02/68:29-68:47/68:51-68:70/68:72-68:102/68:104:01-68:128/68:130:01/68:131-68:133/68:135-68:136/68:138-68:171:01N/68:172-68:183/68:185-68:203N/68:205N-68:217/68:219-68:229/68:231-68:233/68:235/68:237-68:248/68:250/68:252-68:257/68:259-68:261/68:263Q-68:285/68:287-68:292N/68:295-68:298, B*44:497	-, A68(28), A28, Null																					21			
A*02:135/02:309/02:454	-								8	9	10														23
A*02:237/02:580/02:809	-			3	w																			23	
A*02:609	-		w	3																				23	
A*02:662	-														14									23	
A*02:741	-			3	w							11											22	w	
A*02:829, A*32:01:01:01-32:01:06/32:01:08-32:01:11/32:01:13-32:01:47/32:01:49-32:01:54/32:03:01:01-32:03:01:02/32:05-32:12/32:14/32:16/32:18-32:21/32:23/32:25-32:27N/32:29/32:31/32:33:01/32:33:03/32:35-32:51/32:53-32:65/32:67-32:70/32:72-32:74/32:76-32:78/32:80-32:92N/32:95-32:100/32:102-32:108/32:110-32:157/32:159-32:169/32:171-32:173	-, A32(19), Null																		18						
A*02:927	-				w																			23	
A*02:959/02:968/02:982/02:1010	-			3																				?	
A*02:1021/02:1023/02:1033/02:1087	-			?																				23	
A*02:1022	-		2	3																				23	
A*03:01:01:01-03:01:38/03:01:40-03:01:118/03:03N-03:09/03:11N-03:23:02/03:25-03:29/03:33-03:64/03:66-03:68N/03:70-03:72/03:74-03:75/03:77-03:81/03:83-03:88/03:89:02/03:91N-03:102/03:104-03:105/03:108-03:112/03:114-03:151/03:153-03:159/03:161N-03:163/03:165-03:166/03:168N-03:170/03:172-03:193/03:195-03:197:02N/03:199/03:201-03:217/03:220-03:222/03:224/03:226-03:235/03:238-03:240/03:243/03:245-03:252/03:254-03:259/03:261-03:270/03:272/03:276-03:280:02/03:283N-03:284N/03:286N-03:308/03:310-03:313/03:315-03:317/03:320-03:332/03:334N-03:342N/03:344:01-03:347/03:349-03:355/03:357N-03:372/03:374N-03:384/03:386-03:412/03:414-03:427/03:429-03:441/03:443-03:452/03:454-03:457	A3, Null, -				4																				
A*03:01:39	-				4											15					20				

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*03:02:01:01-03:02:06/03:31-03:32/03:65/03:69N/03:73/03:76/03:82/03:90/03:106/03:113/03:160/03:198/03:218/03:223/03:236-03:237/03:242/03:244/03:253/03:274-03:275N/03:281/03:285/03:309/03:314/03:318-03:319/03:333/03:348/03:356/03:385/03:453/03:458Q, A*11:130	A3, -, Null				4										w										
A*03:10/03:167/03:225, A*11:199:01/11:222	-				4										14										
A*03:24	A3				4									13											
A*03:30/03:273	-				4		6																		
A*03:89:01	-				4																			23	
A*03:103:01-03:103:02/03:282	-				4								12												
A*03:152	-				4		6											18							
A*03:171/03:271/03:373/03:428/03:442, A*68:71	-				4																21				
A*03:219, A*32:04/32:52	-				4													18							
A*03:241/03:343, A*69:07	-				4																	22			
A*11:01:01:01-11:01:34/11:01:36-11:01:56/11:01:58-11:01:103/11:01:105-11:09/11:11-11:16/11:18-11:22/11:26/11:29-11:30/11:32:01-11:34/11:36-11:37/11:41-11:44/11:46-11:49/11:51-11:52Q/11:54:01-11:56/11:58-11:77/11:79-11:93/11:95-11:97/11:99N-11:110/11:113-11:115N/11:117-11:129/11:131-11:139/11:141-11:157/11:159-11:175/11:177/11:179-11:182Q/11:184-11:189/11:192-11:198/11:200-11:208:02N/11:210N/11:212-11:221/11:224-11:225/11:227-11:228/11:230-11:249/11:251N-11:257/11:259-11:268/11:270/11:272Q-11:289/11:291-11:309/11:311-11:324/11:327-11:345/11:347N-11:358/11:360-11:371/11:373-11:379/11:381-11:388N/11:390-11:391/11:393-11:395/11:398-11:431/11:433N-11:434/11:436N-11:438, A*24:89/24:290, A*68:129	A11, -, Null														14										
A*11:01:104/11:23/11:98/11:229/11:250/11:258/11:346, A*33:19	-														w										
A*11:10	A11													13	14										
A*11:27/11:38-11:39/11:209	-	1													14										
A*11:57/11:359, A*68:09/68:26/68:28/68:134/68:171:02N/68:204	-, A28, Null														14						21				
A*11:116/11:140	-												12		14										
A*11:178/11:190, A*24:02:17/24:02:49/24:02:89/24:08/24:21:01-24:21:02/24:29/24:31/24:42/24:67/24:77/24:156/24:158N/24:191/24:200/24:237/24:273/24:309-24:310:01/24:326/24:342/24:359N/24:370N/24:389N/24:489/24:506/24:562, B*51:55, B*57:129	-, A24(9), Null							7																	
A*11:199:02-11:199:03	-				w										14										

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*11:325, A*24:19/24:28/24:44/24:379/24:406/24:424	-, A9, A24(9)					6									14										
A*11:372, A*31:89, A*34:19, A*66:23/66:30	-												12												
A*11:389/11:392/11:396-11:397N/11:432N/11:435	-, Null														?										
A*23:01:01:01-23:01:20/23:01:22-23:02/23:04-23:08N/23:10:01-23:17:02/ 23:18-23:22/23:25-23:27/23:29-23:30/23:32-23:33/23:35-23:36/23:38N- 23:44/23:46-23:50/23:52-23:54/23:56/23:58-23:60/23:62-23:65/23:67- 23:68/23:70-23:79/23:81-23:82/23:84N-23:99/23:101-23:103N/23:106N- 23:112/23:114-23:127/23:129, A*24:315	A23(9), -, Null					5	6																		
A*23:01:21	-			3		5	6																		
A*23:03:01/23:83, A*32:170	-					5													18						
A*23:03:02:01-23:03:02:02/23:09/23:23/23:28/23:37:01-23:37:02/23:51/ 23:61/23:80/23:100, A*24:24/24:392, B*18:27	-					5																			
A*23:17:03	-					?	6																		
A*23:24/23:34/23:57/23:66/23:104, A*24:02:74/24:02:146/24:04/24:06/ 24:13:01-24:13:02/24:22/24:87/24:107/24:109/24:167/24:207:01- 24:207:02/24:227/24:285/24:289/24:299/24:339/24:355/24:373/24:381/ 24:405/24:446/24:478/24:569/24:572, A*33:21/33:53	-, A24(9), A9						6																		
A*23:31/23:45/23:55	-	1				5	6																		
A*23:105, A*24:71/24:527	-					5	6	7																	

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*24:02:01:01-24:02:05/24:02:07-24:02:16/24:02:18-24:02:48/24:02:50-24:02:63/24:02:65/24:02:67-24:02:73/24:02:75/24:02:77-24:02:88/24:02:90-24:02:145/24:02:147-24:03:04/24:05:01-24:05:02/24:07:01:01-24:07:04/24:09N-24:11N/24:14:01:01-24:14:01:04/24:17:01:01-24:17:01:02/24:20:01:01-24:20:02/24:23/24:25-24:27/24:30/24:32-24:40N/24:43/24:45N-24:50/24:52-24:64/24:66/24:68-24:70/24:72-24:76/24:78-24:86N/24:88/24:90:01N-24:91/24:93/24:95-24:106/24:108/24:110-24:128/24:130-24:137/24:139-24:155N/24:157/24:159-24:166/24:168-24:187/24:189-24:190/24:192-24:199/24:201-24:203/24:205-24:206/24:208:02:01-24:210/24:212/24:214-24:226:02/24:229-24:234/24:236/24:238-24:272/24:274-24:284/24:286-24:288/24:291-24:292/24:294Q-24:298/24:300-24:308/24:311-24:314/24:316-24:325/24:327-24:338/24:341/24:344-24:354/24:356-24:358/24:360-24:369/24:371-24:372/24:374-24:376/24:378Q/24:380/24:382-24:388N/24:390-24:391/24:393-24:404/24:407-24:414/24:416-24:423/24:425N-24:427N/24:429N-24:445N/24:447Q-24:470/24:472-24:477/24:479Q-24:488/24:490-24:498/24:500-24:503/24:505/24:507-24:514N/24:516-24:521/24:523-24:526/24:528-24:550/24:552-24:555/24:558-24:559/24:561:01-24:561:02/24:563-24:568N/24:570-24:571/24:573-24:589	A24(9), Low A24(9), -, A2403, A9, Null						6	7																	
A*24:02:06/24:340	A24(9), -						6	7																23	
A*24:02:64/24:02:66/24:343/24:415/24:428N	-, Null						w	7																	
A*24:02:76	-			3			6	7																	
A*24:15/24:41/24:51/24:92/24:235/24:499/24:504/24:551/24:557	-	1					6	7																	
A*24:18/24:204/24:213, A*31:08, A*32:79	A24(9)/A3, -						6												18						
A*24:21:03, A*32:02/32:22	-, A32(19)							7											18						
A*24:94/24:138/24:188/24:228/24:515, A*68:249	-						6														21				
A*24:208:01	-						6	7											18						
A*24:293	-						6	7													21				
A*24:310:02	-			3				7																	
A*24:377	-						6	w																	
A*24:522	-						?	7																	
A*24:556N/24:560	Null, -						6	?																	

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*25:01:01:01-25:01:01:12/25:01:03-25:01:07/25:01:10-25:01:16/25:01:18-25:01:21/25:03/25:05-25:10/25:12N-25:13/25:15-25:23/25:25-25:32/25:34-25:35/25:37-25:38/25:40-25:42N/25:45-25:58/25:60-25:64/25:66-25:74/25:76:01-25:76:02/25:78-25:80/25:83-25:84	A25(10), -, Null								8	9															
A*25:01:02, B*07:02:74	A25(10), -								w	9															
A*25:01:08/25:39	-								8	w															
A*25:02	A25(10)								8	9			12												
A*25:04/25:44	-							7	8	9															
A*25:11	-								w	9	10														
A*25:14, A*30:12:01/30:18/30:41/30:166, A*32:101Q, 04:01:49, 12:328, 15:222	-, A30(19)								8																
A*25:24/25:33, A*26:01:02/26:01:13/26:01:63/26:110/26:227, A*66:03:01:01-66:03:01:03/66:05/66:15/66:28N/66:39N, B*18:02:02, B*38:02:04	-, A26(10), A10, Null									9															
A*25:36/25:65	-	1							8	9															
A*25:43	-												13												
A*25:75	-								w	9			12												
A*25:77	-					6			8																
A*25:81	-										?														
A*25:82	-					6					?														
A*26:01:01:01-26:01:01:53/26:01:03-26:01:12/26:01:14-26:01:26/26:01:28-26:01:37/26:01:39-26:01:58/26:01:60-26:01:62/26:01:64/26:01:66-26:03:01:02/26:03:03-26:12:01:02/26:14-26:18/26:20:01-26:28/26:30/26:32/26:34-26:43:02/26:45-26:47/26:50-26:59/26:61-26:68/26:70-26:71N/26:73-26:75/26:77-26:83/26:86-26:109/26:111-26:124/26:126-26:130/26:132-26:142/26:145N-26:169/26:171-26:183/26:185-26:203/26:206:01N-26:224/26:228-26:229	A26(10), -, Null, A10									9	10														
A*26:01:27/26:84/26:184	-										10														
A*26:01:38/26:31/26:76/26:170	-									w	10														
A*26:13/26:226	-									9	10		12												
A*26:19	-										10		12												
A*26:29/26:49/26:205	A26(10), -	1								9	10														

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*26:33, A*34:06, A*66:01:01:01-66:02:01:03/66:04/66:06-66:09/66:11-66:14/66:16-66:18/66:20-66:22/66:24-66:27N/66:29/66:31-66:38/66:40/66:42-66:47	-, A66(10), Null									9			12												
A*26:48/26:69	-									9	10	11													
A*26:72	-		2								10														
A*26:143	-									w	w														
A*26:225	-									9	?														
A*29:01:01:01-29:01:06/29:01:08-29:02:35/29:02:37-29:06/29:08N-29:12/29:15-29:18/29:20-29:31/29:33-29:38/29:40-29:47/29:50/29:52-29:59/29:61-29:68/29:70-29:72/29:74-29:104/29:106-29:110/29:112N-29:158/29:160-29:163/29:165-29:172, A*32:30:02, A*74:10/74:41	A29(19), Null, -															15									
A*29:07/29:49	-					5										15									
A*29:13, A*32:30:01/32:32	-															15			18						
A*29:14	-															15		17							
A*29:19/29:39/29:159	-															15					20				
A*29:32	-												12			15									
A*29:48	-															15				19	20				
A*29:51/29:69/29:73	-															15					21				
A*29:105/29:164	-															15				19					
A*29:111	-															w									
A*30:01:01:01-30:01:07/30:01:09-30:02:09/30:02:11-30:04:01:03/30:04:03-30:04:04/30:06/30:08:01-30:11:02/30:14L-30:15/30:17/30:19-30:20/30:23-30:25/30:27N-30:40/30:42-30:43/30:45/30:47-30:54/30:56-30:81/30:83/30:85-30:88/30:91-30:95/30:97-30:123N/30:125-30:151/30:153-30:162/30:165/30:167-30:183/30:185/30:187-30:188/30:190-30:205	A30(19), -, Null								8								16								
A*30:01:08	-								8							w									
A*30:02:10	-								8		10														
A*30:04:02/30:90/30:96/30:124	-															16									
A*30:07/30:13/30:16/30:22/30:44/30:46	-								8							16							w		
A*30:55	-				4				8																
A*30:89/30:186	-				4				8							16									

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*30:152, A*31:01:02:01-31:01:06/31:01:08-31:01:46/31:01:48-31:06/31:09/31:11:01-31:20/31:23/31:25-31:28/31:30-31:74/31:76-31:88/31:91-31:96/31:98-31:103/31:105-31:108/31:110-31:114/31:116-31:135/31:137-31:201N/31:203/31:205-31:209/31:211-31:219	-, A31(19), Null																	17							
A*30:163	-								8					13			16								
A*30:164	-								8						14		16								
A*30:189/30:206	-								?								?								
A*31:01:47/31:202N/31:204	-, Null																	?							
A*31:07/31:10	-, A31(19)																	17	18						
A*31:24/31:97/31:136	-														w			17							
A*31:29	-					5												17							
A*31:104	-																	w							
A*31:115	-												12					17							
A*32:13	-	1																	18						
A*32:15/32:93	-																		18	19					
A*32:17	-				w								12						18						
A*32:24	-		2																18						
A*32:75/32:94	-																		w						
A*32:109	-																		18			21			
A*33:01:01:01-33:01:04/33:01:06-33:01:18/33:03:01:01-33:03:25/33:03:28-33:03:30/33:03:32-33:03:55/33:03:57/33:03:59-33:18:01/33:20/33:22-33:29/33:31/33:33-33:34/33:36-33:37/33:39-33:48/33:50-33:52/33:54-33:58/33:60-33:65/33:67-33:80N/33:82-33:88/33:90-33:101/33:103-33:121/33:123N-33:124/33:126-33:129N/33:132-33:134/33:136-33:139/33:141-33:143N/33:145-33:149/33:151/33:153-33:157N/33:160-33:174N/33:176N/33:178-33:182/33:184-33:186/33:189/33:191-33:214/33:216-33:234, B*07:02:40, B*08:01:07, B*15:02:07, B*55:02:10	A33(19), -, A19, Null, B8																		19						
A*33:03:26/33:03:31/33:66/33:140N/33:158/33:188	-, Null																			w					
A*33:03:56/33:215	-																			?					
A*33:49/33:144	-													13						19					
A*33:59/33:102/33:150/33:177/33:190	-														w					19					
A*33:125/33:131	-												12							19					
A*33:183	-											11								19					

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24
A*33:187	-				5															19					
A*34:01:01:01-34:01:07/34:05/34:11/34:14/34:16-34:18/34:27	A34(10), -									9		11	12												
A*34:02:01:01-34:04/34:07/34:09-34:10N/34:13/34:15/34:20-34:22/34:25-34:26N/34:29N-34:30N	A34(10), -, Null									9				13											
A*34:08	-				4					9				13											
A*34:12	-											11	12												
A*34:23	-									w		11	12												
A*34:24	-									w															
A*34:28	-									9				?											
A*43:01-43:03	A43, Null, -										10														24
A*66:10	-	1								9			12												
A*66:19	-									9			12		14										
A*66:41	-									w			12												
A*68:01:28/68:01:36/68:130:02	-															15					20	21			
A*68:103:01-68:103:02	-													13								21			
A*68:258	-																					?			
A*68:262	-												12									21			
A*68:294	-									9			12									21			
A*69:01:01:01-69:02/69:04-69:06/69:08-69:09	A69(28), -																					21	22		
A*69:03	-														14							21	22		
A*74:01:01:01-74:03:03/74:05-74:09/74:11-74:12N/74:14N-74:15/74:16:02-74:20/74:22-74:40/74:44	A74(19), -, Null																				20				
A*74:04/74:21/74:43	-																			19	20				
A*74:13	-				4								12								20				
A*74:42	-				?								12			15									
A*80:01:01:01-80:01:02/80:03/80:05-80:09N	A80, -, Null																								24
B*07:225/07:260	-								w																
B*35:01:60, B*51:01:34, 06:247, 07:02:27/07:577	-																							w	
B*35:08:14/35:14:03	-													w											
B*57:148	-								?																
*07:1034	-																						?		

w = weak ; ? = nucleotide sequence information not available for the primer matching sequence or alleles with unknown reactivities

Table 5: Amplification patterns for all **HLA-B*** alleles detected by the HLA-B* CTS-PCR-SSP reagents (**Lot No. B16-2**) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023

Allele	Mixes	Serology
B*07:02:01:01-07:02:05/07:02:07-07:02:10/07:02:12-07:02:27/07:02:29-07:02:32/07:02:34-07:02:42/07:02:44-07:02:84/07:02:86/07:02:88-07:02:97/07:07:01:01-07:07:03/07:10/07:15/07:21-07:23/07:26/07:30-07:31:02/07:35-07:36/07:39/07:41-07:44N/07:46-07:49N/07:51-07:52/07:55-07:59/07:61-07:64/07:66-07:68:03/07:70-07:76/07:79/07:81-07:84/07:86-07:89/07:91-07:93/07:95/07:98-07:99/07:101-07:104/07:106-07:111N/07:113-07:114/07:116-07:117/07:119-07:121/07:124-07:125/07:127-07:130/07:135N-07:136:02/07:141-07:142/07:144-07:145/07:147-07:152/07:154-07:161N/07:164/07:166-07:173/07:177-07:179/07:181N/07:184-07:185/07:187-07:192/07:194-07:196/07:198-07:200/07:203-07:205/07:208/07:211-07:212/07:215-07:218/07:221/07:224/07:226-07:227/07:229:01-07:234/07:238-07:242/07:244-07:245/07:250-07:254/07:256-07:257/07:261/07:263/07:265-07:268/07:271-07:277/07:280-07:282/07:285N/07:288-07:289/07:291-07:292/07:294-07:301/07:303:01-07:303:02/07:305/07:307-07:313/07:315N/07:318N-07:323/07:325N-07:326/07:328-07:331/07:334-07:336/07:338-07:339/07:341-07:342/07:344-07:348/07:350/07:354-07:355/07:360Q-07:367/07:369-07:375/07:377/07:379-07:382/07:384-07:385/07:390-07:394/07:396-07:398/07:401-07:406/07:410-07:419/07:421-07:423/07:425N/07:427/07:429-07:433/07:435N-07:441/07:444-07:452/07:454-07:456/07:458-07:459/07:461-07:463/07:465-07:467Q/07:469N-07:470/07:472N	1+24w+28	B7, -, B7/ B42, Null
B*07:02:11	1+24w+28+39+41	B7
B*07:02:28/07:118/07:126/07:132/07:302/07:337	1+24w+28w	-
B*07:02:33/07:02:43	1w+24w+28	-
B*07:02:85/07:33:01-07:33:03/07:54/07:94/07:122/07:139:01-07:139:02/07:165/07:243/07:290/07:343N/07:386N/07:388/07:442/07:468N, B*55:110	24w+28	-, Null
B*07:02:87	1+24w+28?	-
B*07:03/07:08:01-07:08:02/07:16/07:28/07:37:01-07:38/07:85:01-07:85:02/07:96:01-07:96:02/07:180/07:228:01-07:228:02/07:259/07:286/07:306/07:314/07:351N/07:428	1+24w	B703, -, Null
B*07:04:01-07:04:02/07:25/07:146	1+23+24w+28	B7, -
B*07:05:01:01-07:05:02/07:05:04-07:05:05/07:05:07-07:06:02/07:06:04-07:06:05/07:34/07:80/07:90/07:97/07:105/07:112/07:138/07:175-07:176/07:201N/07:206/07:209-07:210/07:213-07:214/07:247/07:249/07:258/07:262:01-07:262:02/07:264/07:269-07:270/07:278/07:283/07:293/07:304/07:316N-07:317/07:324/07:332-07:333/07:340/07:352-07:353/07:356-07:357/07:368/07:378/07:387/07:389/07:395/07:399/07:407-07:408/07:426Q/07:434/07:443/07:457/07:460N/07:474	1+28	B7, -, Null
B*07:05:03	1+28+39+41	B7
B*07:05:06/07:06:03	1w+28	-
B*07:09:01-07:09:02/07:11/07:255	1+21+24w+28w	B7, -
B*07:12/07:18:01	1+7+19+24w+28	B7, -
B*07:13/07:32/07:69/07:115/07:197/07:207/07:225/07:260/07:284, C*02:23, C*03:260, C*04:77, C*07:38:01-07:38:02	1	-, Cw7
B*07:14:01-07:14:02	1+7w+19w+24w	B7, -
B*07:17	1+21+24w	B7
B*07:18:02	1+7w+19w+24w+28	-
B*07:19, B*42:16:01	23+24w+28	-

Allele	Mixes	Serology
B*07:20	1+24w+25+28+30w	-
B*07:24/07:131/07:220/07:223/07:464	1+24w+28+30	B7, -
B*07:27	1+10+15w+24w	-
B*07:29/07:186/07:383	1+2+24w+28	-
B*07:40/07:53/07:123, B*42:19, B*55:01:13/55:04/55:13/55:17/55:28:01/55:49/55:81/55:86/55:90	28	B7, -, B55(22)
B*07:45/07:420	1+12w+13w+24w+28	-
B*07:50	1+15w+24w	-
B*07:60/07:100	24w+28+30	-
B*07:65/07:134	1+5+28	-
B*07:77/07:163/07:193	1+7+24w+28	-
B*07:78/07:182N	1+27+28	-, Null
B*07:133	1+14+17+24w+26+28	-
B*07:137	1+7+19+28	-
B*07:140	1+25+28	-
B*07:143	1+24w+28+32	-
B*07:153	1+24w+25+28	-
B*07:162	1+7+19+21+28	-
B*07:174/07:202/07:279	1+19+24w+28	-
B*07:183	1+15w+24w+28	-
B*07:219	24w+45	-
B*07:222, B*81:01:01:01-81:10	1+19+28	-, B81, Null
B*07:235	7+19+21+24w+28w+30	-
B*07:236	1+10+24w+28	-
B*07:237	1+21+28w	-
B*07:246/07:327/07:376	1+28w	-
B*07:248	1+7+19w+24w+28	-
B*07:287	1+21+28	-
B*07:349, B*42:05:01	1+23+28	-
B*07:358	1+28+30	-
B*07:359	7w+19w+24w+25+30w	-

Allele	Mixes	Serology
B*07:400, B*35:534/35:541, B*40:209/40:350/40:365	5+15	-
B*07:409	1+28+34+43	-
B*07:424/07:473	1+24w+27+28	-
B*07:453, B*48:08	15+24w	-
B*07:471	1?+24w+28	-
B*08:01:01:01-08:01:31/08:01:33-08:01:64/08:01:66-08:03/08:05/08:08N-08:12:02/08:13-08:16/08:18-08:25/08:27/08:29-08:31/08:33-08:34/08:36/08:38:01-08:48/08:51-08:53:02/08:56:01-08:61/08:63-08:78/08:80-08:86N/08:88/08:90-08:93/08:95-08:106/08:108-08:109/08:113-08:122/08:124-08:125/08:127-08:128/08:130/08:134-08:138/08:140-08:145/08:147-08:154/08:157-08:169/08:171-08:179/08:182-08:186/08:188-08:212/08:214N-08:225/08:227-08:231/08:233-08:238/08:240-08:241/08:243-08:254/08:256-08:262/08:264-08:271/08:273/08:275-08:284Q/08:286-08:290/08:292-08:296/08:298/08:300Q-08:305, C*02:14:01-02:14:02/02:107/02:164, C*06:02:37/06:106:01/06:127:02/06:264, C*07:01:07/07:04:03/07:14/07:27:02/07:31:02/07:43:02/07:184/07:199:02/07:257:02/07:257:04/07:294/07:526:01/07:681/07:850, C*12:16:01-12:16:02/12:217	2	B8, -, Null, Cw7
B*08:01:65	2+6	-
B*08:04:01-08:04:02/08:110	2+15w	-
B*08:07/08:87/08:111	2+24w	B8, -
B*08:17/08:54:01-08:54:02/08:285	2+15	-
B*08:28/08:37/08:89	2+21+24w	-
B*08:32/08:133/08:180	2+5	-
B*08:35/08:291, B*15:180	2+21	-
B*08:49, B*39:112, B*53:09, C*15:27	7	-
B*08:55, B*39:12:01:01-39:12:01:02/39:178	2+12	-, B39(16)
B*08:62/08:226/08:242, B*42:26	22	-
B*08:79:01-08:79:02/08:187	1+2	-
B*08:107	6+21+24w	-
B*08:123/08:132	2+23+28	-
B*08:126/08:131/08:139	2+10	-
B*08:129	2+7w+19w	-
B*08:146	2+23	-
B*08:155	2+17	-
B*08:156	1+2+24w	-
B*08:213:01-08:213:02	2+19	-

Allele	Mixes	Serology
B*08:239, B*15:239/15:376, B*18:01:09/18:01:12/18:01:24/18:07:02-18:07:03:02/18:12:02/18:26/18:29/18:41/18:55/18:82/18:97/18:122, B*35:10/35:13/35:16/35:22:01:01-35:22:01:02/35:28:01:01-35:28:01:03/35:31/35:32:02-35:32:03/35:43:01/35:43:03-35:44/35:51:01/35:67/35:69/35:79-35:80/35:86-35:87/35:99/35:102/35:117-35:118/35:176/35:185/35:199/35:208:01-35:208:02/35:226/35:234/35:242/35:251:02/35:265/35:270/35:318/35:475-35:477/35:549, B*37:11/37:14/37:37, C*07:1011	5	-, B35
B*08:255	2+7w+19w+24w	-
B*08:274	2+38+39+41	-
B*08:306	2?	-
B*13:01:01:01-13:01:21/13:06-13:07N/13:11-13:13:02/13:17/13:20/13:23/13:25/13:28-13:29/13:43/13:50/13:52/13:57/13:60-13:61/13:63N/13:73/13:76N-13:78/13:80/13:83/13:87/13:90/13:92/13:95/13:106-13:107/13:109-13:110/13:122/13:129/13:135/13:138/13:141-13:142/13:144/13:146/13:148/13:161N-13:163/13:168/13:170-13:171/13:173	3	B13, -, Null
B*13:02:01:01-13:02:32/13:02:34-13:02:37/13:08/13:10/13:14-13:16/13:18-13:19/13:26:01-13:27/13:30-13:35/13:38/13:41-13:42/13:44/13:47/13:49N/13:53-13:56:02N/13:59/13:64-13:66/13:68-13:70/13:72/13:74-13:75/13:81-13:82/13:84-13:85/13:88-13:89/13:91/13:94/13:96/13:98-13:105/13:111-13:121/13:123Q-13:127/13:130-13:134/13:136-13:137N/13:139N/13:147/13:149-13:152/13:154-13:158/13:160/13:164-13:165/13:167/13:169/13:175-13:176	3+34	B13, -, Null
B*13:02:33	3?+34	-
B*13:02:38	3+34+40	-
B*13:03-13:04/13:40/13:45/13:48/13:128/13:145/13:153/13:159, B*15:01:49/15:01:75/15:24:01-15:24:02/15:36/15:43/15:57/15:87/15:157/15:256/15:289/15:318/15:442/15:579/15:586/15:594, B*46:01:06/46:01:26/46:41N/46:45/46:54/46:79N, B*57:11/57:83, C*01:73	34	-, B15/B21, B62(15), Null
B*13:09/13:67	3+27+34+43	-
B*13:21/13:71/13:108	3+15	-
B*13:36	3+7+20	-
B*13:37/13:143	3w+34	-
B*13:39	3+27	-
B*13:46	11+34	-
B*13:51/13:79	3+11	-
B*13:62	34+40	-
B*13:93	2+3+34	-
B*13:140	3+11+34	-
B*13:166	3+15+34	-
B*13:172	3+34w	-
B*13:174, B*15:360	15+34	-

Allele	Mixes	Serology
B*14:01:01:01-14:01:11/14:01:13/14:04/14:07N/14:12/14:14/14:26/14:28/14:32/14:40/14:46-14:47/14:49/14:54/14:65/14:71/14:77-14:78/14:80-14:83/14:86-14:87/14:95/14:97/14:103-14:105/14:108Q/14:110/14:112-14:113N	46	B64(14), -, Null
B*14:02:01:01-14:02:02/14:02:04-14:02:05/14:02:07-14:03/14:09/14:11/14:15-14:18/14:20/14:22/14:25/14:29-14:31/14:33-14:36/14:38-14:39/14:41N/14:42:02-14:45/14:48/14:50-14:52/14:57-14:64/14:66-14:70Q/14:72N-14:76N/14:79N/14:84-14:85N/14:88-14:92/14:94/14:98-14:101N/14:106-14:107/14:109/14:111/14:114-14:116	46+47	B65(14), -, B14, Null
B*14:05-14:06:02/14:13/14:96	47	-, B14
B*14:10	7w+19w	-
B*14:19	42+46	-
B*14:23	24w+46+47	-
B*14:24	46+47w	-
B*14:27	42+46+47	-
B*14:37	27+47	-
B*14:56	21+24w+47	-
B*14:93	11+46+47	-
B*14:102	7+19	-
B*15:01:01:01-15:01:01:65Q/15:01:03-15:01:04/15:01:06-15:01:46/15:01:48/15:01:51-15:01:70/15:01:72-15:01:74/15:01:76-15:01:78/15:04:01:01-15:06/15:20/15:25:01:01-15:28/15:32:01-15:32:02/15:34:01-15:35:02/15:38:01-15:40:02/15:50/15:56/15:60/15:63/15:65-15:66/15:70/15:75/15:77/15:79N/15:81-15:85/15:92/15:94N/15:96-15:97/15:102/15:104-15:105/15:109-15:111N/15:113/15:117-15:118/15:120/15:122/15:125:01-15:125:02/15:128-15:129/15:135/15:137-15:138/15:140/15:145-15:147/15:149N/15:152/15:154-15:155/15:159-15:160:01:02/15:163-15:167/15:169/15:171-15:172/15:174/15:178-15:179:02/15:181N-15:185/15:187/15:190N/15:192-15:193/15:195/15:199/15:201/15:203/15:205-15:206/15:211/15:217-15:218Q/15:227-15:228/15:231-15:234/15:236-15:237/15:240-15:241/15:244-15:247/15:250-15:251:03/15:257-15:258N/15:260-15:262:02N/15:264/15:267/15:269/15:271-15:272N/15:276/15:279-15:280/15:285-15:286/15:288/15:295-15:296/15:299-15:300/15:309-15:310/15:315-15:317/15:320-15:322/15:326-15:327/15:332-15:333/15:336/15:340/15:343/15:346-15:349:02/15:353/15:355/15:359/15:363:01-15:366/15:368/15:370-15:373/15:375N/15:377Q/15:385-15:387/15:390-15:393/15:395/15:398/15:400N-15:402/15:406-15:407/15:409-15:410/15:412-15:413/15:415/15:417/15:419/15:422/15:428/15:430/15:432/15:440-15:441/15:443/15:452/15:456-15:459/15:461/15:463N-15:474/15:476-15:478/15:480-15:484/15:489-15:490/15:492-15:493/15:496N-15:497/15:499/15:501/15:507/15:510/15:512/15:514-15:515/15:519-15:522/15:528N-15:529/15:531/15:533:01-15:536/15:538/15:540N/15:542-15:545/15:547/15:549N/15:554/15:556/15:558/15:561-15:562N/15:564-15:566/15:568/15:572/15:576-15:577/15:580-15:581/15:583-15:585/15:587/15:589/15:591/15:593/15:595/15:597-15:601/15:604N-15:605/15:608-15:609/15:618/15:621/15:624-15:626/15:631-15:633/15:635-15:636/15:638-15:641/15:645/15:647/15:651-15:652, B*18:178, B*35:207	34+43	B62(15), Null, -
B*15:01:02/15:277-15:278	34+40+41+43	B62(15), -
B*15:01:47/15:01:50, C*07:336	34w	-
B*15:01:71	34+43?	-

Allele	Mixes	Serology
B*15:02:01:01-15:02:04/15:02:06-15:02:12/15:44/15:88/15:121/15:139/15:144/15:194/15:213-15:214/15:223/15:283/15:291/15:301-15:302N/15:319/15:325/15:328/15:330/15:341/15:345/15:357-15:358/15:374/15:378/15:399/15:425/15:437-15:438/15:454N/15:491/15:495/15:537/15:557/15:567/15:569/15:574/15:596N/15:627/15:630/15:648	7+20+34+43	B75(15), -, Null
B*15:02:05/15:112	7w+20+34+43	-, B15
B*15:03:01:01-15:03:13/15:29:01:01-15:29:01:02/15:47:01-15:47:02/15:49/15:52/15:54/15:61/15:64:01/15:74/15:93/15:98/15:103/15:114/15:123/15:127/15:132/15:151/15:156/15:158/15:173/15:220:01:01-15:220:04/15:238/15:242:01-15:242:02/15:259/15:266/15:274/15:281-15:282/15:369/15:379/15:389/15:397/15:444/15:448/15:475/15:502/15:511/15:517/15:527/15:539/15:560/15:582/15:588/15:610/15:616Q-15:617/15:628/15:650, B*18:37:01, B*38:52/38:88, B*39:43, C*07:46/07:356/07:546	21	B72(70), -, B15, B71(70), B70
B*15:07:01:01-15:07:03/15:33/15:45/15:78:01:01-15:78:04/15:107/15:116/15:126/15:141/15:207/15:248/15:287/15:331/15:352/15:383/15:405/15:421/15:431/15:450/15:524	15+34+43	B62(15), -, B15
B*15:08:01:01-15:08:02/15:11:01-15:11:07/15:15:01:01-15:15:01:02/15:31:01:01-15:31:01:03/15:148/15:191/15:209N/15:215/15:284/15:305/15:350/15:367/15:384:01-15:384:02/15:479/15:508/15:513/15:622/15:629	20+34+43	B75(15), -, B62(15), Null
B*15:09:01-15:09:02/15:435	39+41+42	B70, -
B*15:10:01:01-15:10:08/15:37/15:99/15:312-15:313/15:338/15:445/15:455/15:503/15:506/15:563	24w+42	B71(70), -, B70
B*15:12:01/15:12:03/15:19/15:270/15:298/15:304N/15:649	34+43+44	B76(15), -, Null
B*15:12:02	34+43+44w	-
B*15:13:01-15:13:03/15:418	7+8+20+34+37	B77(15), -
B*15:14/15:518/15:551/15:578/15:634	25+34+43+44	B76(15), -
B*15:16:01:01-15:16:05/15:67/15:95/15:254/15:362/15:408/15:550/15:602/15:613	34+37w+45	B63(15), -
B*15:17:01:01-15:17:07/15:168/15:177/15:196/15:208/15:216/15:230/15:273/15:356/15:361/15:403/15:411/15:423-15:424/15:446/15:462/15:500/15:516/15:523/15:532/15:546Q/15:555/15:575N/15:603/15:619/15:642/15:644	45	B63(15), -, Null
B*15:18:01:01-15:18:02/15:18:04-15:18:10/15:23/15:51/15:72/15:115/15:119/15:134/15:153/15:176/15:197:01-15:198/15:200/15:221/15:226N/15:263/15:275:01-15:275:02/15:290/15:292-15:294N/15:306-15:307/15:311/15:314/15:335/15:337/15:351/15:354/15:380N/15:382/15:388/15:414/15:433/15:447/15:451/15:487N/15:498/15:504:01-15:504:02/15:526/15:548/15:570/15:573/15:607/15:612/15:643	21+42	B71(70), -, B70, Null
B*15:18:03, B*56:33	27+42	B71(70), -
B*15:21:01:01-15:21:01:02/15:394/15:404/15:420/15:615/15:623	7+20+34+42+43	B75(15), -
B*15:30:01:01-15:30:01:03/15:58/15:219/15:334/15:416/15:453/15:488/15:541/15:606	43	B62(15), -
B*15:42/15:485/15:552	27+34+43	-
B*15:46/15:53/15:106/15:143/15:525	17+26	B72(70), -
B*15:48	12w+43	B62(15)
B*15:55	15w+20+34+43	B15

Allele	Mixes	Serology
B*15:62/15:124, B*48:02:01-48:02:03/48:25	19+21	-, B48
B*15:64:02, B*18:37:02/18:125	4+21	-
B*15:68/15:614:01:01-15:614:01:02, B*48:14	15+21	B35, -
B*15:69/15:186/15:235, B*39:07/39:91/39:134/39:169	12+21	-
B*15:71/15:175/15:204/15:225/15:344/15:381/15:559	32+34+43	B62(15), -
B*15:73:01:01-15:73:01:02/15:86/15:224	27+43	B62(15), -
B*15:76	20+28w+30+34+43	-
B*15:80	7w+19w+21+42	B70
B*15:89:01-15:89:02/15:505	7+20+34	-
B*15:90/15:133/15:229/15:329/15:427/15:434/15:436/15:439/15:530/15:611/15:620, B*39:33	42	-
B*15:91/15:131	21+25	-
B*15:101	28w+30+34+43	-
B*15:108	12w+21+42	B71(70)
B*15:136/15:249	12w+34+43	-
B*15:142/15:429	2+34+43	-
B*15:150/15:188/15:324, B*40:21/40:137	15+43	-
B*15:161	21+25+42w	-
B*15:170	7+20w+34w+43w	-
B*15:189:01-15:189:02	20+34+42+43+47	-
B*15:202	5+34w+43w	-
B*15:210/15:571	21+32	-
B*15:212, B*49:02/49:70, B*50:01:18/50:09/50:16/50:33/50:41	17+26+32	-, B49(21)
B*15:222	34+37w	-
B*15:243, C*12:03:31, C*14:02:37	41	-
B*15:252	39+41	-
B*15:255, B*56:03/56:50	21+28w+30	-, B22
B*15:265	7+12+20+34+43	-
B*15:268	34+45	-
B*15:297	1+7+20+34+43	-
B*15:303	20+27+34+43	-

Allele	Mixes	Serology
B*15:308, B*35:01:17/35:01:67/35:03:07/35:04:03/35:12:02/35:20:01/35:59:02/35:83/35:88/35:95/35:112:01-35:112:02/35:164/35:191/35:212/35:215:01-35:215:02/35:230/35:261/35:296/35:382/35:414/35:439:01:01-35:439:01:02/35:463/35:474/35:547, B*37:67	5+7	-, B35
B*15:323	24w	-
B*15:339	18	-
B*15:342	18+43	-
B*15:426	21+32+42	-
B*15:449/15:494	12+34+43	-
B*15:460	7w+20+43	-
B*15:486	20+27+43	-
B*15:509	21+30	-
B*15:590	34?+43w	-
B*15:592	7+20+25+34+43+44	-
B*15:637	20+43	-
B*15:646	7+20+42+43	-
B*18:01:01:01-18:01:08/18:01:10-18:01:11/18:01:13-18:01:23/18:01:25-18:03:01:03/18:04:01-18:07:01/18:08/18:10-18:12:01/18:13/18:17N/18:20/18:23N-18:25/18:27-18:28/18:31-18:34/18:36/18:38-18:39/18:42-18:47/18:50-18:53/18:58-18:60/18:62-18:65/18:68-18:78/18:80-18:81/18:83-18:96/18:98:01-18:98:02/18:100/18:102-18:104/18:106Q/18:109/18:111-18:121/18:123-18:124/18:126-18:131:01:02/18:133-18:135:01:02/18:137-18:141/18:144-18:148/18:151/18:153-18:163/18:166-18:167/18:170-18:177/18:179-18:188/18:190/18:192-18:193/18:196-18:197/18:199-18:203/18:206-18:216/18:218-18:226	4+5	B18, -, Null
B*18:03:02/18:09/18:49/18:66-18:67/18:79/18:108/18:132/18:136/18:142/18:194:01:01-18:194:01:02/18:198	4	-, B18
B*18:14/18:22/18:40/18:217	4+5+15w	-
B*18:15/18:19/18:30/18:57:01-18:57:02/18:101/18:152/18:164	4+5+32	-
B*18:18:01:01-18:18:01:02/18:105/18:204, B*35:20:02	4+5+7	-, B35
B*18:21	4+5+15w+32	-
B*18:35/18:110	1w+4+5	-
B*18:48, B*35:63/35:317, B*44:10/44:14-44:15:01:02/44:17/44:43:01-44:43:02/44:77/44:82/44:140/44:144/44:162/44:226/44:230/44:256/44:276/44:308/44:339-44:340/44:541	17	-, B44(12), B12
B*18:54/18:195	4+10	-
B*18:56/18:169	4+5+25	-
B*18:61/18:189, B*35:72	5+15w	-
B*18:99	4w+5	-
B*18:107	4w+17	-

Allele	Mixes	Serology
B*18:143	4+5+20w+34w+43w	-
B*18:149-18:150	4+5+7w	-
B*18:165	5+18w	-
B*18:168, B*35:276, B*39:19:01-39:19:02/39:186	5+12	-, B39(16)
B*18:191, B*35:60	5+27	-, B35
B*18:205	4+5+39+40w+41	-
B*27:01-27:03/27:05:02:01-27:05:22/27:05:24-27:05:59/27:07:01-27:09/27:11-27:14/27:16-27:17/27:23/27:26-27:28/27:31-27:33/27:35/27:37-27:39/27:42-27:53/27:55-27:62/27:64N-27:65N/27:67/27:70-27:74/27:76-27:77/27:80-27:82/27:84-27:85/27:87-27:90:04/27:93-27:99/27:101-27:102/27:104/27:110-27:111/27:116-27:119/27:121-27:124/27:126/27:129-27:135/27:137/27:139/27:141-27:146/27:148/27:150-27:152/27:155-27:164/27:166-27:167/27:169-27:172/27:174-27:182/27:184-27:188/27:190-27:191/27:193-27:207/27:209/27:212N-27:213/27:216-27:221/27:223N-27:225N/27:228/27:230-27:238/27:241-27:244/27:246N/27:248-27:249/27:251-27:256/27:258-27:261N/27:263N-27:265	9	B27, -, B2708, Null
B*27:04:01-27:04:06/27:10/27:15/27:20/27:25/27:40/27:54/27:63/27:66N/27:68-27:69/27:79/27:86/27:92/27:100/27:103/27:105/27:108-27:109/27:112-27:115/27:120/27:147/27:149/27:165/27:168/27:183/27:210-27:211/27:214/27:226/27:229/27:240/27:245/27:247/27:257	9+32	B27, -, Null
B*27:06:01:01-27:06:02/27:21:01-27:21:02/27:91/27:106/27:138/27:192/27:208/27:222/27:262	9+24w+32	B27, -
B*27:18	9+18+32	-
B*27:19:01:01-27:19:01:02	9+19w	B27, -
B*27:24/27:173	9+16+19+32	-
B*27:29	9+18	B27
B*27:30	8w+9+19w+37w	B27
B*27:34/27:136/27:154/27:227	9+24w	-
B*27:36/27:125	9+16+19	-
B*27:41/27:78	9+21	-
B*27:83/27:128	9+22	-
B*27:107	9+21+32	-
B*27:127	9+19w+24w	-
B*27:153	32+41	-

Allele	Mixes	Serology
B*35:01:01:01-35:01:03/35:01:05-35:01:09/35:01:11-35:01:16/35:01:18-35:01:55/35:01:58-35:01:63/35:01:65-35:01:66/35:01:68-35:01:72/35:02:01:01-35:02:05/35:02:07-35:02:16/35:02:18-35:03:02/35:03:04-35:03:06/35:03:08-35:03:16/35:03:18-35:03:19/35:03:21-35:03:31/35:03:33-35:04:01:02/35:06-35:08:01:13/35:08:04/35:08:06-35:09:04/35:11:01-35:11:03/35:14:01-35:15:02/35:18/35:21/35:24:01/35:29:01-35:29:02/35:33-35:34/35:36/35:38/35:40N-35:42:02/35:48/35:50/35:52-35:55/35:57/35:62/35:64:01-35:65Q/35:70/35:74/35:77-35:78/35:81/35:84/35:90-35:94/35:96/35:98/35:100-35:101:02/35:103-35:111/35:115-35:116/35:119/35:121-35:124/35:126-35:134N/35:136-35:148/35:150:01-35:150:02/35:152-35:153/35:155-35:163/35:165N-35:168/35:170-35:175/35:177-35:184/35:187-35:190/35:192-35:197/35:200/35:204-35:206/35:209/35:211/35:213-35:214/35:216N/35:218-35:220/35:222-35:225/35:227-35:228/35:231/35:237-35:241/35:243-35:250/35:252-35:253/35:255-35:257/35:260/35:262-35:264/35:267-35:269/35:272/35:275/35:277-35:285/35:287/35:289-35:290/35:293/35:295/35:297-35:298/35:300-35:308/35:310/35:312-35:316/35:319-35:354/35:356-35:360/35:362-35:366/35:369-35:372/35:374-35:381N/35:383-35:388/35:390N-35:394/35:396-35:413/35:415/35:417-35:424/35:426/35:428Q-35:431/35:433-35:438/35:441-35:445/35:447-35:449/35:451/35:453N-35:458/35:460Q-35:462/35:464-35:467/35:469-35:470/35:473/35:478-35:480/35:482-35:497/35:500-35:508N/35:510-35:518/35:520/35:524/35:526-35:533/35:536-35:540/35:542N-35:544/35:546/35:548/35:552-35:559/35:562-35:569/35:573-35:574N	5+6+7	B35, -, Null
B*35:01:04/35:01:73/35:02:17/35:03:20/35:08:03/35:12:04/35:23/35:24:02/35:61:02/35:120/35:258/35:288/35:395	5+6+7w	B35, -
B*35:01:10/35:03:32/35:04:02	5+6+7+38+39+41	B35, -
B*35:01:56/35:02:06	5+6w+7	-
B*35:01:57	4+5+6+7	-
B*35:01:64/35:03:17/35:29:03	5+6+7+47	-
B*35:03:03/35:08:02/35:08:05/35:27/35:56/35:203/35:210/35:271/35:299/35:355/35:550-35:551N, B*53:03/53:11-53:13/53:36/53:55/53:65	6+7	B35, -, Null, B53
B*35:05:01:01-35:05:06/35:17:01:01-35:17:05/35:30:01:01-35:30:01:02/35:32:01/35:37/35:51:02/35:58/35:68:01-35:68:02/35:75/35:89/35:97/35:113-35:114/35:125/35:151:01-35:151:02/35:221/35:232/35:235/35:273/35:286/35:292/35:367/35:389/35:416/35:427N/35:440/35:459N/35:472/35:481/35:498-35:499/35:509/35:519/35:521/35:523/35:535/35:572	5+6	B35, -, Null
B*35:12:01:01-35:12:01:04/35:12:03/35:39/35:59:01/35:61:01/35:149/35:201/35:233/35:266/35:425/35:571	5+7w	B35, -
B*35:19/35:47/35:154/35:217/35:274/35:450	6+7+17	B35, -
B*35:25/35:49:01-35:49:02/35:254	6+7+19+21	B35, -
B*35:26	7+19+21+42+47	-
B*35:35/35:468	5+6+7+12	B35, -
B*35:45/35:71/35:169/35:309/35:452	5+6+7+25	-
B*35:46	6+7+20+34+43	B35
B*35:66	1+5+6	-
B*35:76/35:561	5+7+28w+30	B35/22, -
B*35:82/35:85/35:294	5+7+42	-
B*35:135	5+42	-
B*35:186, B*46:40:01-46:40:02	5+33	-

Allele	Mixes	Serology
B*35:198	6+7+19	-
B*35:202/35:373	5+6w	-
B*35:229	5+6+7+28w+30w	-
B*35:236/35:259	5+6+7+12w	-
B*35:251:01, B*78:01:01:02-78:02:02/78:04/78:07-78:08	5+38+39+41	-, B78
B*35:291/35:446	5+6+7+27	-
B*35:311	6+7+20+43	-
B*35:361, B*48:31, C*15:99	19	-
B*35:368	5+6+28w+30	-
B*35:432	5+7+27	-
B*35:471, B*37:10/37:34, B*38:01:01:01-38:02:03/38:02:05-38:02:09/38:02:11-38:04/38:06-38:16/38:18-38:19/38:21-38:25/38:27-38:29/38:31-38:32/38:34N-38:51/38:53-38:62/38:64/38:66-38:72/38:74-38:83N/38:85-38:87/38:89-38:104/38:160-38:173Q, B*39:104, B*57:118, A*24:96/24:256	11	-, B37, B38(16), B16, Null
B*35:522	5?+6+7	-
B*35:525	5?+6	-
B*35:545	5+6+7+25+27	-
B*35:560/35:575	5+6?+7	-
B*35:570, B*46:04/46:47/46:61/46:70/46:74/46:77/46:86, C*01:32:01-01:32:02, C*02:56, C*03:263:01-03:263:03/03:514, C*04:180:01-04:180:02, C*05:217, C*06:20, C*07:81/07:168/07:896, C*08:123/08:139, C*14:82, C*15:126, 16:98/16:102/16:110	33	-, B46
B*37:01:01:01-37:01:11/37:01:13-37:01:29/37:03N-37:04:02/37:13/37:15-37:21/37:23-37:28/37:30N-37:33N/37:35-37:36/37:38-37:66/37:68-37:80/37:82N-37:94/37:96/37:98-37:101/37:103-37:104/37:106-37:107, B*38:17, A*11:133, A*33:60, A*68:234	10	B37, -, Null
B*37:02/37:102	10+18	-
B*37:06:01-37:06:02	10+17	-
B*37:07	10+15+32	-
B*37:08	7+10	B37
B*37:09/37:12/37:81	10+15	B37, -
B*37:97	10+32	-
B*37:105	10?	-
B*38:02:04/38:02:10	11w	-
B*38:05/38:33	11+47	B38(16), -
B*38:20:01	7w+8w+11+19w+37w	-
B*38:20:02	7+8+11+19+37	-

Allele	Mixes	Serology
B*38:26	11+13+28	-
B*38:30/38:65	11+42	-
B*38:63	1+11+19	-
B*38:84	2+11	-
B*39:01:01:01-39:01:01:29/39:01:03:01-39:01:10/39:01:12-39:03:02/39:05:01:01-39:06:04/39:06:06-39:06:08/39:08/39:09:02/39:10:01-39:11/39:13:01-39:16/39:18/39:20/39:22-39:24:01/39:24:03-39:28/39:31:01:01-39:31:01:02/39:35/39:37:01:01-39:41/39:44-39:46/39:48-39:49/39:51-39:68/39:70-39:73/39:75-39:79/39:81/39:83:01-39:87N/39:89-39:90/39:92-39:101/39:103/39:105-39:106/39:108-39:111/39:113-39:126/39:130-39:133N/39:135/39:137-39:154/39:157N/39:159-39:168/39:172/39:175N-39:177/39:179-39:184/39:187-39:194Q, C*07:137:01-07:137:02/07:868, C*14:63	12	B3901, Low B39, -, B3902, B39(16), Null
B*39:04/39:129	12+47	B39(16), -
B*39:09:01:01-39:09:01:06/39:09:03-39:09:05/39:88/39:127/39:155/39:158/39:170-39:171	7+12	B39(16), -
B*39:17	12+27	-
B*39:29	12+24w	-
B*39:32/39:74	12w+42	-
B*39:42	7w+12+19w	-
B*39:47/39:50/39:102/39:174, C*03:215	12w	-
B*39:69, B*40:184	12+17	-
B*39:173	12+25	-
B*39:185	12?	-
B*40:01:01-40:01:05/40:01:07-40:01:08/40:01:10-40:01:18/40:01:20-40:01:55/40:01:57-40:01:80/40:22N/40:30/40:34/40:36/40:38/40:43/40:45/40:47/40:52/40:54-40:55/40:59-40:60/40:62/40:65-40:67/40:69/40:73-40:74/40:76/40:79-40:81/40:84/40:88/40:100-40:102/40:108/40:112-40:114:02/40:116-40:118N/40:123/40:125:01-40:126/40:128/40:130:01-40:130:02/40:134-40:135/40:138-40:141/40:146-40:147/40:149-40:156/40:160:01-40:160:02/40:168/40:170/40:172/40:175/40:179/40:182/40:186:01-40:188/40:192-40:193/40:195/40:204/40:207/40:212-40:213/40:215-40:218/40:221/40:223/40:227-40:228/40:233-40:236/40:238-40:242/40:245:01-40:245:02/40:247/40:249-40:253/40:257-40:264/40:272-40:273/40:277-40:282/40:285-40:286N/40:288/40:298:01-40:299/40:301/40:310/40:312/40:316/40:319/40:321/40:323:01-40:323:02/40:327/40:329/40:332-40:333/40:336/40:338N-40:339/40:344/40:346/40:352-40:353:02/40:357-40:358/40:360/40:367/40:370/40:372N/40:379:01:01-40:379:02/40:381:01-40:383/40:385-40:386/40:391-40:398/40:400-40:404/40:406-40:407/40:416-40:419/40:421Q/40:426N/40:428N/40:431/40:433-40:434/40:436/40:439/40:443/40:446-40:448/40:450/40:452/40:456/40:461/40:463/40:469/40:474-40:475:02/40:479-40:480/40:483N/40:485-40:486/40:490-40:491/40:495/40:497/40:499/40:501-40:503/40:506N-40:510N/40:512/40:516/40:524-40:525/40:528	14+15+16	B60(40), -, Null, B40
B*40:01:06	14+15+16+40+41	B60(40)
B*40:01:09/40:07/40:61/40:87:01-40:87:02/40:196/40:222/40:300/40:368/40:422/40:517	14+16	-, B60(40)
B*40:01:19/40:191/40:197	14+15+16w	-
B*40:01:56/40:171/40:315	14w+15+16	-

Allele	Mixes	Serology
B*40:02:01:01-40:04:02/40:06:01:01-40:06:29/40:06:31/40:11:01:01-40:11:02/40:14:01-40:14:03/40:19-40:20:01:02/40:27:01-40:27:03/40:29/40:35:01-40:35:03/40:37/40:40/40:50:01:01-40:50:01:02/40:53/40:56-40:58/40:64:01:01-40:64:02/40:70:01-40:70:02/40:72:01-40:72:03/40:78/40:82-40:83/40:86/40:89-40:91/40:94/40:96-40:97/40:99/40:103-40:105/40:107/40:110-40:111/40:115/40:119/40:122/40:127/40:129/40:131/40:133Q/40:136/40:142N-40:145/40:157/40:159/40:164-40:165/40:167/40:169/40:173/40:176/40:181/40:189-40:190/40:200-40:203/40:205-40:206/40:211/40:214/40:219/40:224-40:226/40:229-40:230/40:243-40:244/40:246/40:248/40:254-40:256N/40:266-40:271/40:274-40:276/40:283-40:284/40:287/40:289-40:293/40:295-40:297/40:302-40:307/40:311/40:314/40:317-40:318/40:320/40:330-40:331/40:334:01-40:335/40:337N/40:341-40:342/40:345N/40:347-40:349/40:351/40:354/40:356/40:359/40:361N-40:362/40:364/40:366/40:373-40:376/40:384/40:387-40:390/40:399N/40:405/40:408-40:413/40:415/40:420/40:423/40:425/40:427/40:432/40:435/40:437-40:438N/40:440-40:442/40:444-40:445/40:449/40:451/40:453-40:455/40:457/40:459-40:460/40:464-40:465/40:467-40:468/40:470-40:473/40:476/40:478/40:482/40:484/40:487N-40:489/40:492-40:494/40:500/40:504-40:505/40:511N/40:514N-40:515/40:518-40:519/40:521/40:526-40:527	14+15+17	B61(40), -, B60(40), B40, Null
B*40:05:01:01-40:05:02/40:174/40:496	15+17+26+32	B4005, -
B*40:06:30	14+15+17+40+41	-
B*40:08:01:01-40:08:01:02/40:13/40:75/40:77/40:121/40:162/40:220/40:232/40:237/40:322/40:363/40:369/40:481N/40:520	14+17	-, Null
B*40:09/40:18/40:24/40:44/40:120	14+15+17+24w	B61(40), -
B*40:10:01:01-40:10:02/40:49/40:85/40:194/40:210/40:231/40:325/40:380/40:462/40:477, B*44:31, B*48:41	14+15	B60(40), -, B44(12)
B*40:12:01-40:12:02, B*48:01:01:01-48:01:07/48:01:09-48:01:12/48:03:01:01-48:04:02/48:07/48:09/48:11/48:15-48:21/48:24/48:27-48:29/48:32-48:40/48:42-48:45/48:47-48:57	15+19	-, B48, B48/40
B*40:15-40:16:01:02/40:148/40:161/40:309	14+15+17+26+32	-, B61(40)
B*40:23/40:183/40:198/40:466	14+15+16+26+32	-
B*40:25	1+14+15w+16	-
B*40:26/40:28, B*50:14	17+26+32+40+41	B21, -
B*40:31/40:33/40:42	14+15+16+24w	B60(40), -
B*40:32/40:98	14+15+17+24w+26+32	-
B*40:39/40:343, B*41:08-41:09/41:33/41:47/41:62, B*44:166/44:184/44:213,C*16:85	15+17	B41, -
B*40:46/40:430, B*48:22/48:30	15+16+19	-
B*40:48/40:377	14+15+16+27	B60(40), -
B*40:51	15+16+26+32	-
B*40:63/40:92/40:163/40:199/40:208/40:265N/40:324	15+16	-, Null
B*40:68	7+14+17	-
B*40:71	15+17+27	-
B*40:93, B*42:11, B*48:05/48:10/48:12-48:13, C*07:937	15	-, B48
B*40:95	14+17+26+32+40+41	-

Allele	Mixes	Serology
B*40:106	1+14+16	-
B*40:109/40:340	14+15+17+37	-
B*40:124:01-40:124:02	9+14+15+16+18	-
B*40:132/40:178	14+15+16+25	-
B*40:158	14+17+26+32	-
B*40:166	1+14+16+28	-
B*40:177/40:294/40:378	14w+15+17w	-
B*40:180	9+14+15+17+18	-
B*40:185/40:308	5+15+16	-
B*40:313	14+15w+17	-
B*40:326	14+15w+16	-
B*40:328	2+14+15	-
B*40:355	14+15+26+32	-
B*40:371	14+15+17w	-
B*40:414	14+15+17+24	-
B*40:424	12+15+16	-
B*40:429, B*47:01:01:02-47:07/47:09/47:11-47:12	9+14+17+18	-, B47
B*40:458, B*44:02:10/44:09/44:34:01/44:46:01-44:46:02/44:129/44:131/44:254/44:277/44:438N	17+24+25	-, B44(12), B45(12), Null
B*40:498	14+15+16+24	-
B*40:513	14+16+28w	-
B*40:522	14+15+17+27	-
B*40:523	15+16+25	-
B*41:01:01:01-41:02:08/41:02:11-41:03:02/41:05-41:07/41:10-41:18/41:22-41:32/41:34-41:45N/41:48-41:59/41:61/41:64-41:79	15+17+22	B41, -, Null
B*41:02:09	15+17+22+40+41	-
B*41:02:10	15+17?+22	-
B*41:04/41:21/41:63	15+17+22+24w	-
B*41:19	14+15+17+22	-
B*41:20	15+17+22w	-
B*41:46	17+22	-

Allele	Mixes	Serology
B*41:60	15+17+22+37w	-
B*42:01:01:01-42:01:06/42:04/42:06/42:08/42:10/42:12/42:14-42:15/42:21-42:25/42:28-42:31/42:33, B*55:20	23+28	B42, -
B*42:02:01:01-42:02:02/42:17-42:18/42:27	22+23+28	B42, -
B*42:05:02	1w+23+28	-
B*42:07/42:20	23	-
B*42:09	21+22+23+28w	-
B*42:13	7+19+23+28	-
B*42:16:02	23+24w+28w	-
B*42:32	21+23+28w	-
B*44:02:01:01-44:02:09/44:02:11-44:02:12/44:02:14-44:02:69/44:02:71-44:03:11/44:03:13-44:03:36/44:03:38-44:03:39/44:03:41-44:03:50/44:03:52-44:04/44:07/44:11/44:13/44:19N-44:20/44:22-44:24/44:26-44:30/44:32-44:33/44:34:02-44:39/44:41:01-44:42/44:45/44:47-44:48/44:50:02/44:51/44:55-44:56N/44:58N-44:59:02/44:61N/44:63/44:65-44:69:02/44:71-44:74/44:80-44:81/44:84:01-44:89/44:91-44:94/44:96/44:98-44:105/44:108N-44:109/44:111-44:116/44:118-44:122/44:125-44:128:02/44:132-44:133/44:137-44:138Q/44:141-44:143/44:147-44:149N/44:151-44:155/44:157/44:159-44:161/44:163:01-44:165/44:167-44:168/44:170-44:181/44:183/44:185-44:191/44:193-44:197/44:199-44:200/44:202-44:209/44:212/44:214-44:216/44:218-44:225/44:227-44:229/44:231/44:234-44:235/44:237N-44:245:02/44:247-44:249/44:252-44:253/44:255/44:257-44:262/44:265-44:275/44:278-44:287/44:289-44:307/44:309N-44:329/44:331-44:338/44:341N-44:342/44:344-44:346/44:348-44:354N/44:357-44:364/44:437/44:439-44:444/44:446-44:458/44:460-44:464:01:02/44:466N-44:505/44:507-44:524/44:526-44:529/44:531-44:540/44:542-44:543/44:545-44:546	17+24	B44(12), -, Null, B21
B*44:02:70/44:05:01:01-44:05:06/44:70/44:78/44:107/44:123/44:136/44:139/44:198N/44:217N/44:236/44:250/44:264/44:288/44:330/44:465	17+24w	-, B44(12), Null
B*44:03:12	4w+17+24	-
B*44:03:40	17+24+40	-
B*44:06, B*51:42/51:359	7w+8w+25+37w	B44(12), -
B*44:08/44:57/44:60	34+44	B44(12), -
B*44:12/44:232/44:365	7w+17+24	B44(12), -
B*44:16/44:21/44:64:01-44:64:02/44:124/44:347	14+17+24	B47, -
B*44:18:01:01-44:18:01:02/44:355, B*49:01:02/49:11/49:45	17+31+37w	-, B49(21)
B*44:25/44:50:01/44:95	8w+17+24+37w	-
B*44:40/44:44/44:130	19w+24	B44(12), -
B*44:49	2+17+24	-
B*44:50:03	8+17+24+37	-
B*44:52N/44:169/44:233	24	Null, -
B*44:53:01-44:54/44:106/44:110/44:135/44:192:01-44:192:04/44:211/44:343/44:506	15+17+24	B44(12), -

Allele	Mixes	Serology
B*44:62	17+40	-
B*44:75	17+24w+25	-
B*44:76/44:79:01-44:79:02/44:146	17+24+26	-
B*44:90	17+23+24+25+30	-
B*44:97	9+17+24	-
B*44:117/44:251	11+17	-
B*44:145/44:530	17w+24w	-
B*44:150	14+17+24+26	-
B*44:156	2+19w+24w	-
B*44:158	15+17+24w	-
B*44:182	15+24	-
B*44:210:01-44:210:02	19w+24w	-
B*44:263	7+17+24w+42w	-
B*44:356	7w+17+23?+24+25+30w	-
B*44:459	17+26+40	-
B*44:544N	17?+24	Null
B*45:01:01:01-45:01:12/45:03/45:05/45:07/45:11-45:21/45:23/45:25/45:28N/45:30	17+25+27	B45(12), -, Null
B*45:02/45:09/45:29	17+25	-
B*45:04/45:08/45:24, B*50:70, B*56:67	17+27	-
B*45:06	17+23+25+27+28+30w	-
B*45:10	14+17+27	-
B*45:22	5+25+27	-
B*45:27	15+17+25+27	-
B*46:01:01:01-46:01:04/46:01:07-46:01:25/46:01:28-46:03/46:05/46:07N-46:10/46:12-46:13:03/46:15N-46:16/46:20-46:30/46:32:01-46:39/46:42-46:44/46:46/46:48-46:53/46:55N-46:60/46:62-46:69/46:71/46:73/46:75-46:76/46:78/46:80-46:85/46:88-46:92/46:94-46:98N	33+34	B46, -, Null
B*46:01:05/46:31	33w+34	-
B*46:01:27	33+34w	-
B*46:06	1+34	-
B*46:11/46:18	27+33+34	-
B*46:14	32+33+34	-

Allele	Mixes	Serology
B*46:17	25+33+34+44	-
B*46:19/46:93	12w+33+34	-
B*46:72	33+34+43	-
B*46:87	1+33+34	-
B*47:08	9+17+18	-
B*47:10	9+14+17+18+26	-
B*48:01:08	15+19+40+41	-
B*48:06	1+19	-
B*48:23	15+19+21	-
B*48:26	18+19	-
B*48:46	15+19+24w	-
B*49:01:01:01-49:01:01:03/49:01:01:05-49:01:01:18/49:01:04-49:01:17/49:01:19-49:01:24/49:06/49:08-49:09/49:12-49:14/49:16-49:17/49:19N-49:23/49:27-49:30/49:32-49:39/49:41-49:44/49:46/49:48-49:49/49:51-49:53/49:55-49:60N/49:62-49:66/49:68-49:69/49:74-49:75/49:77-49:78	17+26+31+32+37w	B49(21), -, Null
B*49:01:01:04/49:54	17+26+31+32+37	-
B*49:01:03	26+31+32+37w	-
B*49:01:18	17w+26+31+32+37w	-
B*49:03/49:73	31+32+37w	-
B*49:04:01-49:05:02/49:18:02	17+26+32+37w	-, B49(21)
B*49:07	14+17+26+31+32+37w	-
B*49:10/49:24/49:31	17+26+31+37w	-
B*49:18:01	17+26+32+37w+40	-
B*49:25	2+17+26+31+37w	-
B*49:26	31+32+37w+42+47	-
B*49:40	15+17+26+31+32+37w	-
B*49:47:01-49:47:02	17+26+31+32	-
B*49:50	11+17+26+32+37w	-
B*49:61	8w+17+24+26+32+37w	-
B*49:67	6+17+26+31+32+37w	-
B*49:71-49:72	17+26+31+32+37?	-
B*49:76	17+22+26+31+32+37w	-

Allele	Mixes	Serology
B*50:01:01:01-50:01:17/50:01:19-50:01:20/50:01:22-50:01:23/50:04:01-50:05/50:07-50:08/50:10-50:12/50:15/50:17-50:19/50:22/50:32/50:34-50:40/50:42/50:44-50:48/50:50-50:55/50:57-50:59/50:61-50:69/50:71-50:74/50:76-50:79	17+26+27+32	B50(21), -, Null
B*50:01:21	17+26+27+32+40+41	-
B*50:02:01:01-50:02:01:02/50:21/50:56	17+25+26+27+32	B45(12), -
B*50:06/50:75	27+32	-
B*50:13/50:43	17+26+27	-
B*50:20	9+17+18+26+32	-
B*50:31	14+17+26+27+32	-
B*50:49	12+17+27	-
B*50:60	15+17+26+27+32	-
B*51:01:01:01-51:01:03/51:01:05-51:01:10/51:01:12/51:01:14-51:01:16/51:01:18-51:01:37/51:01:39-51:01:67/51:01:69-51:02:03/51:02:05-51:03/51:05:01:01-51:05:01:02/51:08:01:01-51:08:04/51:08:06-51:09:03/51:11N-51:12/51:14/51:16-51:21/51:24:03-51:24:05/51:26-51:40:02/51:43/51:48-51:53/51:55/51:57-51:58/51:60-51:61:02/51:63:01/51:65-51:67/51:69-51:77/51:79-51:80/51:82-51:91/51:93-51:102/51:105/51:109-51:111/51:113/51:115/51:117/51:119-51:134/51:137-51:138/51:140-51:147/51:149N-51:150/51:152-51:156/51:158:01-51:167/51:169/51:171/51:173Q-51:175/51:177/51:179-51:189/51:191-51:193:01/51:195-51:196/51:198/51:201-51:212/51:215-51:219/51:221/51:224/51:226/51:228-51:230:01:03/51:232:01-51:240/51:242-51:244/51:246/51:248-51:250/51:252-51:257/51:259-51:280/51:282-51:292/51:294/51:296-51:299/51:301-51:309/51:311-51:318N/51:320-51:326/51:328-51:337/51:339/51:342-51:351/51:354-51:357N/51:360-51:370/51:373-51:377N	37+38+39	B51(5), -, B5102, B5103, Null
B*51:01:04/51:01:13/51:02:04/51:08:05/51:10/51:24:01-51:24:02/51:63:02/51:136/51:168/51:197/51:199/51:295/51:327/51:338/51:358, B*52:01:05/52:01:15/52:06:01-52:06:03/52:97, B*59:02/59:11	37	B51(5), B5102, -, B52(5), B59
B*51:01:11/51:01:17	37+38w+39	B51(5), -
B*51:01:38/51:178N/51:231	37+38	-, Null
B*51:01:68	37+38+39+47	-
B*51:04:01/51:56:03/51:139	7+8+37+38+39	B51(5), -
B*51:04:02/51:46/51:56:01/51:227	7w+8w+37w+38+39	-, B51(5)
B*51:06:01:01-51:06:04/51:41N/51:54/51:59/51:64/51:78:01-51:78:02/51:135/51:190/51:194/51:213/51:247/51:258/51:293/51:319/51:341/51:352-51:353/51:371	38+39	B51(5), -, Null
B*51:07:01-51:07:02/51:214/51:222/51:251/51:340	37+39+40w	B51(5), -
B*51:13:01/51:92:01, B*52:14:01, B*59:03	11+37	-
B*51:13:02/51:92:02/51:172	11+37+38+39	-
B*51:15/51:106:02/51:157/51:225/51:300, B*52:25:01, B*59:01:01:01-59:01:01:02/59:04-59:05/59:07/59:09-59:10N/59:12	31+37	-, B59, Null
B*51:22	37+39+42	-

Allele	Mixes	Serology
B*51:23/51:108	25+37+38+39	-
B*51:44N/51:118N/51:245N	37+39	Null
B*51:56:02	6+7+8+37+38+39	-
B*51:62, A*24:90:02N	31	-, Null
B*51:68/51:176	2+37+38+39	-
B*51:103/51:116/51:151	32+37+38+39	-
B*51:104:01-51:104:02	20+34+37+38+39	-
B*51:106:01/51:223	31+37+38+39	-
B*51:107/51:170	37+38w+39w	-
B*51:112	17+26+37+38+39	-
B*51:114	37w+38+39	-
B*51:193:02	37+38?+39?	-
B*51:200	7+37+38+39	-
B*51:241	37+38+39+48	-
B*51:281	7w+8w+31+37w+38+39	-
B*51:310, B*53:28	8+37+40	-
B*51:372	6+37w+38+39	-
B*52:01:01-52:01:04/52:01:06-52:01:14/52:01:16-52:01:52/52:03-52:05/52:07:01-52:13/52:15-52:19/52:21:01-52:24/52:26-52:50/52:52-52:53/52:55-52:67/52:69-52:72/52:74/52:76/52:78-52:87/52:89N-52:96N/52:98-52:111	37+40	B52(5), -, Null
B*52:02:01-52:02:02	34+37+40	-
B*52:14:02/52:75	11+37+40	-
B*52:20/52:68	40	-
B*52:25:02	31+37+40	-
B*52:51	32+37+40	-
B*52:54	15+37+40	-
B*52:73	37+40w	-
B*52:77	2+37+40	-
B*52:88	8w+37w+40	-
B*53:01:01-53:01:25/53:01:27-53:02/53:05:01-53:06/53:08:01-53:08:02/53:10/53:15-53:16/53:18-53:27/53:29/53:32/53:34-53:35/53:37/53:41-53:43:02/53:45-53:46/53:48N-53:54/53:56-53:61/53:63-53:64/53:66/53:68-53:71	6+7+8+37	B53, -, Null
B*53:01:26	4+6+7+8+37	-

Allele	Mixes	Serology
B*53:04/53:33	6+7+8+11+37	-
B*53:07	7+11	B53
B*53:14/53:47/53:67, C*06:147	6	-
B*53:17:01-53:17:02, B*58:01:05/58:01:29/58:77/58:82/58:116/58:119	8+37	-
B*53:30/53:39	6+7+8+17+37	-
B*53:31	6+7+11	-
B*53:40	7+8+37	-
B*53:62	7+8w+37	-
B*53:72	6+7w+8w+37w	-
B*54:01:01:01-54:01:06/54:01:08-54:02/54:05N/54:08N/54:10/54:13/54:16-54:29/54:31-54:32/54:34-54:35/54:37/54:39/54:43-54:45, B*55:07/55:124	27+28+29	B54(22), -, Null
B*54:01:07/54:38, B*55:01:01:01-55:01:03/55:01:05-55:01:12/55:01:14-55:02:10/55:02:12-55:03/55:05/55:09-55:12/55:15-55:16/55:19/55:21-55:22/55:24-55:26/55:29-55:31/55:33-55:34/55:36-55:48/55:50/55:52/55:54-55:55N/55:57/55:59-55:60/55:62-55:63/55:65-55:67/55:69-55:77/55:79-55:80/55:82-55:85/55:87-55:89N/55:91-55:93/55:95-55:102/55:104-55:109/55:111-55:117N/55:119-55:120/55:125N-55:129/55:131-55:133, B*56:01:13/56:10/56:19N/56:66	27+28	-, B55(22), B22, Null
B*54:03/54:33/54:40	27+28+29+30	-
B*54:04/54:11/54:42	28+29	-, B54(22)
B*54:06, B*56:22/56:88	28w+30	-
B*54:07	27+28w+29	-
B*54:09	28w	-
B*54:12	28+29+31+37	-
B*54:14	28w+29	-
B*54:15, C*03:464	29	-
B*54:30	27+29	-
B*54:36	27w+28w+29	-
B*54:41/54:46	27w+29	-
B*55:01:04	27+28+39+41	B55(22)
B*55:02:11	27+39+41	-
B*55:08/55:51, B*56:07/56:15/56:46	28+30	B56(22), -
B*55:14	21+28	-
B*55:18/55:61/55:78/55:94/55:118N, B*56:08/56:17/56:23/56:42/56:44/56:54/56:89/56:91, 01:90/01:136/01:145:01N, 03:19/03:307/03:318N, 07:289/07:536/07:725, 08:207, 12:03:66/12:15/12:208, 15:07:01:01-15:07:01:02/15:43/15:85/15:144	27	-, Null, Cw3

Allele	Mixes	Serology
B*55:23/55:32/55:68, B*56:18/56:31-56:32	21+28w	-
B*55:27	24w+28w	-
B*55:28:02	28+39+41	-
B*55:35	27+28+34+43	-
B*55:53	27w+28w	-
B*55:56	23+27+28	-
B*55:64	27+28w	-
B*55:103	28+31+37	-
B*55:122	25+27+28	-
B*55:123	12+13+27+28	-
B*55:130	5+27+28	-
B*55:134, B*56:01:01:01-56:01:12/56:01:14-56:02:03/56:04:01:01-56:04:04/56:16/56:20:01-56:20:02/56:24-56:30/56:34-56:38N/56:40-56:41/56:43/56:47-56:49/56:51-56:53/56:55:01:02-56:59/56:61-56:64/56:68-56:73/56:75-56:87/56:90	27+28+30	-, B56(22), Null
B*56:05:01-56:05:02	28+30+39+41	B56(22)
B*56:06	5+28+30+39+41	B78
B*56:09	7+19+21+28w+30	-
B*56:11	7+19+28w+30	B56(22)
B*56:12	7w+19w+28	B55(22)
B*56:13, B*82:04	23+27+28+30	B56(22), -
B*56:14, C*03:102	27+33	-
B*56:21	28+30+37+39	-
B*56:39	27+28+30+34+43	-
B*56:45	5+27+28+30	-
B*56:60	7+19+27+28	-
B*56:65	28+30+31+37	-
B*56:74	25+27+28+30w	-
B*57:01:01:01-57:01:01:27/57:01:03-57:01:18/57:01:20-57:01:32/57:01:34-57:01:46/57:04:01-57:04:03/57:06/57:08/57:10/57:14:02-57:15/57:18-57:21/57:23-57:27/57:29-57:31:02/57:33-57:38/57:40-57:41/57:43-57:44/57:47-57:50/57:52/57:54-57:56/57:58-57:62/57:64-57:65/57:67:01-57:68/57:71-57:75/57:77-57:79N/57:81/57:85-57:92/57:95/57:97-57:100/57:102-57:106/57:108-57:117/57:121-57:123/57:125-57:130N/57:132-57:136/57:138-57:147/57:149-57:154/57:156-57:160/57:162N-57:165	34+35	B57(17), -, Null
B*57:01:02/57:01:33/57:02:01:01-57:03:05/57:09/57:16-57:17/57:28N/57:32/57:39/57:42/57:46/57:66/57:70/57:80/57:84/57:94/57:96/57:101/57:107/57:120/57:124/57:137/57:148/57:155	35	B57(17), -, Null

Allele	Mixes	Serology
B*57:01:19/57:14:01	34+35+36	-
B*57:05	8+35+37	-
B*57:07	35+44	-
B*57:12	35+43	-
B*57:13/57:22/57:161	34+35+45	-
B*57:45/57:51/57:69	17+35	-
B*57:53	34w+35	-
B*57:57/57:63	35+45	-
B*57:76	11+34+35	-
B*57:93/57:119	34+35+44	-
B*57:131	8+34+37	-
B*58:01:01:01-58:01:04/58:01:06-58:01:15/58:01:17-58:01:19/58:01:21-58:01:28/58:01:30-58:01:44/58:04-58:05/58:09-58:13/58:17N/58:21-58:24/58:28:01-58:29/58:31N-58:33/58:35/58:37/58:39N-58:42/58:45:01-58:45:02/58:47-58:59:02/58:61-58:63/58:65-58:72N/58:74/58:76/58:79-58:81/58:83-58:90/58:92/58:94N-58:103/58:105-58:110/58:112-58:113/58:115/58:117-58:118/58:120/58:123-58:131/58:133N-58:134/58:136/58:139-58:140	8+36+37	B58(17), -, Null
B*58:01:16/58:02:01:01-58:02:02/58:07/58:15-58:16:02/58:18/58:25-58:26/58:38/58:43/58:60/58:104/58:121/58:132/58:137	36	-, B58(17)
B*58:01:20/58:75/58:93N/58:114	8w+36+37w	-, Null
B*58:06/58:46	36+45	-
B*58:08:01-58:08:02	36+37+45	B17, -
B*58:14	35+36	-
B*58:19/58:78/58:91	8+36+37+45	-
B*58:20	11+36	-
B*58:34	8+11+36+37	-
B*58:36	8+34+36+37	-
B*58:44	8+19+21+36+37	-
B*58:64	5+36	-
B*58:73	8+31+36+37	-
B*58:111	8+17+36+37	-
B*58:122	36+37	-
B*58:135/58:138	8+36?+37	-
B*59:06	29+31+37	-
B*67:01:01-67:01:03/67:03-67:07	12+13+28	B67, -

Allele	Mixes	Serology
B*67:02:01:01-67:02:01:02	12+33	-
B*73:01:01:01-73:03	48	B73, -
B*78:03	5+39+41+42	-
B*78:05	5+40+41	-
B*78:06	40+41	-
B*78:09-78:10	38+39+41	-
B*82:01:01:01-82:02:02	23+25+27+28+30w	B82, -
B*82:03	23+25+27+28+30	-
B*83:01	7w+19w+23+24w+25+30w	-
A*26:89, C*03:278/03:504, C*06:77	32	-
C*03:342	19+24w	-
C*03:542	30	-
C*06:02:101	4?	-
C*07:463	2+29w+48w	-
C*07:758	2+29+48	-
C*08:93, C*14:61	25	-
C*12:301	4+15w	-
C*14:92	29+33+41	-
C*15:71	8+31+37	-

w = weak

? = nucleotide sequence information not available for the primer matching sequence or alleles with unknown reactivities

Table 6: Amplification patterns of **HLA-DRB1*** and other HLA-DRB* alleles detected by the HLA-DRB1* CTS-PCR-SSP primer mixes (**Lot No. DR43-3**) based on IMGT/HLA Sequence Database Release 3.51.0, January 2023

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*01:01:01:01-01:01:01:30/01:01:32-01:01:44/01:01:46-01:02:11/01:02:13-01:17/01:19-01:20:02/01:23-01:32/01:35-01:38/01:41-01:43/01:45-01:48/01:50-01:59/01:61-01:64/01:66-01:74/01:77-01:133/01:135-01:141	DR1, -, DR103, Null	1																						
DRB1*01:01:45	-	?																						
DRB1*01:21	-	1							w															
DRB1*01:34	-	1										w												
DRB1*01:39N-01:40N/01:49/01:65:01-01:65:02	Null, -	w																						
DRB1*01:75	-	1											12											
DRB1*03:01:01:01-03:01:01:25/03:01:03/03:01:05-03:01:06/03:01:08-03:01:13/03:01:15-03:01:41/03:04:01-03:04:02/03:06/03:11:01/03:11:03/03:13:01/03:15:01-03:15:01:02/03:16/03:18-03:20/03:22-03:23/03:25:01-03:26/03:28/03:30-03:34/03:36-03:37/03:39/03:43-03:45/03:47-03:48/03:50-03:52:02/03:55-03:56/03:58/03:60-03:62/03:64/03:66-03:73/03:75/03:77-03:78/03:80/03:82-03:83/03:85-03:86/03:89/03:91/03:93-03:96/03:98-03:99/03:104/03:106-03:110/03:112/03:114/03:116-03:118/03:121/03:123-03:124/03:128-03:129/03:132-03:134/03:136-03:139/03:141-03:144/03:146-03:153/03:159-03:161/03:163-03:166/03:168-03:170/03:172-03:173/03:175/03:180-03:184/03:186/03:189N/03:192-03:204N	DR17(3), -, DR3, Null				4	5												17						
DRB1*03:01:02:01-03:01:02:02/03:13:02/03:15:02/03:145	DR17(3), -				w	5												17						
DRB1*03:01:04/03:76/03:84/03:92/03:174N, DRB1*13:129/13:326	DR17(3), -, Null					5												17						
DRB1*03:01:07/03:46/03:59/03:79	DR17(3), -				4	5												w						
DRB1*03:01:14/03:127/03:135/03:190	-				w	5												w						
DRB1*03:02:01:01-03:02:01:02/03:02:03-03:02:06/03:05:01-03:05:02/03:09/03:14/03:27/03:29/03:53/03:74/03:81/03:88/03:90/03:102-03:103/03:115/03:122/03:131/03:154/03:167/03:176/03:178-03:179/03:188	DR18(3), -, DR3				4		6											17						
DRB1*03:02:02/03:05:03	DR18(3), DR3				w		6											17						
DRB1*03:03/03:07:01:01-03:07:02/03:17/03:21/03:24/03:35/03:40-03:41:02/03:57/03:63/03:101/03:105/03:113/03:119-03:120/03:171/03:185	DR18(3), DR3, -				4													17						
DRB1*03:08/03:65/03:140	-				4	5							12											
DRB1*03:10/03:191	DR17(3), -				4	5												17	18					
DRB1*03:12	DR3				4	?										15		17						
DRB1*03:38	-				4		6											w						

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*03:42/03:87, DRB3*01:34, DRB3*03:33	-					5																21		
DRB1*03:49	-				4			7										17						
DRB1*03:54/03:100:01-03:100:02/03:205	-				4	w												17						
DRB1*03:97, DRB1*11:293, DRB1*13:45/13:76/13:197/13:206, DRB1*14:01:01:01-14:01:05/14:05:01:01-14:05:05/14:07:01-14:08/14:14/ 14:16/14:22-14:23:04/14:25:01-14:26/14:32:01-14:32:03/14:35-14:36/ 14:38:01-14:38:02/14:42-14:45/14:49/14:53-14:54:01:17/14:54:03-14:54:07/ 14:54:09-14:56/14:58/14:60/14:62/14:69-14:70/14:72/14:74-14:75/14:86- 14:88/14:90-14:92N/14:96-14:97/14:99/14:101/14:103/14:105/14:110/ 14:112-14:114/14:117-14:119/14:122/14:124/14:126:01-14:126:02/14:128- 14:131/14:137N/14:139-14:140/14:142/14:146-14:147/14:149-14:151/14:153/ 14:155/14:157-14:158/14:160/14:162-14:164/14:166N-14:168/14:172-14:173/ 14:175/14:178/14:184-14:193/14:197N/14:202-14:204/14:205:02/14:207- 14:208/14:211/14:213-14:217/14:219/14:221/14:225/14:227-14:228/14:230/ 14:232/14:234-14:238/14:240/14:243-14:245/14:251	-, DR14(6), DR6, DR13(6), Null																	17	18					
DRB1*03:111/03:162, DRB1*09:06	-					5																		
DRB1*03:125, DRB1*13:27:01-13:27:02/13:176	-, DR13(6)					5									14			17						
DRB1*03:126, DRB1*14:95/14:132	-					5												17	18	19				
DRB1*03:130	-				w		6											w						
DRB1*03:155	-				4											15		17						
DRB1*03:156N	Null				4	5																		
DRB1*03:157	-					5												17		w				
DRB1*03:158	-				4	5										15		17						
DRB1*03:177	-					5												17		19				
DRB1*03:187	-				4	5												?						
DRB1*04:01:01:01-04:05:11/04:05:13-04:11:03/04:11:05-04:14/04:16-04:60/ 04:62-04:137/04:139-04:144/04:146-04:153/04:155/04:157N-04:169/04:171- 04:210/04:212N-04:264N/04:266N-04:290/04:292-04:322/04:324-04:338/ 04:340-04:361, DRB1*14:10:01:01-14:10:01:02/14:57/14:143/14:182	DR4, -, Null, DR14(6)							7																
DRB1*04:15/04:61/04:138/04:145/04:154/04:170/04:265/04:291/04:323/ 04:339, DRB1*11:102:01-11:102:02/11:175/11:275/11:278	DR4, -							7					12											
DRB1*04:211	-							7												19				
DRB1*07:01:01:01-07:01:01:02/07:01:01:04-07:01:30/07:03-07:139/07:141- 07:146	DR7, -, Null								8															
DRB1*07:140	-								?															

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*08:01:01:01-08:01:02:08:01:04-08:01:05:08:01:07-08:02:04/ 08:03:02:01-08:04:01:05:08:04:03-08:07:08:09:01-08:19:08:21-08:24:02/ 08:26-08:28:08:30:01-08:30:03:08:33:08:35-08:37:08:39-08:40:08:42-08:46/ 08:48-08:51:08:54-08:56:08:58-08:62:08:64-08:66:08:68:01-08:68:02:08:70- 08:74:08:76-08:79:08:81-08:86:08:88-08:108:08:110-08:111:08:114-08:119, DRB1*14:15	DR8, -, Null									9								17						
DRB1*08:02:05	-									?								?						
DRB1*08:08/08:29, DRB1*14:04:01:01-14:04:09/14:28/14:31/14:50:01- 14:50:02/14:68:01-14:68:02/14:71/14:73/14:76/14:79/14:93/14:107/14:120/ 14:138/14:145/14:152N/14:180/14:196/14:201/14:205:01/14:206/14:210Q/ 14:220/14:224/14:229/14:233/14:241/14:250	-, DR1404, DR14(6), Null									9								17	18					
DRB1*08:20/08:67, DRB1*11:264, DRB1*13:07:01-13:09/13:11:01-13:11:02/ 13:14:01-13:14:03/13:17/13:19/13:22:01:01-13:23:02/13:25/13:37/13:44/ 13:46-13:47/13:50:01/13:52/13:57/13:60/13:62/13:64/13:68/13:70/13:72/ 13:82-13:84/13:86/13:97:01-13:97:02/13:106-13:107/13:114/13:116/13:132/ 13:135/13:150/13:172/13:175/13:178/13:180/13:182/13:192/13:204/13:224/ 13:229/13:231-13:232/13:246/13:250/13:278Q/13:284/13:286/13:314/ 13:316/13:330/13:335, DRB1*14:20-14:21/14:34/14:37/14:40-14:41/14:77/ 14:84/14:94/14:100/14:102/14:116/14:123/14:127:01-14:127:02/14:136/ 14:156/14:159/14:181/14:183	-, DR13(6), DR14(6)																	17						
DRB1*08:25/08:47, DRB1*12:16:01-12:16:03/12:22/12:39:01-12:39:02	-									9				13				w						
DRB1*08:31, DRB1*11:01:01:01/11:01:01:03-11:02:02/11:02:04-11:02:07/ 11:04:01:01-11:04:12/11:04:14-11:06:03/11:08:01-11:08:03/11:10:01- 11:10:02/11:12:01-11:15:02/11:18-11:19:03/11:21/11:23:01-11:25/11:27:01- 11:28:01:02/11:28:03-11:35/11:37:01-11:39/11:42:01-11:51/11:56-11:57/ 11:58:02/11:60-11:62:02/11:64-11:66:02/11:69-11:70/11:72-11:75/11:77- 11:82/11:84:01-11:84:03/11:88/11:90-11:93/11:95-11:101:03/11:104/11:106/ 11:108-11:109/11:111-11:112/11:114/11:116-11:121/11:124:01-11:124:02/ 11:126-11:134/11:137/11:140-11:149/11:152-11:155/11:157-11:158/11:160/ 11:162-11:166/11:169N-11:171/11:174/11:177-11:181/11:183-11:184/11:186/ 11:188-11:191/11:194:01-11:201/11:204-11:207/11:209-11:210/11:212/ 11:214-11:230/11:232/11:234-11:236/11:238-11:240/11:242-11:253/11:255- 11:259/11:261-11:263/11:265-11:266/11:268-11:274/11:276-11:277/11:279/ 11:281-11:282/11:284-11:292/11:295-11:300/11:302-11:312	DR11(5), -, Null												12											
DRB1*08:32/08:53/08:87, DRB1*12:46	-									9				13				17						
DRB1*08:34/08:52/08:109/08:112-08:113	-									9								w						
DRB1*08:38	-							7		9								17						

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*08:41/08:63/08:75, DRB1*11:67/11:193:01-11:193:02	-									9			12											
DRB1*08:69	-									9														
DRB1*09:01:02:01-09:05/09:08/09:10-09:22/09:24-09:52N	DR9, -, Null										10													
DRB1*09:07	-										10										20			
DRB1*10:01:01:01/10:01:01:03-10:11/10:13-10:22/10:24-10:30/10:32-10:45	DR10, -											11												
DRB1*10:23	-											11	12											
DRB1*10:31	-											w												
DRB1*11:03:01-11:03:04/11:11:01/11:11:03/11:36/11:41/11:55/11:63:01-11:63:02/11:68/11:76/11:85/11:138-11:139/11:151/11:161/11:168/11:176/11:213/11:231/11:233/11:237/11:241/11:260/11:301N	DR11(5), -, Null												12			15								
DRB1*11:07:01/11:103:01-11:103:02/11:105/11:107/11:125/11:136/11:173	DR11(5), -				4								12											
DRB1*11:07:02	-				w								12											
DRB1*11:09:01:01-11:09:01:02/11:87	DR11(5), -					6							12		14		16			19				
DRB1*11:16:01:01-11:16:01:02	DR11(5)/ DR13 (6), -												12		14									
DRB1*11:17/11:52/11:54:01-11:54:02/11:89/11:150/11:159/11:202/11:267/11:283, DRB1*14:11:01:01-14:11:01:02	-, DR11(5), DR14(6)												12						18					
DRB1*11:20	DR11(5)					6							12		14									
DRB1*11:22	-							7					?											
DRB1*11:26	DR11(5)												?											
DRB1*11:28:02	-					6							12		14		16			w				
DRB1*11:40	-												?		w	15	w							
DRB1*11:53	-				4								12						18					
DRB1*11:58:01/11:115	-												12		14		16			w				
DRB1*11:59	-					5							12		w	15	w							
DRB1*11:83/11:182	-					5							12		14		16			19				
DRB1*11:113	-												12		14		16			19				
DRB1*11:122	-						6						12		w	15	w							
DRB1*11:135	-					5							12											
DRB1*11:203/11:211/11:254	-												w											
DRB1*11:280	-												12							19				

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*12:01:01:01-12:01:04/12:01:06-12:01:07/12:01:09-12:01:11/ 12:02:01:01-12:02:05/12:02:07-12:03:03/12:05-12:08/12:10-12:15/12:17/ 12:19:01-12:19:02/12:21/12:23/12:26-12:30/12:32/12:34-12:38/12:40-12:41/ 12:43-12:45/12:50-12:53/12:55-12:56/12:59-12:61/12:63-12:65/12:68-12:75/ 12:77-12:82/12:84-12:88/12:90-12:95/12:97-12:103	DR12(5), -, Null													13				w						
DRB1*12:01:05/12:20	-													w				w						
DRB1*12:01:08/12:02:06/12:18/12:33/12:47	-													13										
DRB1*12:01:12	-													13				?						
DRB1*12:04/12:49	-												12	13										
DRB1*12:09/12:24N-12:25/12:48/12:54/12:62/12:67/12:76	-, Null													13				17						
DRB1*12:58	-													13				w		19				
DRB1*12:66	-													13				w	18					
DRB1*12:89	-													?				?						
DRB1*12:96	-									9				13				17	18					
DRB1*13:01:01:01-13:01:40/13:02:04/13:10/13:15-13:16:01:02/13:28:01- 13:28:02/13:35/13:40/13:51/13:53/13:59/13:61:01-13:61:02/13:69/13:79- 13:80/13:87/13:91-13:92/13:98/13:102/13:105/13:110/13:112-13:113N/ 13:117:01:01-13:117:01:02/13:121/13:125/13:127/13:130-13:131/13:137N- 13:138/13:141/13:148/13:153/13:166/13:173/13:177/13:184-13:187/13:190/ 13:201/13:205/13:213-13:215/13:221-13:222/13:226/13:233/13:238/ 13:242:01-13:242:02/13:245/13:251-13:252N/13:254/13:256/13:258-13:263/ 13:266-13:267/13:270-13:274/13:279/13:281/13:283/13:287/13:290-13:291/ 13:294/13:300/13:303-13:305/13:308-13:309/13:311/13:313/13:315/13:322N- 13:324N/13:327-13:328	DR13(6), -, Null														14			17						
DRB1*13:01:41/13:160/13:218	-														14			w						
DRB1*13:02:01:01-13:02:03/13:02:05-13:02:11/13:02:13-13:02:22/13:31/ 13:34/13:36/13:39/13:41/13:73-13:74/13:85/13:96:01-13:96:02/13:99/ 13:104/13:109/13:123-13:124/13:126/13:128/13:140/13:142N-13:143/ 13:145/13:147/13:155/13:168/13:170/13:171:02/13:198/13:202/13:207- 13:212/13:220/13:225/13:236-13:237/13:239-13:241/13:243-13:244/13:255N/ 13:257/13:264-13:265/13:268N-13:269/13:275/13:280/13:282/13:285/13:288- 13:289N/13:292-13:293/13:297-13:298N/13:301/13:306/13:310N/13:318/ 13:320/13:325/13:332/13:334	DR13(6), -, Null					6									14			17						
DRB1*13:02:12/13:103/13:165	-					6									14			w						

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*13:03:01:01-13:04/13:12:01-13:13/13:20:01:01-13:21:02/13:24:01:01-13:24:01:03/13:30/13:33:01-13:33:03/13:38/13:48-13:49:02/13:55/13:58/13:66:01-13:66:02/13:75/13:78/13:81/13:88-13:90/13:94:01-13:95/13:101/13:108/13:115/13:118/13:122/13:133-13:134/13:149/13:151-13:152/13:154/13:161/13:164/13:167/13:169/13:174/13:194/13:216/13:219/13:227/13:230/13:247/13:253/13:276/13:295N/13:307/13:317/13:319N/13:321/13:329N/13:333, DRB1*14:65	DR13(6), -, Null															15		17						
DRB1*13:05:01:01-13:05:03/13:26:01-13:26:02/13:136/13:158:01-13:158:02/13:203/13:277, DRB1*14:27:01-14:27:02	DR13(6), -, DR14(6)					6									14		16	17	19					
DRB1*13:06/13:146/13:156	DR13(6), -														14			17	19					
DRB1*13:18/13:42:01-13:42:02	DR13(6), -														14		16	17	19					
DRB1*13:29:01:01-13:29:01:02	DR13(6), -					6										15		17						
DRB1*13:32/13:93	-														14	15		17						
DRB1*13:43/13:188, DRB1*14:169	-															15		17	18					
DRB1*13:50:02-13:50:03/13:189	-					6									14		16	17	w					
DRB1*13:54	DR14(6)															15		w						
DRB1*13:56/13:119/13:235, DRB1*14:02:01:01-14:02:07/14:02:09-14:03:02/14:09/14:30/14:47/14:89/14:115/14:135/14:144/14:174/14:179/14:194/14:198/14:209/14:218/14:223/14:239/14:242:01-14:242:02	-, DR14(6), DR1403					6												17	19					
DRB1*13:63	-					6									w	15	w	17						
DRB1*13:65/13:120/13:139	-					6									14	15		17						
DRB1*13:67	-					6									14									
DRB1*13:71	-				5											15		17						
DRB1*13:77, DRB1*14:177/14:195N	-, Null																	w	19					
DRB1*13:100/13:162	-																	w						
DRB1*13:111	-														w	15	w	17						
DRB1*13:144	-				5										14		16	17	19					
DRB1*13:157/13:331, DRB1*14:67	-					6									14			17	19					
DRB1*13:159/13:171:01/13:179	-					6									14			17	18					
DRB1*13:163/13:181/13:193, DRB1*14:02:08/14:06:01-14:06:06/14:12:01-14:12:02/14:17/14:29/14:52/14:80-14:81/14:83/14:108/14:121/14:165/14:176/14:212	-, DR14(6), DR6																	17	19					
DRB1*13:183/13:200N/13:248-13:249N/13:299	-, Null														14									
DRB1*13:191/13:217	-														14			17	18					

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*13:196/13:312, DRB1*14:19/14:24/14:109/14:226	-, DR14(6)						6											17						
DRB1*13:199	-						6								14		16	w		19				
DRB1*13:223/13:234	-														w			17						
DRB1*13:228	-															15								
DRB1*13:296	-						6								w			17						
DRB1*13:302	-						6								14		16	17	18	19				
DRB1*14:13/14:85	DR14(6), -						6									15		17		19				
DRB1*14:18/14:33/14:59/14:64/14:106/14:133/14:154/14:246	DR6, -																	17	18	19				
DRB1*14:39/14:54:02/14:61/14:199	-																		18					
DRB1*14:46/14:222N	-, Null						6													19				
DRB1*14:48	-						6											w		19				
DRB1*14:51/14:134/14:200	-						6											17	18	19				
DRB1*14:54:08	-																	w	18					
DRB1*14:63	-						6								14	15		17		19				
DRB1*14:78	-														14	15		17		19				
DRB1*14:82	-					5												17	18					
DRB1*14:98	-						6											17		w				
DRB1*14:104/14:111	-				w													17	18					
DRB1*14:125	-																	w	w					
DRB1*14:141, DRB3*01:01:02:01-01:07/01:09-01:13/01:15-01:33/01:35-01:41/01:43-01:45/01:47-01:51/01:53-01:67/01:69-01:70/01:72-01:94/01:96-01:103/01:105N-01:106/01:108/01:110-01:112, DRB3*02:01-02:05:02/02:07-02:17/02:19/02:21-02:25/02:28-02:46/02:48-02:57/02:59-02:65/02:67N-02:77/02:79-02:90/02:92-02:108/02:110/02:112-02:116/02:118-02:123/02:125N-02:130/02:132-02:142/02:144-02:159/02:161/02:163-02:189, DRB3*03:01:01-03:01:12/03:01:14-03:21/03:23-03:32/03:34-03:45/03:47-03:61	-, DR52, Null																					21		
DRB1*14:148	-									9									18					
DRB1*14:161	-							7										17	18					
DRB1*14:170	-															15		17		19				
DRB1*14:171	-				4													17	18					
DRB1*14:248	-																	?	?					
DRB1*14:249	-									9								17	?					

Allele	Serology	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23
DRB1*15:01:01:01-15:01:21/15:01:23-15:20/15:22-15:33/15:35-15:98/15:100-15:108/15:110-15:112/15:114-15:216	DR15(2), -, DR2, Null		2																					
DRB1*15:21, DRB1*16:09:01-16:10:02/16:33/16:36-16:37/16:58/16:66	-, DR16(2)		2	3																				
DRB1*15:99/15:109	-		2										12											
DRB1*15:113N	Null		w																					
DRB1*16:01:01:01-16:01:08/16:01:10-16:05:02/16:07/16:11-16:15/16:17-16:26/16:28-16:32/16:34-16:35/16:38:01-16:49/16:51-16:57/16:59Q-16:65/16:67-16:73	DR16(2), -, DR2, Null			3																				
DRB1*16:01:09	-			w																				
DRB1*16:08	-			3		6									14		16			w				
DRB1*16:16/16:27	-			3									12											
DRB3*01:08/01:46/01:52/01:71, DRB3*02:20	-					6																21		
DRB3*01:14	-				w													w						
DRB3*01:42, DRB3*02:18	-												12									21		
DRB3*01:107, DRB3*02:190	-																					?		
DRB3*02:06	-					?																21		
DRB3*02:78/02:109N/02:117/02:124/02:162, DRB3*03:01:13/03:46	-, Null																					w		
DRB4*01:01:01:01/01:01:01:03-01:01:01:04/01:02-01:03:01:01/01:03:01:03-01:03:01:12/01:03:01:14-01:03:02:02/01:152, DRB4*02:01N	DR53, -, Null																						22	
DRB4*01:03:01:02N/01:03:01:13N	Null																							23
DRB5*01:01:01:01-01:01:09/01:01:11-01:02:01/01:03-01:12/01:14-01:25/01:27N-01:40/01:42-01:47/01:49N-01:50/01:52N-01:72/01:75-01:78/01:80-01:86/01:88-01:96/01:98-01:109/01:111-01:112/01:114-01:130, DRB5*02:02:01-02:37	DR51, -, Null																				20			
DRB5*01:02:02/01:26/01:48N/01:74/01:113	-, Null																				w			
DRB5*01:13/01:41	-												12								20			

w = weak

? = nucleotide sequence information not available for the primer matching sequence or alleles with unknown reactivities (HLA-DRB1*08:04:02 is not listed)

14. Certificates of Analysis



HEIDELBERG
UNIVERSITY
HOSPITAL

Heidelberg University Hospital
Institute of Immunology
Transplantation Immunology
Im Neuenheimer Feld 305
69120 Heidelberg – Germany
Phone: +49 6221 564013
Fax: +49 6221 564200
www.ctstransplant.org

CTS Collaborative Transplant Study

Certificate of Analysis

HLA-A* CTS-PCR-SSP Tray Kit

Product number	101
Lot number	A17-2
Number of HLA-specific primer mixes per test	24

Mix specifications

The specificity of each primer pair has been tested against a panel of well characterized DNAs.

Result

No false positive or false negative amplifications were obtained under our test conditions of the bulk reagents.

Date of approval	20.03.2023
Approved by	H. Tran, M.D. Quality Control, Supervisor



HEIDELBERG
UNIVERSITY
HOSPITAL

Heidelberg University Hospital
Institute of Immunology
Transplantation Immunology
Im Neuenheimer Feld 305
69120 Heidelberg – Germany
Phone: +49 6221 564013
Fax: +49 6221 564200
www.ctstransplant.org

CTS Collaborative Transplant Study

Certificate of Analysis

HLA-B* CTS-PCR-SSP Tray Kit

Product number	102
Lot number	B16-2
Number of HLA-specific primer mixes per test	48

Mix specifications

The specificity of each primer pair has been tested against a panel of well characterized DNAs.

Result

No false positive or false negative amplifications were obtained under our test conditions of the bulk reagents.

Date of approval	20.03.2023
Approved by	H. Tran, M.D. Quality Control, Supervisor



HEIDELBERG
UNIVERSITY
HOSPITAL

Heidelberg University Hospital
Institute of Immunology
Transplantation Immunology
Im Neuenheimer Feld 305
69120 Heidelberg – Germany
Phone: +49 6221 564013
Fax: +49 6221 564200
www.ctstransplant.org

CTS Collaborative Transplant Study

Certificate of Analysis

HLA-DRB1* CTS-PCR-SSP Tray Kit

Product number	104
Lot number	DR43-3
Number of HLA-specific primer mixes per test	23

Mix specifications

The specificity of each primer pair has been tested against a panel of well characterized DNAs.

Result

No false positive or false negative amplifications were obtained under our test conditions of the bulk reagents.

Date of approval	24.04.2023
Approved by	H. Tran, M.D. Quality Control, Supervisor